

Excellence in **Digital Healthcare**

DECEMBER 2020

ABOUT **CLINICAL SOLUTIONS APPLICATIONS**

IMAGE-COM **MEASUREMENTS 3D OPTION IMAGE-COM** AutoLV AutoLA AutoIMT STRESS ECHO AutoStrain LV/RV/LA **2D STRAIN FETAL 2D STRAIN 4D LV-ANALYSIS 4D RV-FUNCTION 4D MV ASSESSMENT 4D CARDIO-VIEW**

TOMTEC ACADEMY IMPRESSUM

CATH OPTIONS

TOMTEC

TOMTEC Imaging Systems GmbH is a global leader in medical imaging software solutions. Specializing in ultrasound, TOMTEC offers state-of-the-art solutions for our clinical customers and industry partners.

TOMTEC was founded in 1990 and our Headquarters are located in Munich, Germany. Our products encompass a wide range of 2D and 3D/4D technologies for visualization, automated image analysis, quantification, reporting and image management.

TOMTEC proudly offers solutions for adult and pediatric cardiology, obstetrics, gynecology, radiology and vascular diagnostics. We partner with industry experts to deliver education and best practice sharing via TOMTEC ACADEMY.

At TOMTEC, we are passionate about providing worldclass support and service for all of our customers.

TOMTEC stands for "Excellence in Digital Healthcare". Our mission and core value is a commitment to continuous improvement, delivery of excellent image and analysis quality, ease-of-use, reproducibility through automation, world-class service, support and education.

TOMTEC focuses on emerging, future technologies like Cloud-based Software Solutions, Artificial Intelligence and Virtual Reality to support healthcare providers in delivering best possible patient care.

The company maintains close working relationships with many leading universities and research institutes around the world.

ABOUT **CLINICAL SOLUTIONS APPLICATIONS**

IMAGE-COM MEASUREMENTS 3D OPTION IMAGE-COM AutoLV **AutoLA** AutoIMT STRESS ECHO AutoStrain LV/RV/LA **2D STRAIN** FETAL 2D STRAIN **4D LV-ANALYSIS 4D RV-FUNCTION 4D MV ASSESSMENT 4D CARDIO-VIEW**

CATH OPTIONS



TOMTEC ARENA INCREASE EFFICIENCY & DIAGNOSTIC QUALITY

Solution for all your needs in Cardiovascular Imaging

WORKFLOW INTEGRATION MULTI-MODALITY VENDOR NEUTRAL AUTOMATION CLINICAL USE

ABOUT CLINICAL SOLUTIONS APPLICATIONS

IMAGE-COM MEASUREMENTS 3D OPTION IMAGE-COM AutoLV AutoLA AutoIMT STRESS ECHO AutoStrain LV / RV / LA 2D STRAIN FETAL 2D STRAIN 4D LV-ANALYSIS 4D RV-FUNCTION 4D MV ASSESSMENT 4D CARDIO-VIEW

CATH OPTIONS

DIAGNOSTIC IMAGING

Efficient multi-

modality image

to clinical tools

review with access

O'NOSTIC O'NOSTIC

1.1100.001

L



Review of prior studies

CLINICAL AREAS \checkmark MV Repair/Replacement - ASD/VSD Closure CLINICAL NEEDS \checkmark Clinical guidance to determine time of intervention for asymptomatic vs. symptomatic patients Comprehensive pre-surgical assessment for device selection and communication with surgeons Verify cardiac CT findings Complete ASD/VSD assessment (Size / Shape / Location)

Increase diagnostic confidence

STRUCTURAL HEART DISEASE

RISK MANAGEMENT

Monitor subclinical markers for left ventricular dysfunction

THERAPY MONITORING

Track clinical

markers to assess

patient outcome

CLINICAL AREAS

- Cardio-Oncology
- Diabetes
- Pulmonary Hypertension (early stage)
- Coronary Artery Disease
- Valvular Heart Disease

CLINICAL NEEDS

- Subclinical marker for LV Dysfunction
- Consistent results
- High reproducibility
- Solutions tailored to individual needs

CLINICAL AREAS

- Heart Failure
- AV Optimization
- Pulmonary Hypertension

CLINICAL NEEDS

- Consistent results independent from vendor

 \checkmark

 \checkmark

- Clinical marker to monitor response to treatment
- Changes in EF
- Changes in End-Diastolic Volume
- Changes in GLS
- High reproducibility through automation

ABOUT **CLINICAL SOLUTIONS APPLICATIONS**

IMAGE-COM

MEASUREMENTS

3D OPTION IMAGE-COM

AutoLV

AutoLA

AutoIMT

STRESS ECHO

AutoStrain LV/RV/LA

2D STRAIN

FETAL 2D STRAIN

4D LV-ANALYSIS

4D RV-FUNCTION

4D MV ASSESSMENT

4D CARDIO-VIEW

CATH OPTIONS

DIAGNOSTIC IMAGING



Efficient study review

TOMTEC's solutions have been designed with the needs of a busy cardiovascular laboratory in mind. With the ability to quickly retrieve and navigate through studies, clinicians are now able to efficiently interpret patient data, while having access to all relevant clinical tools. Our software enables clinicians to review images from all diagnostic modalities, including previous studies to cross examine clinically warranted findings. All measurements that are performed on the cart are imported and populated to worksheets or reports.

In addition, the user has the ability to perform any echo, vascular or cath lab measurements offline. Both the tools for clinical practice and advanced analysis are vendor neutral and seamlessly integrated into your workflow.

DIAGNOSTIC QUALITY



IMAGE-COM¹

2D and 3D multi-

with automated

quantification

CARDIAC MEASUREMENTS¹ modality image viewer Complete echocardiographic measurement package with labeled measurements and access to advanced M-Mode, Doppler and 2D measurements



Cath-QCA¹ Clinical solution for quantitative coronary analysis based on an automated contour detection of coronary angiograms



DIGITAL SUBTRACTION ANGIOGRAPHY (DSA)¹ Enhanced display of peripheral vessels with automatic image summation over time and compensation of motion artifacts







AutoLV¹ "Two click" Biplane left ventricular volume quantification with automated contour proposals for





3D Quantification of left ventricular volume and function

function

Monitoring of high risk patients

Patients that are currently treated for conditions such as HTN, CAD, diabetes or undergoing chemotherapy are at risk to develop heart failure. The ability to closely monitor subclinical markers is essential to detect even subtle changes in cardiac function. This allows clinicians to properly manage patients in order to minimize the risk of developing heart failure. Accurate measurements that are not impacted by user variability are needed. Automated strain imaging solutions provide highly eproducible data which can be used as subclinical markers to onitor cardiac function before visual changes are observed. MTEC has a portfolio of automated 2D and 3D strain imaging itions to support clinicians in managing these patients.

RISK MANAGEMENT

AutoStrain LV/RV/LA¹ The application provides measurement tools for Strain analysis. It is powered by Auto View Recognition, Auto Contour Placement EF assessment and speckle tracking.



2D STRAIN^{1,4} Global and regional assessement of myocardial mechanics





VASCULAR **MEASUREMENTS¹** Vascular measurement package with labeled measurements for all major vessels



PEDIATRIC **MEASUREMENTS¹** Measurement package with a comprehensive range of 2D. M-Mode and Doppler pediatric measurements



STRESS ECHO^{1,3} Solution to review and analyze exercise and pharmaceutical stress echo studies



AutoLA¹ Biplane LA volume quantification by just selecting two clips



AutoIMT¹ Solution to analyze intima-media thickness of the carotid artery



4D RV-FUNCTION¹ Quantification of right ventricular volumes.



CARDIAC **MEASUREMENTS**¹ Complete echocardiographic measurement package with labeled M-Mode, Doppler and 2D measurements

ABOUT **CLINICAL SOLUTIONS APPLICATIONS**

IMAGE-COM MEASUREMENTS 3D OPTION IMAGE-COM AutoLV AutoLA AutoIMT STRESS ECHO AutoStrain LV/RV/LA **2D STRAIN FETAL 2D STRAIN 4D LV-ANALYSIS 4D RV-FUNCTION**

4D MV ASSESSMENT

4D CARDIO-VIEW

CATH OPTIONS

THERAPY MONITORING



Close monitoring based on precise data of heart failure patients is essential. Monitoring specific data allows providers to better manage and thus potentially slow down the progression to later stages of CHF.

Response to treatment and the prognosis of these conditions are directly related to the stage of the disease and the time treatment is initiated. Echocardiography is the procedure of choice in the assessment and workup for patients suffering from CHF and/or RV-Failure. The ability to closely monitor changes in volume allows clinicians to assess the trend of how a specific condition progresses and alter treatment plans as required. TOMTEC has products that allow users to monitor changes in volume, based on automated solutions for 2D and 3D echocardiography.

INCREASE EFFICIENCY



Increase diagnostic confidence

With new transcatheter therapies and minimally invasive valve repair procedures the need for advanced imaging prior to the intervention has increased.

TOMTEC's software solutions for the analysis of the mitral valve apparatus and also the aortic valve, give surgeons and interventional cardiologists a clear understanding of even the most complex anatomy. The most accurate and current information increases diagnostic confidence to successfully perform both transcatheter and surgical interventions.



4D LV-ANALYSIS¹ 3D Quantification of left ventricular volume and function function

4D RV-FUNCTION¹ Quantification of right ventricular volumes,





AutoStrain LV/RV/LA¹ AutoLV¹ "Two click" Biplane The application provides measurement tools left ventricular volume for Strain analysis. It quantification with is powered by Auto automated contour View Recognition, Auto proposals for Contour Placement FF assessment and speckle tracking.

AutoLA¹



AutoStrain LV/RV/LA¹

measurement tools

for Strain analysis. It

View Recognition, Auto

is powered by Auto

Contour Placement

and speckle tracking

4D CARDIO-VIEW^{1,2}

Generic three dimensi-

onal software solution.

Distance, area, curve

VSD size and location

The application provides



AutoLA¹ Biplane LA volume quantification by just selecting two clips



Comprehensive assess- 3D Quantification of ment and quantification left ventricular volume of mitral valve anatomy and function and angle measurements and function based on for pathologies like ASD, an automated contour detection of MV annulus and leaflets

STRUCTURAL HEART DISEASE

IMAGE-COM¹

quantification









CARDIAC MEASUREMENTS¹ Complete echocardiographic measurement package with labeled M-Mode, Doppler and 2D measurements

AutoIMT¹

Solution to analyze

of the carotid artery

intima-media thickness



Biplane LA volume quantification by just selecting two clips



2D and 3D multimodality image viewer with automated measurements and access to advanced



1 Is part of TOMTEC-ARENA.

- 2 4D CARDIO-VIEW and TOMTEC-ARENA are trademarks of TOMTEC Imaging Systems GmbH
- 3 STRESS ECHO is called ECHO-COM.
- 4 2D STRAIN is called 2D CARDIAC PERFORMANCE ANALYSIS

ABOUT **CLINICAL SOLUTIONS APPLICATIONS**

IMAGE-COM MEASUREMENTS 3D OPTION IMAGE-COM AutoLV AutoLA AutoIMT STRESS ECHO AutoStrain LV/RV/LA **2D STRAIN FETAL 2D STRAIN 4D LV-ANALYSIS 4D RV-FUNCTION 4D MV ASSESSMENT**

4D CARDIO-VIEW

CATH OPTIONS



CLINICAL APPLICATION PORTFOLIO TOMTEC ARENA TTA2 LOT40

Improve efficiency and diagnostic quality!



ABOUT CLINICAL SOLUTIONS APPLICATIONS

IMAGE-COM MEASUREMENTS 3D OPTION IMAGE-COM AutoLV AutoLA AutoIMT STRESS ECHO AutoStrain LV / RV / LA 2D STRAIN FETAL 2D STRAIN 4D LV-ANALYSIS 4D RV-FUNCTION 4D MV ASSESSMENT 4D CARDIO-VIEW

CATH OPTIONS

BASIC

Your Daily Use Routine reading and reporting, manual measurements and worksheet

CARDIO VASCULAR ULTRASOUND



IMAGE-COM¹ 2D and 3D multimodality image viewer with automated measurements and access to advanced quantification

CARDIO VASCULAR ULTRASOUND



CARDIAC MEASUREMENTS¹ measurement package with a comprehensive range labeled M-Mode, Doppler and 2D measurements



PEDIATRIC MEASUREMENTS¹ Complete echocardiographic Measurement package with of 2D, M-Mode and Doppler nediatric measurements



VASCULAR MEASUREMENTS¹ Vascular measurement package with labeled measurements for all maior vessels

STRESS ECHO^{1,3} Solution to review and analyze exercise and pharmaceutical stress echo studies

ESSENTIALS

Automated, Reproducible, Fast Speed-up your workflow, improve diagnostic quality and reproducibility

AutoStrain LV/RV/LA¹ The application provides measurement tools for Strain analusis. It is powered by Auto View Recognition, Auto Contour Placement and speckle tracking.

AutoLV "Two click" Biplane left

ventricular volume quantification with automated contour proposals for EF assessment



AutoLA¹ Biplane LA volume guantification by just selecting two clips



3D OPTION IMAGE-COM¹ Display of 3D/4D cardiac ultrasound side-by-side intima-media thickness with 2D images of the carotid arteru

Solution to analyze

MEASUREMENT **MAPPING SERVICE**

Standardization of diverse **DICOM SR formats of various** ultrasound vendors





MEASUREMENT MAPPING SERVICE¹ Standardization of different DICOM SR formats, tremendously reducing mapping efforts to information and reporting systems

PREMIUM

4

C

r

Advanced Research and Clinical Tools Tools for expert users in both, research and clinical environments

CARDIO VASCULAR ULTRASOUND



2D STRAIN^{1,4} Global and regional assessement of myocardial mechanics



4D LV-ANALYSIS¹ 3D Quantification of left ventricular volume and function



4D RV-FUNCTION¹ Quantification of right ventricular volumes and function







4D MV-ASSESSMENT¹ Comprehensive assessment and quantification of mitral valve anatomy and function based on an automated contour detection of MV annulus and leaflets

CATHLAB



IMAGE-COM¹ Advanced cathlab viewer with Biplane imaging support, display of ECG and acquisition angles



DIGITAL SUBTRACTION ANGIOGRAPHY (DSA)1 Enhanced display of peripheral vessels with automatic image summation over time and compensation of motion artifacts

CATHLAB



Cath-QCA¹ Clinical solution for quantitative coronary analysis based on an automated contour detection of coronary angiograms



Cath-QLVA¹ Quantification of left ventricular volumes and function

ABOUT **CLINICAL SOLUTIONS APPLICATIONS**

IMAGE-COM MEASUREMENTS 3D OPTION IMAGE-COM AutoLV **AutoLA** AutoIMT STRESS ECHO AutoStrain LV/RV/LA **2D STRAIN FETAL 2D STRAIN 4D LV-ANALYSIS 4D RV-FUNCTION 4D MV ASSESSMENT 4D CARDIO-VIEW**

CATH OPTIONS



IMAGE-COM

The cardiovascular solution for routine and advanced workflows at hand



ABOUT CLINICAL SOLUTIONS APPLICATIONS

IMAGE-COM MEASUREMENTS 3D OPTION IMAGE-COM AutoLV AutoLA AutoIMT STRESS ECHO AutoStrain LV/RV/LA

2D STRAIN

FETAL 2D STRAIN

4D LV-ANALYSIS

4D RV-FUNCTION

4D MV ASSESSMENT

4D CARDIO-VIEW

CATH OPTIONS

IMAGE-COM

Comprehensive image-review and measurements

IMAGE-COM¹ is an established DICOM viewer for ultrasound (US) and X-Ray Angiography (XA) examinations. Prior studies can easily be compared with current examinations and simultanously displayed with US, XA and IVUS, IVOCT/OCT and NM examinations.

IMAGE-COM provides a comprehensive feature set for review of still images and image sequences.

- A study can be reviewed with various functions like different screen tiles, zoom, flip, rotate etc.
- More efficient use of your ultrasound machines by image acquisition with modality and measuring in IMAGE-COM.
- Better traceability and transparency of the measurements through a link between each measurement and the corresponding image.



Imaging review

- Multi modality viewer
- Vendor neutral 3D/4D MPR viewer
- All TOMTEC applications can be launched from IMAGE-COM
- Side-by-side comparison (US, XA, NM, IVUS, IVOCT)
- Complete review of secondary studies and quick acess to prior finalized reports
- AVI, BMP, JPEG or DCM export
- "DSA" digital subtraction angiography
- Filter on thumbnails preview (2D, 3D, SR, PDF, Bookmark)
- Dual monitor setup capability (combined with Reporting)

Ultrasound (US) and X-Ray Angiography (XA)

- Automated US measurements value import
- XA viewer with biplane imaging support, display of ECG, acquisition angles
- Smart link between measurement and image
- Images and worksheet on different monitors
- Workflow support through measurement tools
- Worksheet management
- up to five values for each parameter
- ¬ avg, first, last, min, max selection
- XA analysis for Coronary Stenosis Quantification and Left Ventricle Analysis (auto/manual)
- Measurement export in different formats (auto/manual)
- Cardiovascular measurement packages for 2D echo, M-Mode and Doppler studies
- Control for completeness of measurement data

IMAGE-COM offers a vast number of offline measurements.

The IMAGE-COM measurement modules provide measurement capabilities for echocardiography and vascular examinations, and for X-Ray angiographic images.

IMAGE-COM has been developed to speed up your workflow in the daily routine.

VALUABLE EXTENSIONS

CARDIOVASCULAR ULTRASOUND





VASCULAR **MEASUREMENTS¹** Vascular measurement package with labeled measurements for all major vessels and 2D measurements



CATHLAB



3D OPTION IMAGE-COM¹ Display of 3D/4D cardiac ultrasound images side-by-side with 2D images



AutoIMT¹ Solution to analyze intima-media thickness of the carotid artery

nerinheral vessels with automatic image summation over time and compensation of motion artifacts



Complete echocardiographic measurement package with labeled M-Mode, Doppler







Auto Contour Placement and speckle tracking.



AutoLV¹ "Two click" Biplane left ventricular volume Strain analysis. It is powered quantification with automated contour proposals for EF assessment



AutoLA¹ Biplane LA volume quantification by just selecting two clips



DIGITAL SUBTRACTION ANGIOGRAPHY (DSA)1 Enhanced display of



Clinical solution for quantitative coronary analysis based on an automated contour detection of coronary angiograms

Cath-QLVA¹ Quantification of left ventricular volumes and

function

ABOUT **CLINICAL SOLUTIONS APPLICATIONS**

IMAGE-COM

MEASUREMENTS 3D OPTION IMAGE-COM AutoLV **AutoLA** AutoIMT STRESS ECHO AutoStrain LV/RV/LA **2D STRAIN FETAL 2D STRAIN 4D LV-ANALYSIS**

4D RV-FUNCTION

4D MV ASSESSMENT

4D CARDIO-VIEW

CATH OPTIONS

TOMTEC ACADEMY IMPRESSUM

ECHO



IMAGE-COM OPTIONS MEASUREMENTS

Томтес

IMAGE-COM OPTIONS*

MEASUREMENTS

IMAGE-COM Measurements offer a huge number of offline measurements. DICOM SRs from US carts are automatically imported into a comprehensive worksheet for easy comparison, re-measuring and export. Findings can easily be transferred to report or EMR/HIS and are available for later data mining and further analysis.

Export of multiple findings can be selected individually according to average, lowest, highest, first or last. Conformity to common guidelines: ASE, IAC, ACC and ESC.

Import from modality, re-measure, manage and export all your findings

CARDIAC MEASUREMENTS

VASCULAR MEASUREMENTS





Complete echocardiographic measurement package with labeled M-Mode, Doppler and 2D measurements.

- measurements for all major vessels
- Re-measurement of complete studies
- All measurements, including measurements imported from US systems, are summarized in the Worksheet
- Link between measurement and image: one click to get from the measurement to the image
- Up to five values for each measurement parameter with individual selection of the final value according to average, first, last, minimum or maximum
- Manual or automated measurement export in several formats: DICOM Structured Report: TID5100, TID5200, TID5220, TID3202, TID3213; EPICXML: XML

Extensive solution for cardiovascular offline measurements and multi-purpose processing of findings



				-					
	and			and a second second		and in our set			
					-				
our restaures									1.000.744
	-			a start	THE PART	HEA.	-	and the second	See out
and an	140	1.00		E.e.	1.4	6.0	1.4		323524
	1410								-
			-						Conceptore.
	-		-	-	117	100	100		in .
		-	-	1.00		100	110		1.1
entere-									a lower
ates:	141	4.05		NR.	1.	0.0	140	40.96	· Terrery Terr
at he he	Sec.	-	-		-	-	1	100	
	-		-			1 H	1.		A CONTRACTOR
**	-		-		-		*		-
	-		-			10.0	84.4	100	No.
1741-1981	- 14	1.16	-	1.10	1.	655	1.16	10.0	100
and report.	1910		-	1.00		1.1	144		and the second second

PEDIATRIC MEASUREMENTS

Measurement package with a

comprehensive range of 2D, M-Mode

and Doppler pediatric measurements

Comprehensive vascular measurement package with labeled

- Generic measurements can be assigned to a specific measurement value in a later stage (first measure, then label)
- Generic measurements such as distance, area, time, acceleration/ deceleration, VTI, velocity

ABOUT **CLINICAL SOLUTIONS APPLICATIONS**

IMAGE-COM MEASUREMENTS 3D OPTION IMAGE-COM AutoLV **AutoLA** AutoIMT STRESS ECHO AutoStrain LV/RV/LA **2D STRAIN FETAL 2D STRAIN 4D LV-ANALYSIS 4D RV-FUNCTION**

4D MV ASSESSMENT

4D CARDIO-VIEW

CATH OPTIONS



3D OPTION IMAGE-COM

TOMTEC

3D OPTION IMAGE-COM

Display 3D/4D cardiac ultrasound datasets side-by-side with 2D images. It can read and process proprietary 3D data formats of most vendors. The display and navigation of multiplanar reconstructions (MPRs) and volume rendering (VR) helps to analyze complex morphologies in 3D. With the unique navigation tool D⁺art, 3D views can be cropped easily with just two mouse-clicks.

2D/3D clinical study evaluation - it has never been so easy!

3D imaging navigation and rendering

- 3D imaging navigation and rendering
- Real-time 3D volume rendering
- ¬ D[↑]art navigation tool with `auto cropping`
- ¬ Navigation options via reference images (LOI: Lines of intersection), center point navigation and image navigation (sweep, rotate, pan, orbit)
- Auto-sweep through dataset

Visualization of vendorindependent 3D echo-data

- ¬ Color-coded 3D display for optimal depth display
- Setting of threshold and transparency
- Switch from rendering to cross plan views

3D echo data vendor-independent review



Empower the IMAGE-COM review and measurement workflow

- Empower the IMAGE-COM review and measurement workflow
- Measurements on a 2D or MPR view
- Exporting of AVI/BMP, screenshots and measurements
- Side-by-side comparison 2D and 3D

ABOUT **CLINICAL SOLUTIONS APPLICATIONS**

- IMAGE-COM **MEASUREMENTS 3D OPTION IMAGE-COM** AutoLV **AutoLA** AutoIMT STRESS ECHO AutoStrain LV/RV/LA **2D STRAIN FETAL 2D STRAIN 4D LV-ANALYSIS 4D RV-FUNCTION 4D MV ASSESSMENT**
- **4D CARDIO-VIEW**
- **CATH OPTIONS**



AutoLV

TOMTEC



AutoLV* is an automated measurement of Simpson's biplane method.

By selecting apical 2- and 4-chamber views, AutoLV automatically proposes initial contouring of the left ventricle in end-diastole and end-systole.

Left ventricular volumes and Ejection Fraction results are visualized for immediate assessment of left ventricular function

Analysis of the left ventricle with three clicks - integrated in your clinical workflow



Standard LV volume measurements with minimal user interaction

- Automated contouring in end-diastole and end-systole
- Verification of automated results by manual contour editing possible
- ¬ Volumes and EF are displayed immediatelu



Precise results within seconds

- Guideline conform measurements
- Automated biplane Simpson's method
 - High reproducibility due to automation without user interaction
- Approved algorithm by multicenter study

Biplane left ventricular volume quantification with only three clicks



Simpson	
EDV	99.92 ml
ESV	38.10 ml
SV	61,81 ml
EF	61.9%

Remain in your regular clinical workflow

- Stay in your regular review application
- ¬ No additional software has to be opened
- Directly access your results in the worksheet and report

ABOUT **CLINICAL SOLUTIONS APPLICATIONS**

IMAGE-COM **MEASUREMENTS 3D OPTION IMAGE-COM AutoLV** AutoLA AutoIMT STRESS ECHO AutoStrain LV/RV/LA **2D STRAIN** FETAL 2D STRAIN **4D LV-ANALYSIS 4D RV-FUNCTION 4D MV ASSESSMENT 4D CARDIO-VIEW**

CATH OPTIONS



AutoLA



AutoLA

Left Atrial volume is one of several parameters to identify diastolic dysfunction and a marker for LA pressure elevation.

AutoLA* is a fast and intuitive automation of Simpson's biplane method.

By selecting apical 4- and 2-chamber views, AutoLA automatically finds end-systole and proposes tracings of the left atrial blood tissue interface.



LA volume measurement with three clicks integrated in your clinical workflow



Standard LA volume measurements with minimal user interaction

- Automated contouring of maximal volume in end-systole
- Verification of automated results by manual contour editing possible
- Volume is displayed immediately

Biplane LA volume quantification with only three clicks

* AutoLA is part of TOMTEC-ARENA. TOMTEC-ARENA is a trademark of TOMTEC Imaging Systems GmbH.



Ensures diagnostic quality results

- Guideline-conform measurement
- Automated biplane disk summation
- High reproducibility due to
- automation without user interaction



Remain in your regular clinical workflow

- Stay in your regular review application
- No additional software has to be opened
- Directly access your results in the worksheet and report

ABOUT **CLINICAL SOLUTIONS APPLICATIONS**

IMAGE-COM **MEASUREMENTS 3D OPTION IMAGE-COM AutoLV AutoLA** AutoIMT STRESS ECHO AutoStrain LV/RV/LA **2D STRAIN FETAL 2D STRAIN 4D LV-ANALYSIS 4D RV-FUNCTION 4D MV ASSESSMENT**

4D CARDIO-VIEW

CATH OPTIONS

TOMTEC ACADEMY IMPRESSUM

ECHO



AutoIMT



TOMTEC

AutoIMT

AutoIMT* is an automated TOMTEC solution for the quantification of the intima media thickness. By defining a region of interest on a vascular long axis view AutoIMT automatically calculates the maximal and mean intima media thickness as well as the standard deviation and the quality index.

Automated IMT quantification - integrated in your clinical workflow

- Automated contour detection of lumen-intima and media-adventitia vessel walls
- Calculation of max and mean value of IMT
- Calculation of standard deviation and IMT quality index
- IMT quantification with high precision and reproducibility
- Automated contouring
- ¬ Immediate calculation of max and mean IMT as well as the SD and QI. - Contour overlay for visual result
- verification possible - ACCF/AHA guideline conform

Complete and automated intima media thickness quantification



- Easy measurement in familiar software environment
 - No additional software has to be opened

Integrated in IMAGE-COM Workflow

Results in the worksheet and report

ABOUT **CLINICAL SOLUTIONS APPLICATIONS**

IMAGE-COM **MEASUREMENTS 3D OPTION IMAGE-COM** AutoLV AutoLA AutoIMT STRESS ECHO AutoStrain LV/RV/LA 2D STRAIN FETAL 2D STRAIN **4D LV-ANALYSIS 4D RV-FUNCTION 4D MV ASSESSMENT**

4D CARDIO-VIEW

CATH OPTIONS



STRESS **ECHO**

TOMTEC

STRESS ECHO

STRESS ECHO* is optimized to provide a streamlined workflow for routine stress echo analysis.

The software immediately arranges the automatically synchronized loops for rapid viewing.

Additional shuffle functions are available for the comparison of arbitrary loops. An easy Wall Motion Scoring with quick assignment tools makes your analysis more comfortable.

Your vendor independent offline solution for pharmacological and exercise stress echo analysis



Sorting and scoring of stress echo studies

- Reduced mouse clicks by functions like auto-layout, sorting and synchronisation
- Quick browse through all selected clips in the review by using the dog-ear



Analysis and reporting of stress echo loops

- Easy colour-coded Wall Motion Scoring with quick assignment tools ¬ 16- and 17-segment scoring layouts - Export of Wall Motion Scoring values
 - to reporting

The intuitive workflow solution to review and analyze your 2D stress echo studies

* STRESS ECHO is called ECHO-COM and is part of TOMTEC ARENA. TOMTEC-ARENA is a trademark of TOMTEC Imaging Systems GmbH.







Customization

- User setting options for an optimal customization to your personal preferences
- Different scoring options like 4, 5 or 7 points model
- Different workflow options with or without pre-selection of clips

ABOUT **CLINICAL SOLUTIONS APPLICATIONS**

IMAGE-COM MEASUREMENTS 3D OPTION IMAGE-COM AutoLV AutoLA AutoIMT STRESS ECHO AutoStrain LV/RV/LA **2D STRAIN FETAL 2D STRAIN 4D LV-ANALYSIS 4D RV-FUNCTION**

4D MV ASSESSMENT

4D CARDIO-VIEW

CATH OPTIONS



AutoStrain LV/RV/LA

TOMTEC

AutoStrain LV/RV/LA

The AutoStrain* application provides measurement tools for Strain analysis of the left ventricle (LV), the right ventricle (RV) or the left atrium (LA) to assess cardiac function.

Just by selecting the loops and starting the application, the user gets immediate results based on automated proposals for view labels and endocardial border definitions based on Speckle tracking. This allows the user to focus directly on the interpretation of measurments.

Cardiac function assessment of LV and beyond -Strain Analysis dedicated to LA and RV



Cardiac motion analysis based on 2D clips

- Auto View Recognition: AutoSTRAIN recognizes A4C, A2C and A3C view and assigns labels
- Auto Contour Placement: Automated contour definition and speckle tracking for user review
- Cardiac cycle can be defined manually in case of missing or bad ECG signal
- Support of TTE and TEE data

Simple, fast workflow for robust and reproducible GLS measurements

* AutoStrain LV/RV/LA is part of TOMTEC-ARENA. TOMTEC-ARENA is a trademark of TOMTEC Imaging Systems GmbH



Flexible and standardized LV function assessment

- Robust, reproducible one-button GLS measurements for each view and average
- Waveform analysis for each view or all segments together
- Possibility to exclude segments with bad visibility





Deformation imaging of RV and LA

- Direct measurement of LA reservoir. conduit and contraction strain
- LA reference time point at end-diastole or right before atrial contraction
- Global RV and free wall strain measurement
- Supported strain analysis tool to perform LV, RV or LA function assessment

ABOUT **CLINICAL SOLUTIONS APPLICATIONS**

IMAGE-COM MEASUREMENTS 3D OPTION IMAGE-COM AutoLV **AutoLA** AutoIMT STRESS ECHO AutoStrain LV/RV/LA **2D STRAIN FETAL 2D STRAIN 4D LV-ANALYSIS 4D RV-FUNCTION 4D MV ASSESSMENT**

4D CARDIO-VIEW

CATH OPTIONS



2D STRAIN

TOMTEC

2D STRAIN

2D STRAIN* is a vendor independent, offline solution for the quantification of cardiac deformation. Detailed analysis of myocardial velocity, displacement, strain and strain rate is performed based on 2D speckle tracking in the long or short axis views of the left ventricle, in the right ventricle or in the left atrium. Basic parameter assessment as well as advanced quantifications, together with comprehensive result export options, make 2D STRAIN suitable for research and routine use.



Versatile analysis of left ventricle, right ventricle or left atrium





Versatile analysis options

- Analysis of LV, RV and Left Atrium
- Bull's eye analysis based on long axis or short axis views
- Same analysis independent from ultrasound system
- One application for a variety of indications
- Regional myocardial function analysis
- Assess diastolic function - Early detection of heart function
- Monitor subclinical changes or therapy over time

impairment in cardiomyopathies

Vendor independent speckle tracking analysis to quantify cardiac wall motion in 2D ultrasound data

* 2D STRAIN is called 2D CARDIAC PERFORMANCE ANALYSIS and is part of TOMTEC-ARENA. TOMTEC-ARENA is a trademark of TOMTEC Imaging Systems GmbH.







Advanced cardiac motion assessment

- Semi-automatic initial contour proposal
- Extensive export options for further post-processing
- Regional analysis of velocity, displacement, strain and strain rate
- Advanced motion visualization

ABOUT **CLINICAL SOLUTIONS APPLICATIONS**

IMAGE-COM MEASUREMENTS 3D OPTION IMAGE-COM AutoLV AutoLA AutoIMT STRESS ECHO AutoStrain LV/RV/LA **2D STRAIN FETAL 2D STRAIN 4D LV-ANALYSIS 4D RV-FUNCTION 4D MV ASSESSMENT**

4D CARDIO-VIEW

CATH OPTIONS



FETAL **2D STRAIN**

TOMTEC

FETAL 2D STRAIN

FETAL 2D STRAIN* is a vendor independent, offline solution for the quantification of cardiac deformation of the fetal heart. Detailed analysis of myocardial velocity, displacement, strain and strain rate is performed based on 2D speckle tracking in the long axis views of the left ventricle and right ventricle. Basic parameter assessment as well as advanced quantifications, together with comprehensive result export options, make FETAL 2D STRAIN suitable for research and routine use.

Vendor independent speckle tracking analysis to quantify cardiac wall motion in 2D ultrasound data







Versatile analysis options

- Analysis of LV, RV on the same 4CH clip
- Possibility to set different ES and ED for LV and RV for severe dyssynchrony cases
- Same analysis independent from ultrasound system

of indications

- Regional myocardial function analysis
- Assess diastolic function
- Early detection of heart function impairment in cardiomyopathies and congenital pathologies
- Monitor subclinical changes or therapy over time

Versatile analysis of the fetal heart





Advanced cardiac motion assessment

- Semi-automatic initial contour proposa
- ¬ Extensive export options for further post-processing
- Regional analysis of velocity displacement, strain and strain rate
- Advanced motion visualization

ABOUT **CLINICAL SOLUTIONS APPLICATIONS**

IMAGE-COM **MEASUREMENTS 3D OPTION IMAGE-COM** AutoLV AutoLA AutoIMT STRESS ECHO AutoStrain LV/RV/LA **2D STRAIN FETAL 2D STRAIN 4D LV-ANALYSIS**

4D RV-FUNCTION

4D MV ASSESSMENT

4D CARDIO-VIEW

CATH OPTIONS



40 LV-ANALYSIS

TOMTEC

4D LV-ANALYSIS

4D LV-ANALYSIS* is a vendor independent offline solution for LV assessment based on 3D speckle tracking. It provides an automated workflow for quantitative and reproducible analysis of left ventricular function and global strain values.

Regional strain analysis is mapped onto the LV model and shown in a bull's eye plot for clear visualization.

All results can be stored and exported for routine, advanced and research applications.

Quick and reproducible analysis of LV function that includes volumes, EF and strain





Fast and easy analysis of global and regional deformation

- Automated adjustment of view planes
- Easy editing of Beutel* surface in end-diastole and end-systole
- ¬ Volumes and GLS are displayed in tracking review

Ready for clinical use

systems of multiple vendors (TEE and TTE) - Fast, easy to use and accurate quantification (no geometrical

assumptions) - Reproducible results through

automated workflow

3D Quantification of left ventricular global and regional function





- Supports 3D data from ultrasound



Global and regional clinical parameters

- Accurate and reproducible volumes, EF and LV-mass
- Global/regional longitudinal and circumferential strain
- Twist and torsion
- Displacement

ABOUT **CLINICAL SOLUTIONS APPLICATIONS**

- **IMAGE-COM MEASUREMENTS 3D OPTION IMAGE-COM** AutoLV **AutoLA** AutoIMT STRESS ECHO AutoStrain LV/RV/LA **2D STRAIN FETAL 2D STRAIN 4D LV-ANALYSIS**
- **4D RV-FUNCTION**
- **4D MV ASSESSMENT**
- **4D CARDIO-VIEW**

CATH OPTIONS



4D RV-FUNCTION

TOMTEC

4D RV-FUNCTION

The right ventricular quantification tool for clinical routine work, pulmonary hypertension and right sided heart failure.

4D RV-FUNCTION* allows complete evaluation of the right ventricle and combines 3D and 2D values including EDV, EDVi, ESV, ESVi, EF and SV, RVLS, TA PSE and FAC.

4D RV-FUNCTION helps to overcome the RV complexity by calculating standard values based on a 3D surface model. A detailed and easy to follow workflow allows users to adapt the proposed contour quick and easy.

Overcome the challenges of right ventricle geometry.





Automated 2D measurements with every 3D analysis

- 3D measurement results: RV EDV, EDVi, ESV, ESVi, EF, SV
- Guideline conform 2D measurement results: TAPSE, Distances, FAC

The two minutes tool for your right ventricle analysis

- Automated contour proposal no laborious contour drawing
- Quick and reproducible results

Comprehensive right ventricular analysis, the right way





Visualize the complex RV shape in 3D

- Calculating right ventricular data based on a complex 3D surface model
- Validated model of the right ventricle
- Display your model in combination with 2D planes for a better control

ABOUT **CLINICAL SOLUTIONS APPLICATIONS**

IMAGE-COM MEASUREMENTS **3D OPTION IMAGE-COM** AutoLV **AutoLA** AutoIMT STRESS ECHO AutoStrain LV/RV/LA **2D STRAIN FETAL 2D STRAIN 4D LV-ANALYSIS 4D RV-FUNCTION**

4D MV ASSESSMENT

4D CARDIO-VIEW

CATH OPTIONS



P1

P2

4D MV-ASSESSMENT

TOMTEC



4D MV-ASSESSMENT* helps to analyze the complex anatomy of the mitral valve in 3D as well as its dynamic mechanics in systole. The mitral valve anatomy and topology are visualized with a comprehensive static and dynamic model.

Geometric measurements such as annular dimensions, leaflet morphology and coaptation descriptions can be used from the initial discovery of MV disease/pathology, support in device planning up to the monitoring of pre- and postoperative cases. Additionally, advanced navigation options and flexible manual measurements in 2D and 3D allow for individual planning of new devices.

Automated alignment and initialization proposals facilitate the workflow. Advanced edit options allow for definition and quantification of Open Coaptation Regions.

Visualization and quantification of the mitral valve in 4D



Supports decision making and monitoring

- Pre-/postoperative monitoring
- Generic manual measurements allow for individual assessment in special cases
- Highlight findings with landmarks



Facilitate communication within the heart team

Increase your diagnostic confidence for mitral valve therapy

* 4D MV-ASSESSMENT is part of TOMTEC-ARENA. TOMTEC-ARENA is a trademark of TOMTEC Imaging Systems GmbH



- Schematic visualization of MV model together with 3D volume rendering - Visualization from surgical or ventricular view direction
- Evaluation of flail/prolapse topology

Mitral valve analysis based on 4D echo data

- Streamlined workflow with semiautomated annulus detection
- Manual review and the possibility to edit the automated model proposal
- Comprehensive automatic measurements for annulus, leaflets and coaptation

ABOUT **CLINICAL SOLUTIONS APPLICATIONS**

IMAGE-COM MEASUREMENTS 3D OPTION IMAGE-COM AutoLV **AutoLA** AutoIMT STRESS ECHO AutoStrain LV/RV/LA **2D STRAIN FETAL 2D STRAIN**

4D LV-ANALYSIS

4D RV-FUNCTION

4D MV ASSESSMENT

4D CARDIO-VIEW

CATH OPTIONS

TOMTEC ACADEMY IMPRESSUM

ECHC



4D CARDIO-VIEW

TOMTEC

4D CARDIO-VIEW

4D CARDIO-VIEW* is a vendor independent offline solution to review and analyze 3D echo data.

It offers easy and fast navigation to get the perfect 3D view with just two clicks by using the unique D⁺art tool. Features like the multi-slice (multiple 2D slices), basic measurements and workflow based volume measurements make 4D CARDIO-VIEW an all-purpose solution for any cardiac structure.

All measurements and views can be stored as bookmarks for easy retrieval at any time.





No laborious navigation - Get your perfect view of complex structures in a few seconds

- Unique D[^]art navigation displays any 3D view with only two clicks
- ¬ D[↑]art multiple slice offers multiple 2D slices of a region of interest
- Smart region navigation reduces mouse miles substantially
- Landmark navigation shows your points of interest in one view
- Restore your views and measurements from bookmarks anytime

Your vendor independent all-purpose solution to review and analyze 3D echo data

3D data - Enhanced image rendering

- Adapted presets for different vendors
- Store your own render preset adjusted to your preference







- One generic tool for any cardiac structure
- Improved workflow for easier LV and generic volume measurements including LV myocardial mass
 - Distance, area, curve and angle measurements for TAVR planning and pathologies like ASD, VSD size and location
 - Measurements in volume rendered and 2D views

ABOUT **CLINICAL SOLUTIONS APPLICATIONS**

IMAGE-COM MEASUREMENTS 3D OPTION IMAGE-COM AutoLV **AutoLA AutoIMT** STRESS ECHO AutoStrain LV/RV/LA **2D STRAIN FETAL 2D STRAIN 4D LV-ANALYSIS 4D RV-FUNCTION**

4D MV ASSESSMENT

4D CARDIO-VIEW

CATH OPTIONS



IMAGE-COM OPTIONS CATH

TOMTEC



IMAGE-COM OPTIONS*

CATHLAB

Coronary angiography is routinely used for diagnosis and treatment of patients who have or are at risk of coronary artery disease. In clinical practice visual estimation has limitations such as high variability associated with operator experience and only general estimation of the vessel dimensions. Quantitative Coronary Analysis being highly accurate and reproducible can refine the visual estimate and provide important measurements of the coronary anatomy.

Integrated cath solution for QCA and QLVA

DSA





DSA is a clinical analysis package to visualize digital subtraction angiograms. It provides enhanced display of XA vessel images using automated image masking and DSA summation over time (vascular trace).

- DSA with default or user selected reference frame
- Compensation of motion artifacts using automated pixel shift

Clinical solution for quantitative coronary analysis based on an auto-

mated contour detection of coronary angiograms

- ¬ Two-clicks one border proposal ¬ Total control – verification of final
- results by user ¬ Stenosis quantification on calibrated
- and uncalibrated images Calculation of stenosis diameter
- and area, obstruction and reference diameters and obstruction length
- Graphical display: Diameter curves, Stenosis Overlays

Cath tools embedded into your daily cath routine



CATH-OLVA



Cath-QLVA is a clinical analysis package to quantify left ventricular volumes and function.

- Left ventricular volumes like EDV, ESV
- Ejection Fraction (EF), Stroke Volume (SV) and Cardiac Output (CO)
- Various Single Plane and Biplane calculation methods like Area-Length and Simpson method of disk for all results
- Review, editing and verification of final results

ABOUT **CLINICAL SOLUTIONS APPLICATIONS**

IMAGE-COM

MEASUREMENTS

3D OPTION IMAGE-COM

AutoLV

AutoLA

AutoIMT

STRESS ECHO

AutoStrain LV/RV/LA

2D STRAIN

FETAL 2D STRAIN

4D LV-ANALYSIS

4D RV-FUNCTION

4D MV ASSESSMENT

4D CARDIO-VIEW

CATH OPTIONS



TOMTEC ACADEMY



THE TRAINING AND EDUCATION **PROGRAM FOR OUR VALUED** CUSTOMERS

Detailed product trainings are offered exclusively for our up-to-date maintenance contract customers.

Use these educational programs as a chance to increase your level of expertise or simply stay up-to-date with the newest TOMTEC software.

INFO SESSION

- Complimentary for all customers
- Gain insight into our product
- portfolio

- Customer ID is required

Keeping you up-to-date on software advancements and clinical trends

- Live chat - Open discussion round
 - Demonstration of different
 - case studies

TOMTEC





TRAINING SESSION

- Exclusively for our up-to-date maintenance contract customers Refresh your knowledge and skills

Register today for unlimited complimentary sessions at www.tomtec.de/ academy



ABOUT **CLINICAL SOLUTIONS APPLICATIONS**

IMAGE-COM MEASUREMENTS 3D OPTION IMAGE-COM AutoLV AutoLA **AutoIMT** STRESS ECHO AutoStrain LV/RV/LA **2D STRAIN FETAL 2D STRAIN 4D LV-ANALYSIS 4D RV-FUNCTION 4D MV ASSESSMENT 4D CARDIO-VIEW**

CATH OPTIONS

Excellence in Digital Healthcare

The TOMTEC portfolio in this eflyer is part of TOMTEC-ARENA TTA2.40. TOMTEC-ARENA is a trademark of TOMTEC Imaging Systems GmbH.

TOMTEC IMAGING SYSTEMS GM Freisinger Strasse 9 85716 Unterschleissheim Deutschland/Germany • +49 89 32175 500 • info@tomtec.de

w www.tomtec.de

TOMTEC CORPORATION 525 W Monroe Street, Suite 2360, Chicago IL 60661, USA • + 1866 735 3740 • info@tomtecusa.com w www.tomtecusa.com

M.30.0001-06 © 2020 TOMTEC Imaging Systems GmbH. All rights reserved.



ABOUT CLINICAL SOLUTIONS APPLICATIONS

IMAGE-COM MEASUREMENTS 3D OPTION IMAGE-COM AutoLV AutoLA AutoIMT STRESS ECHO AutoStrain LV / RV / LA 2D STRAIN FETAL 2D STRAIN 4D LV-ANALYSIS 4D RV-FUNCTION 4D MV ASSESSMENT

ECHO

CATH OPTIONS

4D CARDIO-VIEW

TOMTEC ACADEMY