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
**Mesures acquises**

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Application	Abdominal
Abdomen General	Abdomen Arterial
Abdomen Renal	Abdomen Venous
ElastPQ Stiffness	<b>EQI Stiffness</b>
Left	Right

Measurements	
Liver	
Liver EQI 1	6.34 kPa
Liver EQI 2	6.57 kPa
Liver EQI 3	6.21 kPa
Liver EQI 4	6.27 kPa
Liver EQI 5	6.36 kPa
Liver EQI 6	5.41 kPa
Liver EQI 7	
Liver EQI 8	
Liver EQI 9	

Calculations	
Liver EQI Med	6.31 kPa
Liver EQI IQR/Med	2 %
Liver EQI IQR	0.15 kPa
Liver EQI Std	0.37 kPa

## Onglet Compte-rendu

Patient Demographics				
Temp ID: 20180918154741	Study Date: 09/18/2018			
Patient ID: 4141520180918	Age:	Accession #: _____	Gender:	MR ID: _____
Institution: Test Hospital		HC: _____	WR: _____	BSA: _____
Referring Physician: _____		Physician of Record: _____		Performed By: _____
Comments: _____				

Abdominal Measurements and Calculations	
EQI Liver Stiffness Calculations	
Liver EQI Avg	6.19 kPa
Liver EQI Med	6.31 kPa
Liver EQI IQR/Med	8 %
Liver EQI IQR	0.50 kPa
Liver EQI Std	0.40 kPa
Liver EQI Avg Vel	1.43 m/s
Liver EQI Med Vel	1.44 m/s
Liver EQI IQR/Med Vel	4 %
Liver EQI IQR Vel	0.06 m/s
Liver EQI Std Vel	0.04 m/s

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### Update to the Society of Radiologists in Ultrasound Liver Elastography Consensus Statement

Richard G. Barr, MD, PhD • Stephanie R. Wilson, MD • Deborah Rubens, MD • Guadalupe Garcia-Tsao, MD • Giustina Ferraroli, MD

From the Department of Radiology, Northwestern Ohio Medical University, Rootstown, Ohio (R.G.B.); Department of Radiology, University of Calgary, Calgary, Canada (S.R.W.); Department of Imaging Science, Oncology and Biomedical Engineering, University of Rochester Medical Center, Rochester, NY (D.R.); Section of Digestive Diseases, Department of Medicine, Yale University, New Haven, Conn (G.G.T.); and Ultrasound Unit, Department of Clinical Sciences and Infectious Diseases, Fondazione IRCCS Policlinico San Matteo, University of Pavia, Pavia, Italy (G.F.). Received October 31, 2019; revision requested December 11; revision received April 2, 2020; accepted April 23. Address correspondence to R.G.B., Southwoods Imaging, 7623 Market St, Youngstown, OH 44512 (e-mail: rpbarr52@gmail.com).

Conflicts of interest are listed at the end of this article.

Radiology 2020; 00:1-12 • <https://doi.org/10.1148/radiol.201912437> • Content codes [ ] [ ]

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© 2019 Koninklijke Philips N.V. Tous droits réservés. Philips se réserve le droit d'apporter des modifications aux spécifications et/ou d'abandonner tout produit à tout moment sans préavis ni obligation et ne sera pas responsable des conséquences résultant de l'utilisation de cette publication. Les marques déposées sont la propriété de Koninklijke Philips N.V. ou de leurs propriétaires respectifs.