PHILIPS

General care

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Bibliografia

- Rowan KM, Kerr JH, Major E, McPherson K, Short A, Vessey MP. Intensive Care Society's APACHE II study in Britain and Ireland - I. Variations in case mix of adult admissions to general intensive care units and impact on outcome. *British Medical Journal* 1993; 307: 972-7.
- Murphy A, Cronin J, Whelan R, Drummond FJ, Savage E, Hegarty J. Economics of Early Warning Scores for identifying clinical deterioration - a systematic review. *Ir J Med Sci*. 2018; 184: 193-205.
- American Hospital Association. Trend Watch Chartbook 2018, disponibile su: https://www.aha.org/system/files/2018-05/2018-AHA-Chartbook_0.pdf.
- HealthGrades, Inc. *Third Annual Patient Safety in American Hospitals Study*. Aprile 2006. <http://www.proctmasspatients.org/docs/PatientSafetyInAmericanHospitalsStudy2006.pdf>
- Young MP, Gooder VJ, McBride K, James B, Fisher ES. Inpatient transfers to the intensive care unit: delays are associated with increased mortality and morbidity. *J Gen Intern Med*. 2003; 18: 77-83.
- Rutherford P, Lee B, Greiner A. *Transforming Care at the Bedside*. IHI Innovation Series white paper. Boston: Institute for Healthcare Improvement; 2004. (Disponibile su: www.IHI.org).
- Fuhrmann L, Lippert A, Perner A, Østergaard D. Incidence, staff awareness and mortality of patients at risk on general wards. *Resuscitation*. 2008; 77(3):325-30.
- Monteiro C, Ferreira Machado Avelar A, da Luz Gonçalves Pedreira M. Interruptions of nurses' activities and patient safety: an integrative literature review. *Rev Lat Am Enfermagem*. 2015; 23(1): 169-179.
- Slawomirski L, Auroaen A, Klazinga N. The economics of patient safety. Report dell'OCSE. 2017. Disponibile su: <https://www.oecd.org/els/health-systems/The-economics-of-patient-safety-March-2017.pdf>.
- Healy MA, Mullard AJ, Campbell DA Jr, Dimick JB. Hospital and Payer Costs Associated With Surgical Complications. *JAMA Surg*. 2016; 151(9):823-830.
- Moerer O, Burchardi H. The Cost of Sepsis. *ICU Management and Practice*. 2004; 4(2). Disponibile su: <https://healthmanagement.org/c/icu/issuearticle/the-cost-of-sepsis>, vedere i dati di base
- Kazure H, Roman S, Sosa J. A. Population-Level Analysis of 5620 Recipients of Multiple In-Hospital Cardiopulmonary Resuscitation Attempts. *Journal of Hospital Medicine*. 2014; 9(1): 29-34.
- IHI. *Early Warning Systems: Scorecards That Save Lives*. Disponibile su: <http://www.ihl.org/resources/Pages/ImprovementStories/EarlyWarningSystemsScorecardsThatSaveLives.aspx>
- Smith MEB, Chiovaro J, O'Neil M, Kansagara D, Quinones A, Freeman M, Motu'apuaka M, Slatore CG. *Early Warning Scoring Systems: A Systematic Review*. VA-ESP Project #05-225; 2013
- Royal College of Physicians. National Early Warning Score (NEWS) 2. Disponibile su: <https://www.rcplondon.ac.uk/projects/outputs/national-early-warning-score-news-2>
- NCEC. Guideline: National Early Warning Score. 2013. <http://hdl.handle.net/10147/317078>.
- IHI. 5 Million Lives Campaign. Disponibile su: <http://www.ihl.org/Engage/Initiatives/Completed/5MillionLivesCampaign/Pages/default.aspx>
- Fernando SM, Fox-Robichaud AE, Rochweg B, Cardinal P, Seely AJE, Perry JJ, McIsaac DI, Tran A, Skitch S, Tam B, Hickey M, Reardon PM, Tanuseputro P, Kyeremanteng K. Prognostic accuracy of the Hamilton Early Warning Score (HEWS) and the National Early Warning Score 2 (NEWS2) among hospitalized patients assessed by a rapid response team. *Crit Care*. 2019; 23(1):60.
- Ravikirti. Early Warning Scoring System for Early Recognition of and Timely Intervention in Deteriorating Patients in the Hospital. *JAPI* 2016; 64.
- Royal College of Physicians. National Early Warning Score (NEWS) 2: Standardising the assessment of acute-illness severity in the NHS. *Updated report of a working party*. London: RCP, 2017.
- Green N, Lander H, Snyder A, Hudson P, Churpek M, Edelson D. Comparison of the Between the Flags calling criteria to the MEWS, NEWS and the electronic Cardiac Arrest Risk Triage (eCART) score for the identification of deteriorating ward patients. *Resuscitation*. 2018; 123:86-91.
- Jones DA, DeVita MA, Bellomo R. Rapid-response teams. *N Engl J Med*. 2011; 365:139-46 [PMID: 21751906].
- Churpek MM, Yuen TC, Edelson DP. Risk stratification of hospitalized patients on the wards. *Chest*. 2013; 143(6):1758-1765.
- Subbe CP, Kruger M, Rutherford P, Gemmel L. MEWS: Validation of a modified Early Warning Score in medical admissions. *QJM*. 2001; 94(10):521-6.
- Smith GB, Prytherch DR, Meredith P, Schmidt PE, Featherstone PI. The ability of the National Early Warning Score (NEWS) to discriminate patients at risk of early cardiac arrest, unanticipated intensive care unit admission, and death. *Resuscitation*. 2013; 84(4):465-70.
- Churpek MM, Yuen TC, Winslow C, Robicsek AA, Meltzer DO, Gibbons RD, Edelson DP. Multicenter development and validation of a risk stratification tool for ward patients. *Am J Respir Crit Care Med*. 2014; 190(6):649-55.
- Kipnis P, Turk BJ, Wulf DA, LaGuardia JC, Liu V, Churpek MM, Romero-Brufau S, Escobar GJ. Development and validation of an electronic medical record-based alert score for detection of inpatient deterioration outside the ICU. *J Biomed Inform*. 2016; 64:10-19.
- Bein B, Seewald S, Gräsner JT. How to avoid catastrophic events on the ward. Best practice & research. *Clinical Anaesthesiology*. 2016; 30(2), 237-245.
- Fox A, Elliott N. Early warning scores: a sign of deterioration in patients and systems. *Nurs Manag (Harrow)*. 2015; 22(1):26-31.
- PeraHealth. Leveraging the Power of EHR Data to Improve Patient Safety, a Comparison of the Rothman Index with MEWS and NEWS. 18 settembre 2017. <https://www.perahealth.com/resources/white-papers/intelligent-early-warning-a-comparison-of-the-rothman-index-with-mews-and-news/>.
- Ghosh E, Eshelman L, Yang L, Carlson E, Lord B. Early deterioration Indication: data-driven approach to detecting deterioration in general ward. *Resuscitation*. 2018; 122:99-105.
- van Galen LS, Struik PW, Driesen BE, Merten H, Ludikhuizen J, van der Spoel JI, Kramer MH, Nanayakkara PW. Delayed Recognition of Deterioration of Patients in General Wards Is Mostly Caused by Human Related Monitoring Failures: A Root Cause Analysis of Unplanned ICU Admissions. *PLoS One*. 2016; 11(8): e0161393.
- Fletcher GS, Aaronson BA, White AA, Julka R. Effect of a Real-Time Electronic Dashboard on a Rapid Response System. *J Med Syst*. 2017; 42: 5.
- Simmes FM, Schoonhoven L, Mintjes J, Fikkers BG, van der Hoeven JG. Incidence of cardiac arrests and unexpected deaths in surgical patients before and after implementation of a rapid response system. *Annals of intensive care*. 2012; 2(1):20
- Ludikhuizen J, Smorenburg SM, de Rooij SE, de Jonge E. Identification of deteriorating patients on general wards; measurement of vital parameters and potential effectiveness of the Modified Early Warning Score. *Journal of Critical Care*. 2012; 27(4):424.e7-13.
- Pedersen NE, Rasmussen LS, Petersen JA, Gerds TA, Østergaard D, Lippert A. A critical assessment of early warning score records in 168,000 patients. *Journal of Clinical Monitoring and Computing*. 2018; 32(1): 109-16.
- Subbe CP, Duller B. Continuous Monitoring of Respiratory Rate on General Wards What might the implications be for Clinical Practice? *Acute Med*. 2018; 17(1):5-9.
- Jensen JK, Skår R, Tveit B. The impact of Early Warning Score and Rapid Response Systems on nurses' competence: An integrative literature review and synthesis. *J Clin Nurs*. 2018; (7-8):e1256-e1274.
- Burns KA, Reber T, Theodore K, Welch B, Roy D, Siedlecki SL. Enhanced early warning system impact on nursing practice: A phenomenological study. *J Adv Nurs*. 2018; 74(5):1150-1156.
- Saab MM, McCarthy B, Andrews T, Savage E, Drummond FJ, Walshe N, Forde M, Breen D, Henn P, Drennan J, Hegarty J. The effect of adult Early Warning Systems education on nurses' knowledge, confidence and clinical performance: A systematic review. *J Adv Nurs*. 2017; 73(11):2506-2521.
- Downey CL, Tahir W, Randell R, Brown JM, Jayne DG. Strengths and limitations of Early Warning Scores: a systematic review and narrative synthesis. *International Journal of Nursing Studies*. 2017; 76:106-19.
- Prgomet M, Cardona-Morrell M, Nicholson M, Lake 1, Long J, Westbrook J, Braithwaite J, Hillman K. Vital signs monitoring on general wards: clinical staff perceptions of current practices and the planned introduction of continuous monitoring technology. *Int J Qual Health Care*. 2016; 28(4): 515-21.
- Douw G, Huisman-de Waal G, van Zanten ARH, van der Hoeven JG, Schoonhoven L. Capturing early signs of deterioration: the dutch-early-nurse-worry-indicator-score and its value in the Rapid Response System. *J Clin Nurs*. 2017; 26(17-18):2605-2613.
- Subbe CP, Duller B, Bellomo R. Effect of an automated notification system for deteriorating ward patients on clinical outcomes. *Critical Care*. 2017; 21(1):52.
- Michard F, Bellomo R, Taenzer A. The rise of ward monitoring: opportunities and challenges for critical care specialists. *Intensive Care Med*. 2019; 45(5):671-673.
- Bellomo R, Ackerman M, Bailey M, Beale R, Clancy G, Danesh V, Hvarfner A, Jimenez E, Konrad D, Lecardo M, Pattee KS, Ritchie J, Sherman K, Tangkau P. Vital Signs to Identify, Target, and Assess Level of Care Study (VITAL Care Study) Investigators. A controlled trial of automated advisory vital signs monitoring in general hospital wards. *Crit Care Med*. 2012; 40(8).
- Fletcher GS, Aaronson BA, White AA, Julka R. Effect of a Real-Time Electronic Dashboard on a Rapid Response System. *J Med Syst*. 2017; 42: 5.
- Danesh V, Neff D, Jones TL, Aroian K, Unruh L, Andrews D, Guerrier L, Venus SJ, Jimenez E. Can proactive rapid response team rounding improve surveillance and reduce unplanned escalations in care? A controlled before and after study. *Int J Nurs Stud*. 2019; 91:128-133.
- Bein B, Seewald S, Gräsner JT. How to avoid catastrophic events on the ward. Best practice & research. *Clinical Anaesthesiology*. 2016; 30(2): 237-245.
- Capan M, Ivy JS, Wilson JR, Huddleston JM. A stochastic model of acute-care decisions based on patient and provider heterogeneity. *Health Care Management Science*. 2017; 20: 187-206.
- Michard F, Gan T, Bellomo R. Protecting ward patients the case for continuous monitoring. *ICU Management & Practice*. 2019.
- Brown H, Terrence J, Vasquez P, Bates DW, Zimlichman E. Continuous monitoring in an inpatient medical-surgical unit: a controlled clinical trial. *Am J Med*. 2014; 127:226-232.
- Winters BD. Rapid Response Systems. A brief review of the evidence. In: DeVita MA, Hillman K, Bellomo R (Eds). *Textbook of Rapid Response Systems*. 2nd edition. Springer 2017.
- McGaughy J, Alderdice F, Fowler R, Kapila A, Mayhew A, Moutray M. Outreach and Early Warning Systems (EWS) for the prevention of intensive care admission and death of critically ill adult patients on general hospital wards. *Cochrane Database Syst Rev*. 2007; (3):CD005529.
- Winters BD, Pham JC, Hunt EA, Guallar E, Berenholtz S, Pronovost PJ. Rapid response systems: a systematic review. *Crit Care Med*. 2007; 35(5):1238-43.
- Ranji SR, Auerbach AD, Hurd CJ, O'Rourke K, Shojanian KG. Effects of rapid response systems on clinical outcomes: systematic review and meta-analysis. *J Hosp Med*. 2007; 2(6):422-32.
- Jones DA, DeVita MA, Bellomo R. Rapid-response teams. *N Engl J Med*. 2011; 365:139-46.
- Chan PS, Jain R, Nallmothu BK, Berg RA, Sasson C. Rapid response teams: a systematic review and metaanalysis. *Arch Intern Med*. 2010; 170:18-26.
- Alan M, Hobbelenk EL, van Tienhoven AJ, van de Ven PM, Jansma EP, Nanayakkara PW. The impact of the use of the Early Warning Score (EWS) on patient outcomes: A systematic review. *Resuscitation*. 2014; 85(5):587-94.
- Smith ME, Chiovaro JC, O'Neil M, Kansagara D, Quinones AR, Freeman M, Motu'apuaka ML, Slatore CG. Early warning system scores for clinical deterioration in hospitalized patients: a systematic review. *Ann Am Thorac Soc*. 2014; 11(9):1454-65.
- Maharaj R, Raffaele I, Wendon J. Rapid response systems: a systematic review and meta-analysis. *Critical Care*. 2015; 19(1): 254.
- Winters BD, Weaver SJ, Pfoh ER, Yang T, Pham JC, Dy SM. Rapid-response systems as a patient safety strategy: a systematic review. *Ann Intern Med*. 2013; 158(5 Pt 2):417-25.
- Mathukia C, Fan WQ, Vadyak K, Biege C, Krishnamurthy M. Modified Early Warning System improves patient safety and clinical outcomes in an academic community hospital. *Journal of Community Hospital Internal Medicine Perspectives*. 2015; 5:2, 26716.
- Ludikhuizen J, Brunsvelde-Reinders AH, Dijkgraaf MG, Smorenburg SM, de Rooij SE, Adams R, de Maaijer PF, Fikkers BG, Tangkau P, de Jonge E. Cost and Outcomes of Medical Emergency Teams Study Group. Outcomes Associated With the Nationwide Introduction of Rapid Response Systems in The Netherlands. *Crit Care Med*. 2015; 43(12):2544-51.
- Heller AR, Mees ST, Lauterwald B, Reeps C, Koch T, Weitz J. Detection of Deteriorating Patients on Surgical Wards Outside the ICU by an Automated MEWS-Based Early Warning System With Paging Functionality. *Ann Surg*. 2018; Pubblicazione elettronica prima della stampa. DOI: 10.1097/SLA.0000000000002830.
- Moore D, Poynton M. Business case for investing in a quality improvement programme to reduce harm caused by clinical deterioration. Sapere Report 2015; disponibile su: <https://www.hqsc.govt.nz/assets/Deteriorating-Patient/PR/Deteriorating-patient-investment-case-Sapere-report-Oct-2015.pdf>.
- HIQA (Autorità per la qualità e l'informazione sanitaria) (2015) Health technology assessment of the use of information technology for early warning and clinical handover systems. Health Information and Quality Authority, Irlanda
- Bailey JG, Davis PB, Levy AR, Molinari M, Johnson PM. The impact of adverse events on health care costs for older adults undergoing non-elective abdominal surgery. *Canadian Journal of Surgery*. 2016; 59(3), 172-179.
- Simmes F, Schoonhoven L, Mintjes J, Adang E, van der Hoeven JG. Financial consequences of the implementation of a rapid response system on a surgical ward. *J Eval Clin Pract*. 2014; 20(4):342-7.
- Bonafide CP, Localio AR, Song L, Roberts KE, Nadkarni VM, Priestley M, Paine CW, Zander M, Lutts M, Brady PW, Keren R. Cost-Benefit Analysis of a Medical Emergency Team in a Children's Hospital. *Pediatrics*. 2014; 134(2):235-241.
- Ward W. The Business Case for Implementing Rapid Response Teams. Presentation for VHA INC. 2006; disponibile su: <http://www.ihl.org/resources/Pages/Presentations/BusinessCaseforImplementingRRTsPresentation.aspx>