A healthcare professional, likely a nurse or technician, is shown in a clinical setting. She is wearing blue scrubs and has a name tag. She is holding a drip chamber, a component used in medical procedures, and looking upwards with a focused expression. The background is slightly blurred, showing what appears to be a patient bed and other medical equipment.

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
Fig.1 Beat to beat variable delay algorithm

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Reference

1. SCOT-HEART Investigators, Newby, D. E., et al. (2018). Coronary CT Angiography and 5-Year Risk of Myocardial Infarction. *New England Journal of Medicine*, 379(10), 924–933.
2. Knuuti, J., et al. (2020). 2019 ESC Guidelines for the diagnosis and management of chronic coronary syndromes. *Eur Heart J*, 41(3), 407–477.
3. Vembar, M. (2003). A dynamic approach to identifying desired physiological phases for cardiac imaging using multislice spiral CT. *Med Phys*, 30(7), 1683–1693.
4. Garcia, M. J., CATSCAN Study Investigators., et al. (2006). Accuracy of 16-row multidetector computed tomography for the assessment of coronary artery stenosis. *Jama*, 296(4), 403–411.
5. Lehmkühl, L., et al. (2011). Diagnostic performance of prospectively ECG triggered versus retrospectively ECG gated 64-slice computed tomography coronary angiography in a heterogeneous patient population. *Eur J Radiol*, 80(2), 342–348.
6. Chao, S.-P., et al. (2010). The diagnostic accuracy of 256-row computed tomographic angiography compared with invasive coronary angiography in patients with suspected coronary artery disease. *Eur Heart J*, 31(15), 1916–1923.
7. Kim, J., et al. (2019). Diagnostic performance and image quality of iterative model-based reconstruction of coronary CT angiography using 100 kVp for heavily calcified coronary vessels. *PLoS One*, 14(9), e0222315–.
8. Ecabert, O., et al. (2011). Segmentation of the heart and great vessels in CT images using a model-based adaptation framework. *Med Image Anal*, 15(6), 863–876.
9. Lee, M. S., et al. (2010). Reproducibility in the assessment of noncalcified coronary plaque with 256-slice multi-detector CT and automated plaque analysis software. *Int J Cardiovasc Imaging*, 26(Suppl 2), 237–244.

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©2020 Philips Japan, Ltd.