


# Optimize stents with IVUS

### Optimize short- and long-term patient outcomes using a defined pre- and post-stent criteria with Philips IVUS.<sup>1-2</sup>

The treatment of coronary lesions with DES is complicated by the challenges of suboptimal post stent results, such as stent under expansion, incomplete stent apposition, edge dissection and geographic miss, which may contribute to the risk of stent failure.<sup>3</sup>


Recent publications suggest angiography alone is not enough and that the use of IVUS to guide stent planning and optimization delivers enduring reduction in MACE and TVF.



### ULTIMATE trial<sup>1</sup>

Compared to angiography-guided PCI alone, in all comer patients, ULTIMATE shows that IVUS-guided PCI significantly reduces clinically driven Target Lesion Revascularization (TLR). This benefit was sustained over three years.

**Only 1.6% TVF\* at 12 months and 4.2% TVF at 3 years when optimal IVUS-guided PCI criteria was met.**



\*TVF categorized as cardiac death, target vessel MI and clinically driven TVR


Study outcomes were strongest when specified criteria for an optimal PCI procedure were met.

### IVUS-XPL 5-yr follow up<sup>2</sup>

Patients in the IVUS-guided stent group who did not meet IVUS criteria for stent optimization had a significantly higher incidence of MACE at 5 years compared with those who met IVUS criteria for stent optimization.

**MACE 7.4% vs 4.0% (p=0.048)**

**Only 1.5% MACE at 12 months and 4.0% at 5 years when optimal IVUS-guided PCI criteria was met.**




Number at risk					
315	285	276	270	261	253
363	334	326	320	313	292

# Plan, land and expand

Plan the procedure with pre-stent IVUS to size stents and identify the optimal landing zone. Perform post-stent IVUS to confirm good stent expansion, apposition and no edge dissection.


## Plan

Pre-stent IVUS to guide procedural planning




### Stent/vessel diameter

- Measure distal reference and determine stent diameter
- Measure proximal reference to guide post stent dilation



### Landing zone/stent length

Identify most normal segments adjacent to target lesion, with plaque burden <50% at 5mm. Measure distance between distal and proximal reference point, round up to the nearest stent length.




### Lesion

Assess lesion characteristics to guide plaque modification strategy.

## Land and expand

Post-stent IVUS to guide and confirm stent optimization



Criteria used for IMPact on Revascularization Outcomes of intraVascular Ultrasound Guided Treatment of Complex Lesions and Economic Impact (IMPROVE) trial. Shlofmitz et al.<sup>4</sup>

### Landing zone

Plaque burden <50% at 5 mm proximal and distal to stent edge.

### Expansion and apposition

MSA  $\geq$ 90% of the distal reference lumen area and full stent apposition throughout.

### Stent edges

No edge dissection involving media with length > 3mm and arc  $\geq$  60°.

1. Zhang J et al. Three year outcomes of the The ULTIMATE trial. Intravascular Ultrasound-Guided Versus Angiography-Guided Implantation of Drug-Eluting Stent in All-Comers. Presented at TCT Connect, 15th October 21.
2. Hong et al. 5-Year Follow-Up of the IVUS-XPL Randomized Trial. JACC; 13 Jan, 2020:62 -71.
3. N. Pal, J. Din. Et al Contemporary Management of Stent Failure: Part One/Interventional Cardiology Review 2019;14(1):10-6. DOI:<https://doi.org/10.15420/icr.2018.39.1>
4. Shlofmitz et al. IMPROVE trial: Study design and rationale. AHJ; (2020) Oct, Vol 228, doi.org/10.1016/j.ahj.2020.08.002

