VisiSheath™

Dilator Sheath

Instructions for Use

Spectranetics®

Instructions for Use

Instructions for Use - Sections by Language

Pg	Language	Instructions For Use
3	English	Instructions For Use
7	Bulgarian	Инструкции за употреба
12	Croatian	Upute za uporabu
16	Czech	Návod k použití
20	Danish	Brugsanvisninger
24	Dutch	Gebruiksaanwijzing
28	Estonian	Kasutusjuhend
32	Finnish	Käyttöohjeet
36	French	Mode d'emploi
40	German	Gebrauchsanweisung
44	Greek	Οδηγίες χρήσης
48	Hungarian	Használati utasítások
52	Icelandic	Leiðbeiningar um notkun
56	Italian	Istruzioni per l'uso
60	Latvian	Lietošanas instrukcija
64	Lithuanian	Naudojimo instrukcija
68	Norwegian	Bruksinstruksjoner
72	Polish	Instrukcja użytkowania
76	Portuguese	Instruções de utilização
80	Romanian	Instrucțiuni de utilizare
84	Russian	Инструкции по применению
88	Serbian	Upustva za upotrebu
92	Slovak	Návod na použitie
96	Slovenian	Navodila za uporabo
100	Spanish	Instrucciones de uso
104	Swedish	Bruksanvisning
108	Turkish	Kullanım Talimatları

Instructions for Use

English / English

Table of Contents

1.	DESCRIPTION	3
2.	INDICATIONS FOR USE	3
3.	CONTRAINdications	3
4.	WARNINGS	3
5.	PRECAUTIONS	3
6.	ADVERSE EVENTS	4
7.	HOW SUPPLIED	4
8.	COMPATIBILITY	4
9.	DIRECTIONS FOR USE	5
10.	MANUFACTURER'S LIMITED WARRANTY	6
11.	NON-STANDARD SYMBOLS	6

CAUTION: Federal law restricts this device to sale by or on the order of a physician.

1. DESCRIPTION

The VisiSheath Dilator Sheath is a single lumen sheath used independently or as a support for an inner sheath to facilitate tissue dilation. One end is terminated with a 45° angle cut, while the other end is blunt, see Figure 1. Both ends contain a radiopaque marker band to enable fluoroscopic identification of tip location and orientation. An additional exterior mark aligned with the tip of the 45° angle cut permits visual identification of sheath orientation. There are multiple diameter and length options available, see Table 1.

Figure 1: VisiSheath Dilator Sheath Tip Configuration

2. INDICATIONS FOR USE

The VisiSheath Dilator Sheath is intended for use in patients requiring the percutaneous dilation of tissue surrounding cardiac leads. The device is also intended for use in the introduction and support of intravascular catheters.

3. CONTRAINDICATIONS

None known.

4. WARNINGS

- Dilator sheaths should be used only at institutions with thoracic surgical capabilities by physicians knowledgeable in the techniques and devices for lead removal.
- When using dilator sheaths, do not insert sheaths over more than one lead at a time. Severe vessel damage, including venous wall laceration requiring surgical repair, may occur.
- Do not maintain a stationary position with the VisiSheath tip at the Superior Vena Cava (SVC)-right atrial (RA) junction as it may result in damage to this delicate area during subsequent lead extraction and reinsertion procedures (e.g., manipulating the dilator sheath or implanting a new lead).
- Weigh the relative risks and benefits of intravascular lead dilation procedures especially in cases when:
 - the object to be dilated away from adherent tissue is of a dangerous shape or configuration,
 - the likelihood of lead disintegration may result in increased risk of fragment embolization,
 - vegetations are attached directly to the lead body.

5. PRECAUTIONS

- For single use only.** The VisiSheath Dilator Sheath must not be resterilized and/or reused.
- DO NOT resterilize or reuse this device, as these actions can compromise device performance or increase the risk of cross-contamination due to inappropriate reprocessing.
- Reuse of this single use device could lead to serious patient injury or death and voids manufacturer warranties.
- Do not alter the sheath from its original state prior to use.
- When the VisiSheath Dilator Sheath is in the body, it should be manipulated only under fluoroscopic observation with radiographic equipment that provides high quality images.
- Prior to the procedure, evaluate the physical dimensions of the lead, catheter, or inner sheath in relation to the specifications of the dilator sheath to determine possible incompatibility.
- If selectively removing leads with the intent to leave one or more chronic leads implanted intact, the non-targeted leads must be subsequently tested to ensure that they were not damaged or dislodged during the procedure.
- When advancing dilator sheaths, use proper sheath technique. Maintain adequate tension and coaxial alignment on the lead to minimize the risk of vessel wall damage.
- If excessive scar tissue or calcification prevents safe advancement of dilator sheaths, consider an alternate approach.
- Excessive force with dilator sheaths used intravascularly may result in damage to the vascular system requiring emergency surgical repair.

Instructions for Use

English / English

- If the lead breaks, evaluate fragment for retrieval.
- If hypotension develops, rapidly evaluate; treat as appropriate.
- Due to rapidly evolving lead technology, this device may not be suitable for dilation of tissue around all types of leads. If there are questions or concerns regarding compatibility of this device with particular leads, contact the lead manufacturer.
- Do not pull on the lead because it may stretch, distort, or break, making subsequent removal more difficult. Damage to a lead may prevent passage of a lead locking device through the lumen and/or make dilation of scar tissue more difficult.
- If a lead locking device has not been used, be aware that damage to the lead caused by pulling on it may prevent subsequent passage of a lead locking device through the lumen and/or make dilation of scar tissue more difficult.
- If removing a chronic pacing lead, be aware that if it is freed spontaneously during the extraction procedure, the lead tip may become trapped in the upper vasculature. Dilator sheaths, advanced at least to the innominate vein, are often necessary to extract the lead tip through the scar tissue at the site of venous entry, and to avoid a venotomy.
- If the dilator sheath fails to progress after initial success, or if advancing the sheath was difficult, remove the sheath to inspect the tip. If the tip is distorted or frayed, exchange the damaged sheath for a new sheath before continuing treatment.
- When advancing a sheath around a bend, keep the point of the sheath's beveled tip oriented toward the inside of the bend.

6. ADVERSE EVENTS

Potential Adverse Events

Potential adverse events related to the procedure of intravascular removal of leads include (listed generally in order of increasing potential effect):

- dislodging or damaging non-targeted lead
- chest wall hematoma
- thrombosis
- arrhythmias
- bacteremia
- hypotension
- pneumothorax
- migrating fragment from lead
- migration of vegetation from lead
- pulmonary embolism
- laceration or tearing of vascular structures or the myocardium
- hemopericardium
- cardiac tamponade
- hemothorax
- stroke
- death

7. HOW SUPPLIED

7.1 Sterilization

For single use only. Do not re-sterilize and/or reuse. The VisiSheath Dilator Sheath is supplied sterile and non-pyrogenic. Sterility is guaranteed only if the packaged is unopened and undamaged.

7.2 Device / Packaging Care

Store devices in a dry cool place (below 60° C / 140° F) until use.

7.3 Inspection Prior to Use

Prior to use, visually inspect the sterile package to ensure that seals have not been broken. The sheath should be examined carefully for defects or damage. Do not use if it has apparent defects or is damaged.

8. COMPATIBILITY

Use with other devices

The VisiSheath Dilator Sheath may be used as a support sheath for compatibly sized inner sheaths, including the Spectranetics Laser Sheath (SLS™ II) used for removal of cardiac leads. It may also be used in conjunction with a Spectranetics Lead Locking Device (LLD™).

Be sure to closely follow the "Instructions for Use" for each device used.

Table 1 provides physical dimensions and specifications for use in determining VisiSheath compatibility with other devices.

Table 1: VisiSheath Dilator Sheath Size Compatibility Specifications

Model (Ref)	Length (cm)	Sheath Size (Diameter)	Minimum Inner Diameter (mm/in/F)	Maximum Outer Diameter (mm/in/F)	SLS™ II Compatibility
501-012	43	S	4.2 / 0.168 / 12.8	5.5 / 0.215 / 16.4	12 F
501-112	33				
501-212	23				
501-014	43	M	5.0 / 0.198 / 15.0	6.5 / 0.253 / 19.3	14 F
501-114	33				
501-214	23				
501-016	43	L	5.9 / 0.236 / 17.9	7.5 / 0.293 / 22.4	16 F
501-116	33				
501-216	23				

Instructions for Use

English / English

9. DIRECTIONS FOR USE

9.1 Procedure Setup

- Obtain a thorough patient history, including blood type. Appropriate blood products should be readily available.
- Determine the manufacturer, model number, physical dimensions, and implant date of the target lead. Perform radiographic evaluation of position, type and condition of target lead.
- Use a procedure room that has high quality fluoroscopy, pacing equipment, defibrillator, and thoracotomy and pericardiocentesis trays.
- Prep and drape the patient's chest for possible thoracotomy.
- Arrange for immediate surgical back-up.
- Establish back-up pacing as needed.
- Using sterile technique, open the sterile package.
- Make available, and open as needed, any other adjunct devices such as inner sheaths (laser or non-laser), lead locking devices, or related accessories.

9.2 Clinical Technique

1. Surgically expose the proximal end of the target lead(s) and remove the lead from its connections (if connected).
2. Remove all suture and tie-down materials.
3. Cut off all proximal fittings, if present, using clippers or other cutters. It is important to cut the lead very close to the connector (but past any crimp joints) leaving as long a portion of the target lead(s) to work with as possible. Avoid closing off the interior lumen (or coil) of the lead when cutting it.

PRECAUTION: Do not pull on the lead because it may stretch, distort, or break, making subsequent removal more difficult. Damage to a lead may prevent passage of a lead locking device through the lumen and/or make dilation of scar tissue more difficult.

4. Advance a lead locking device down the lead lumen as distal as possible and deploy the locking mechanism.
5. Unless the lead exterior is damaged, degraded or too thin, tie a suture at the proximal end of the lead to use as a traction element (the suture may be attached to the proximal portion of a lead locking device).

6. For an active fixation lead, attempt to unscrew the lead fixation mechanism via counterclockwise rotation by rotating the lead (and lead locking device if used) counterclockwise.
7. Gently apply traction to the lead to determine if it is still engaged in tissue. If the lead is sufficiently free of binding tissue, gently pull on the lead locking device (if used) and lead to remove it.

PRECAUTION: If a lead locking device has not been used, be aware that damage to the lead caused by pulling on it may prevent subsequent passage of a lead locking device through the lumen and/or make dilation of scar tissue more difficult.

PRECAUTION: If removing a chronic pacing lead, be aware that if it is freed spontaneously during the extraction procedure, the lead tip may become trapped in the upper vasculature. Dilator sheaths, advanced at least to the innominate vein, are often necessary to extract the lead tip through the scar tissue at the site of venous entry, and to avoid a venotomy.

8. If the lead is not removed from the vessel with gentle pulling, dilator sheaths (or other retrieval devices) may help separate the lead from any tissue encapsulation by advancing the dilator sheath over the target lead and any accompanying lead locking device.

9. A size-compatible inner sheath, including a laser sheath, may be used with a VisiSheath Dilator Sheath acting as a support sheath. The VisiSheath Dilator Sheath is loaded over the inner sheath by inserting the distal end of inner sheath through the proximal end of the dilator sheath until it appears beyond the proximal end of the dilator sheath. The inner sheath and VisiSheath Dilator Sheath may be advanced in alternating fashion to progressively dilate tissue along the length of the lead.

10. Always maintain adequate tension and coaxial alignment on the lead to support the maneuvering of the dilator sheaths and any accompanying inner sheaths to properly guide them within the patient anatomy.

With too little tension, the sheaths may damage the vein. Too much tension may cause a myocardial avulsion. The size of the dilator sheath should be large enough that the sheath can be advanced over the lead without causing the lead to buckle or its exterior jacket to distort, but the sheaths should not be overly loose.

Rotating the sheaths during advancement may facilitate progress through resistant scar tissue.

Always use fluoroscopic monitoring when manipulating the dilator sheath within the vascular system. Never use sheaths on more than one lead at a time.

PRECAUTION: If the dilator sheath fails to progress after initial success, or if advancing the sheath was difficult, remove the sheath to inspect the tip. If the tip is distorted or frayed, exchange the damaged sheath for a new sheath before continuing treatment.

PRECAUTION: When advancing a sheath around a bend, keep the point of the sheath's beveled tip oriented toward the inside of the bend.

Instructions for Use

English / English

11. For cardiac leads, if the lead has not been freed by the time the dilator sheath nears the myocardial lead fixation point, position the sheath end against the myocardium, with any inner sheath retracted several centimeters into the dilator sheath. Countertraction may be applied to dilate the remaining tissue at the lead tip. This is accomplished by holding the dilator sheath about one centimeter from the myocardium while steadily and gently applying tension to the lead. Rotation of the sheath may help dilate remaining tissue at the lead tip.
12. When all binding tissue surrounding the target lead has been successfully dilated, the lead may slide freely out of the body with applied traction.
13. Dispose of the used products according to local biological handling and disposal procedures.

10. MANUFACTURER'S LIMITED WARRANTY

Manufacturer warrants that the VisiSheath is free from defects in material and workmanship when used by the stated "Use By" date and when package is unopened and undamaged immediately before use. Manufacturer's liability under this warranty is limited to replacement or refund of the purchase price of any defective VisiSheath. Manufacturer will not be liable for any incidental, special, or consequential damages resulting from use of the VisiSheath. Damage to the VisiSheath caused by misuse, alteration, improper storage or handling, or any other failure to follow these Instructions for Use will void this limited warranty. **THIS LIMITED WARRANTY IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING THE IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.** No person or entity, including any authorized representative or reseller of Manufacturer, has the authority to extend or expand this limited warranty and any purported attempt to do so will not be enforceable against Manufacturer.

11. NON-STANDARD SYMBOLS

Inner Diameter		Outer Diameter	
Quantity	QTY	Working Length	
Non-Pyrogenic		Importer	
CAUTION: Federal law restricts this device to sale by or on the order of a physician.			Rx ONLY

References

1. Verma A, Wilkoff BL. Intravascular pacemaker and defibrillator lead extraction: A state-of-the-art review. Heart Rhythm, Vol. 1, No. 6, December 2004; 739-745.
2. Smith HJ, et al. Five-years experience with intravascular lead extraction. PACE 1994; 17:2016-2020.
3. Byrd CL, et al. Lead extraction: indications and techniques. Cardiol Clin 1992; 10:735-748.
4. Byrd CL, et al. Intravascular techniques for extraction of permanent pacemaker leads. J Thorac Cardiovasc Surg 1991;101:989-997.
5. Byrd CL, et al. Intravascular lead extraction using Locking Stylets and sheaths. PACE 1990; 13:1871-1875.

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1. Verma A, Wilkoff BL. Intravascular pacemaker and defibrillator lead extraction: A state-of-the-art review. Heart Rhythm, Vol. 1, No. 6, December 2004; 739-745.
2. Smith HJ, et al. Five-years experience with intravascular lead extraction. PACE 1994; 17:2016-2020.
3. Byrd CL, et al. Lead extraction: indications and techniques. Cardiol Clin 1992; 10:735-748.

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4. Byrd CL, et al. Intravascular techniques for extraction of permanent pacemaker leads. J Thorac Cardiovasc Surg 1991;101:989-997.
5. Byrd CL, et al. Intravascular lead extraction using Locking Stylets and sheaths. PACE 1990; 13:1871-1875.

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Rx ONLY

Reference

1. Verma A, Wilkoff BL. Intravascular pacemaker and defibrillator lead extraction: A state-of-the-art review. Heart Rhythm, Vol. 1, No. 6, December 2004; 739-745.
2. Smith HJ, et al. Five-years experience with intravascular lead extraction. PACE 1994; 17:2016-2020.
3. Byrd CL, et al. Lead Extraction: indications and techniques. Cardiol Clin 1992; 10:735-748.
4. Byrd CL, et al. Intravascular techniques for extraction of permanent pacemaker leads. J Thorac Cardiovasc Surg 1991;101:989-997.
5. Byrd CL, et al. Intravascular lead extraction using Locking Stylets and sheaths. PACE 1990; 13:1871-1875.

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Rx ONLY

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- Verma A, Wilkoff BL. Intravascular pacemaker and defibrillator lead extraction: A state-of-the-art review. Heart Rhythm, Vol. 1, No. 6, December 2004; 739-745.
- Smith HJ, et al. Five-years experience with intravascular lead extraction. PACE 1994; 17:2016-2020.
- Byrd CL, et al. Lead extraction: indications and techniques. Cardiol Clin 1992; 10:735-748.
- Byrd CL, et al. Intravascular techniques for extraction of permanent pacemaker leads. J Thorac Cardiovasc Surg 1991;101:989-997.
- Byrd CL, et al. Intravascular lead extraction using Locking Stylets and sheaths. PACE 1990; 13:1871-1875.

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Rx ONLY

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
- Verma A, Wilkoff BL. Intravascular pacemaker and defibrillator lead extraction: A state-of-the-art review. Heart Rhythm, Vol. 1, No. 6, December 2004; 739-745.
- Smith HJ, et al. Five-years experience with intravascular lead extraction. PACE 1994; 17:2016-2020.
- Byrd CL, et al. Lead extraction: indications and techniques. Cardiol Clin 1992; 10:735-748.
- Byrd CL, et al. Intravascular techniques for extraction of permanent pacemaker leads. J Thorac Cardiovasc Surg 1991;101:989-997.
- Byrd CL, et al. Intravascular lead extraction using Locking Styles and sheaths. PACE 1990; 13:1871-1875.

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©2022 Spectranetics Corporation
Spectranetics Corporation