# **Quick-Cross<sup>TM</sup>**

## **EXTREME<sup>TM</sup>**

Support Catheter

## Instructions for Use

**Quick-Cross™  
EXTREME™**  
Support Catheter

## Instructions for Use

English / English

### Instructions for Use - Sections by Language

#### Table of Contents

Pg	Language	Instructions For Use
3	English	Instructions for Use
5	Bulgarian	български език
7	Croatian	Hrvatski
9	Czech	Česky
11	Danish	Dansk
13	Dutch	Nederlands
15	Estonian	Eesti
17	Finnish	Suomi
19	French	Français
21	German	Deutsch
23	Greek	Ελληνικά
25	Hungarian	Magyar
27	Icelandic	Íslenska
29	Italian	Italiano
31	Latvian	Latviešu valoda
33	Lithuanian	lietuvių kalba
35	Norwegian	Norsk bokmål
37	Polish	Polski
39	Portuguese	Português
41	Romanian	Română
43	Russian	Русский язык
45	Serbian	Српски језик
47	Slovak	slovenčina
49	Slovenian	Slovenčina
51	Spanish	Español
53	Swedish	Svenska
55	Turkish	Türkçe

# Instructions for Use

English / English

## Table of Contents

1. DESCRIPTION.....	3
2. INDICATIONS FOR USE.....	3
3. CONTRAINDICATIONS.....	3
4. WARNINGS .....	3
5. PRECAUTIONS.....	3
6. ADVERSE EFFECTS.....	3
7. HOW SUPPLIED.....	3
8. COMPATIBILITY.....	4
9. PROCEDURE SET-UP AND DIRECTIONS FOR USE.....	4
10. MANUFACTURER'S LIMITED WARRANTY.....	4
11. NON-STANDARD SYMBOLS.....	4

**CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician**

### 1. DESCRIPTION

The Spectranetics Quick-Cross Extreme Support Catheters are intravascular catheters. These catheters are available in a variety of lengths and tip configurations. All models have 3 radiopaque markers spaced equally along the distal shaft to aid in estimating geometry within the vascular system. The distal radiopaque marker is positioned within 3 mm of the distal catheter tip. A standard female luer is placed on the proximal end of each model. The catheter is coated with a lubricious, hydrophilic coating.

### 2. INDICATIONS FOR USE

Quick-Cross Extreme Support Catheters are intended to guide and support a guidewire during access of the coronary or peripheral vasculature, allow for wire exchanges and provide a conduit for the delivery of saline solutions or diagnostic contrast agents.

### 3. CONTRAINDICATIONS

None known.

### 4. WARNINGS

- Maximum Infusion Pressure: **300 psi**. for 0.014" & 0.018" catheters and **500 psi** for 0.035" catheters.
- The catheter is designed and intended for one time use only. **Do not re-sterilize and/or reuse.**
- The catheter should only be used by physicians qualified to perform percutaneous vascular interventions.

### 5. PRECAUTIONS

- DO NOT resterilize or reuse this device, as these actions can compromise device performance or increase the risk of cross-contamination due to inappropriate reprocessing. Reuse of this single use device could lead to serious patient injury or death and voids manufacturer warranties.
- Store in a cool, dry place. Protect from direct sunlight and high temperature (greater than 55°C or 131°F).
- Do not use if device or packaging is damaged.
- Use the catheter prior to the "Use By" date (Expiration Date) specified on the package.

### 6. ADVERSE EFFECTS

Vascular catheterization and/or vascular intervention may result in complications including but not limited to:

- Vessel dissection, perforation, rupture or total occlusion
- Infection
- Hematoma
- Unstable angina
- Embolism
- Hypo/hypertension
- Acute myocardial infarction
- Arrhythmia, including ventricular fibrillation
- Death

### 7. HOW SUPPLIED

The Spectranetics Quick-Cross Extreme Support Catheters are supplied **STERILE**. The devices are designated and designed for **SINGLE USE ONLY** and must not be resterilized and/or reused.

#### 7.1 Sterilization

- The sterility of the product is guaranteed only if the package is unopened and undamaged. Before use, visually inspect the sterile package to ensure that the seals have not been broken. Do not use the catheter if the integrity of the package has been compromised. Do not use catheter if its "Use Before Date," found on package labeling, has been passed.

#### 7.2 Inspection Prior to Use

- Before use, examine carefully for defects, all of the equipment to be used. Do not use any equipment if it is damaged.

# Instructions for Use

English / English

## 8. COMPATIBILITY

Product Features & Model Numbers Table

Model (Ref)	Working Length	Tip Shape / Type	Hydrophilic Coating Length	GW Compatibility	Minimum Guide I.D.	Minimum Sheath I.D.	Maximum Catheter O.D.	Marker Band Spacing
518-084	135 cm	Straight	100 cm	0.014" 0.36 mm	5F	4F	3.2F	15 mm
518-086	150 cm	Straight					0.042" 1.07 mm	
518-088	90 cm	Straight	60 cm	0.018" 0.46 mm	5F	4F	3.4F	15 mm
518-090	135 cm	Straight	100 cm				0.044" 1.12 mm	
518-092	150 cm	Straight	0.035" 0.89 mm	N/A	5F	4.5F	50 mm	
518-076	65 cm	Straight				45 cm		0.059"
518-078	90 cm	Straight				70 cm		1.50 mm
518-080	135 cm	Straight	115 cm					
518-082	150 cm	Straight	130 cm					

## 9. PROCEDURE SET-UP AND DIRECTIONS FOR USE

**Note:** Follow instructions for use for all equipment to be used with the Quick-Cross Extreme Support Catheter.

- Preparation: Using sterile technique, open the sterile package. Gently remove the protective hoop with the catheter from the pouch. Fill a sterile standard luer-lock syringe with sterile saline. Before removing the catheter from the hoop, connect the syringe to the catheter proximal luer fitting, flush the catheter and allow the saline to fill the hoop. Set catheter in hoop aside until ready for use.
- Insertion: Through a previously inserted, appropriately sized guiding catheter or introducer sheath, introduce the catheter over an appropriately sized guidewire (see specifications) using standard technique.
- Advancement: Use fluoroscopic guidance when advancing the catheter to the desired location within the vasculature.
- Removal: Gently withdraw the catheter using standard technique, being careful to maintain guidewire position if the guidewire is to remain in place.

**Infusion:** To perform infusion, withdraw the guidewire and reference the specifications for maximum infusion pressure.

**Note:** Do not exceed the maximum infusion pressures.

After use, dispose of all equipment in accordance with applicable requirements relating to hospital waste, and potentially bio-hazardous materials.

## 10. MANUFACTURER'S LIMITED WARRANTY

Manufacturer warrants that the Quick-Cross Extreme is free from defects in material and workmanship when used by the stated "Use By" date and when package is unopened and undamaged immediately before use. Manufacturer's liability under this warranty is limited to replacement or refund of the purchase price of any defective Quick-Cross Extreme. Manufacturer will not be liable for any incidental, special, or consequential damages resulting from use of the Quick-Cross Extreme. Damage to the Quick-Cross Extreme caused by misuse, alteration, improper storage or handling, or any other failure to follow these Instructions for Use will void this limited warranty. **THIS LIMITED WARRANTY IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING THE IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.** No person or entity, including any authorized representative or reseller of Manufacturer, has the authority to extend or expand this limited warranty and any purported attempt to do so will not be enforceable against Manufacturer.

## 11. NON-STANDARD SYMBOLS

Distal Marker Spacing	
Shape	
GW Compatibility	
Maximum OD	
Maximum PSI/kPa	
Sheath Compatibility	
Working Length	
Quantity	QTY
CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician.	Rx ONLY
Importer	

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©2020 Spectranetics Corporation
**Spectranetics Corporation**
**Spectranetics International B.V.**