Quick-Access™
Needle Holder

Instructions for Use

 **Spectranetics®**

CE
0086

Instructions for Use - Sections by Language

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
Table of Contents

1.	DESCRIPTION	3
2.	INDICATIONS FOR USE	3
3.	CONTRAINDICATION	3
4.	WARNINGS	3
5.	PRECAUTIONS	3
6.	POTENTIAL COMPLICATIONS OR ADVERSE EVENTS	3
7.	HOW SUPPLIED	3
8.	COMPATIBILITY.....	3
9.	PROCEDURE SET-UP AND DIRECTIONS FOR USE	3
10.	MANUFACTURER'S LIMITED WARRANTY	4
11.	NON-STANDARDS SYMBOLS	4

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician with appropriate training.

1. DESCRIPTION

The Quick-Access Needle Holder is used to extend the proximal portion of the access needle during standard vascular access technique. The Quick-Access Needle Holder is a sterile, single-use, single lumen rigid polymer tube having a male luer connector at the distal end to provide connection to a standard needle. A wire lock shut-off connector is located at the proximal end to lock a guidewire in position or release it. A side port is located adjacent to the connector.

2. INDICATIONS FOR USE

The Quick-Access Needle Holders are intended to facilitate the placement of guidewires into the vascular system.

3. CONTRAINDICATION

None.

4. WARNINGS

- Do not advance the guidewire or the Quick-Access Needle Holder if resistance is met.
- Use only appropriate guidewires of size 0.035" and below.

5. PRECAUTIONS

- Procedures should be conducted under fluoroscopic guidance with appropriate x-ray equipment.
- The sealed needle container should be inspected prior to opening. If the seal is broken, or the container has been damaged, sterility cannot be assured.
- Careful attention must be paid to maintain tight connection between the Quick-Access Needle Holder and needle.
- Do not re-sterilize or reuse this device, as these actions can compromise device performance or increase the risk of cross-contamination due to inappropriate reprocessing.
- Reuse of this single use device could lead to serious patient injury or death and voids manufacturer warranties.

6. POTENTIAL COMPLICATIONS OR ADVERSE EVENTS

Potential complications related to angioplasty include, but are not limited to clot formation and embolism, nerve damage, vascular perforation requiring surgical repair, damage to the vascular intima, or death.

7. HOW SUPPLIED

7.1 Sterilization

Supplied sterilized and sealed in peel-open package. Intended for single use; do not re-sterilize, re-process or re-use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether packaging has been compromised.

7.2 Inspection prior to use

- Upon removal from the package, thoroughly inspect the product to ensure that no damage has occurred during shipment.
- Do not use Quick-Access Needle Holder if its expiration date, found on the package labeling, has passed.

8. COMPATIBILITY

Guidewire: compatible with guidewire sizes up to and including 0.035".

9. PROCEDURE SET-UP AND DIRECTIONS FOR USE

- Remove the Quick-Access Needle Holder from its packaging.
- Inspect the Quick-Access Needle Holder for damage before use.
- Purge the Quick-Access Needle Holder through lumen, observing for obstructions.
- Attach a needle to the Quick-Access Needle Holder.
- Insert guidewire at least 3 cm into Quick-Access Needle Holder's proximal end.
- Fully engage the shut-off connector slider so that the guidewire is locked in position.
- Insert the needle into the target vessel and advance using standard vascular access technique.
- Confirm vascular access via blood back flow, disengage the shut-off connector slider and advance the guidewire.

**Spectranetics**[®]
Quick-Access[™]
Needle Holder

**Instructions
for Use**

English / English

- i) Maintain guidewire position and remove the needle and Quick-Access Needle Holder.
- j) Proceed with further percutaneous intervention as planned.
- k) Apply pressure to the insertion site according to standard practice or hospital protocol after removing all catheters and guidewires.

10. MANUFACTURER’S LIMITED WARRANTY

Manufacturer warrants that the Quick-Access Needle Holder is free from defects in material and workmanship when used by the stated “Use By” date and when package is unopened and undamaged immediately before use. Manufacturer’s liability under this warranty is limited to replacement or refund of the purchase price of any defective Quick-Access Needle Holder. Manufacturer will not be liable for any incidental, special, or consequential damages resulting from use of the Quick-Access Needle Holder. Damage to the Quick-Access Needle Holder caused by misuse, alteration, improper storage or handling, or any other failure to follow these Instructions for Use will void this limited warranty. **THIS LIMITED WARRANTY IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING THE IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.** No person or entity, including any authorized representative or reseller of Manufacturer, has the authority to extend or expand this limited warranty and any purported attempt to do so will not be enforceable against Manufacturer.

11. NON-STANDARDS SYMBOLS

Guidewire Compatibility		Quantity	QTY
Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.			Rx ONLY

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©2014 Spectranetics Corporation
Manufactured by Spectranetics Corporation