Quick-CrossTM
CaptureTM
Guidewire Retriever

Instructions for Use

CE
0086

Instructions for Use

Instructions for Use - Sections by Language

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Table of Contents


1.	DESCRIPTION	3
2.	INDICATIONS FOR USE	3
3.	CONTRAINDICATION	3
4.	WARNINGS	3
5.	PRECAUTIONS	3
6.	POTENTIAL COMPLICATIONS OR ADVERSE EVENTS	4
7.	HOW SUPPLIED	4
8.	COMPATIBILITY	4
9.	PROCEDURE SET-UP AND DIRECTIONS FOR USE	4
10.	MANUFACTURER'S LIMITED WARRANTY	4
11.	NONSTANDARDS SYMBOLS	4

Caution:Federal (USA) law restricts this device to sale by or on the order of a physician with appropriate training.

1. DESCRIPTION

The Quick-Cross Capture Guidewire Retriever is a coaxial lumen design catheter with a balloon and flexible funnel mounted on its distal tip. The outer lumen is for balloon inflation and is accessed through the balloon port, the side female luer at the proximal end of the Y-connector. The inner through lumen permits the use of guidewires to facilitate the positioning of the catheter. The wire lumen is accessed through the second port of the Y-connector. A balloon protector is provided to facilitate insertion of the funnel and balloon into the introducer sheath. Two radiopaque marker bands identify the distal section of the catheter.

Model	Funnel Size	Working Length
519-106	6 mm	110cm

Pressure (ATM)	0.5	1.0	1.5	2.0	2.5	3.0	3.5	4.0	4.5	5.0
Balloon Diameter (mm)	6.4	6.6	6.8	7.3	7.5	7.8	7.9	8.2	8.5	9.0

Rated Burst Pressure: 5 ATM

2. INDICATIONS FOR USE

The Quick-Cross Capture Guidewire Retriever is intended to be used in conjunction with a steerable guidewire to access discrete regions of the vasculature and for guidewire exchange. The Quick-Cross Capture Guidewire Retriever is not intended for use in the coronary, cerebral or carotid vasculature.

3. CONTRAINDICATION

None.

4. WARNINGS

- Do not advance the guidewire or the Quick-Cross Capture Guidewire Retriever if resistance is met.
- Pressure in excess of 5 ATM can cause balloon rupture.
- Use an inflation device with a pressure gauge for balloon inflation/deflation.
- Use only appropriate balloon inflation medium. Do not use air or gaseous medium to inflate the balloon.

5. PRECAUTIONS

- Conduct procedures under fluoroscopic guidance with appropriate x-ray equipment.
- Inspect sealed catheter container prior to opening. If the seal is broken, or the container has been damaged, sterility cannot be assured.
- Ensure tight catheter connections to avoid air introduction into the system.
- Completely deflate balloon before removing the catheter from the sheath.
- Do not advance the catheter system against resistance. Identify the cause of the resistance and take corrective action. If resistance is felt upon catheter removal, then remove the system and the sheath simultaneously, particularly if balloon rupture or leakage is suspected.
- Care should be used when handling. Damage may result from kinking, stretching, or forcefully wiping the catheter.
- Do not resterilize or reuse this device, as these actions can compromise device performance or increase the risk of cross-contamination due to inappropriate reprocessing.
- Reuse of this single use device could lead to serious patient injury or death and voids manufacturer warranties.

Instructions for Use

English / English

6. POTENTIAL COMPLICATIONS OR ADVERSE EVENTS

- Potential complications related to angioplasty include but are not limited to: clot formation and embolism, nerve damage, vascular perforation requiring surgical repair, damage to the vascular intima, or death.
- Potential complications related to introducing the catheter into the body include but are not limited to: infection, hematoma formation, and air emboli.
- Potential balloon or funnel separation or rupture due to misuse may necessitate the use of a snare or other medical interventional technique to retrieve detached pieces.

7. HOW SUPPLIED

7.1 Sterilization

- Supplied sterilized and sealed in peel-open package. Intended for single use; do not re-sterilize, re-process or re-use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether packaging has been compromised.

7.2 Inspection prior to use

- Upon removal from the package, thoroughly inspect the product to ensure that no damage, bends or kinks have occurred during shipment.
- Do not use catheter if expiration date, found on the package labeling, has passed.

8. COMPATIBILITY

- Sheath: Compatible with 6 French or larger Introducer Sheath.
- Guidewires: Compatible with 0.014", 0.018", or 0.035" guidewires.

9. PROCEDURE SET-UP AND DIRECTIONS FOR USE

- Inspect sealed catheter container prior to opening. If the seal is broken, or the container has been damaged, sterility cannot be assured.
- Remove the Quick-Cross Capture Guidewire Retriever from its packaging and inspect the catheter for damage before use.
- Advance the balloon protector of the catheter so that it covers the balloon and funnel at the catheter's distal tip to facilitate introduction to introducer. Purge the catheter's guidewire lumen, observing for leaks and obstructions.
- When inflating the balloon, use either 50/50 or 75/25 dilution of saline and contrast medium, respectively.
- Attach an inflation device with a pressure gauge, partially filled with contrast solution, to the balloon port of the catheter.
- Point inflation device with pressure gauge nozzle downward, aspirate until all air is removed from the balloon and bubbles no longer appear in the contrast solution.
- Prior to catheter insertion, confirm the appropriate guidewire being used.
- Enter the vessel percutaneously using a 0.014", 0.018", or 0.035" guidewire, using standard vascular access technique. Insert a 6 French compatible sheath or larger. Insert the balloon protector into the hub of the sheath, advance the Quick-Cross Capture Guidewire Retriever approximately 5cm into the sheath, and retract the balloon protector to the proximal end of the catheter.
- Advance the Quick-Cross Capture Guidewire Retriever using accepted angioplasty techniques (see references¹) until gaining confirmation fluoroscopically that the funnel, identified by the distal marker, is proximal to the occlusion.
- Pull the catheter backwards about 1cm to allow the funnel to regain its intended shape, and allow easier guidewire exchange.
- Inflate the balloon with saline and contrast medium to desired size (see chart above).
- Advance the guidewire across the lesion, into the funnel of the Quick-Cross Capture Guidewire Retriever and exit through its proximal port.
- Exchange guidewires as necessary, using the two radiopaque markers to confirm guidewire location and trajectory.
- Deflate the catheter's balloon by drawing a vacuum via an inflation device with a pressure gauge. The greater the vacuum applied and held during the withdrawal, the lower the deflated balloon profile.
- Gently withdraw the catheter, leaving the guidewire in place to complete the interventional procedure. As the balloon exits the vessel, use a smooth, gentle, steady motion. If resistance is felt upon removal, then the catheter should be removed under fluoroscopic guidance, particularly if the balloon ruptures or leakage is suspected.
- Then proceed with further percutaneous intervention such as atherectomy, balloon dilation and/or stent placement.
- Apply pressure to the insertion site according to standard practice or hospital protocol after removing all catheters and guidewires.

10. MANUFACTURER'S LIMITED WARRANTY

Manufacturer warrants that the Quick-Cross Capture Guidewire Retriever is free from defects in material and workmanship when used by the stated "Use By" date and when package is unopened and undamaged immediately before use. Manufacturer's liability under this warranty is limited to replacement or refund of the purchase price of any defective Quick-Cross Capture Guidewire Retriever. Manufacturer will not be liable for any incidental, special, or consequential damages resulting from use of the Quick-Cross Capture Guidewire Retriever. Damage to the Quick-Cross Capture Guidewire Retriever caused by misuse, alteration, improper storage or handling, or any other failure to follow these Instructions for Use will void this limited warranty. **THIS LIMITED WARRANTY IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING THE IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.** No person or entity, including any authorized representative or reseller of Manufacturer, has the authority to extend or expand this limited warranty and any purported attempt to do so will not be enforceable against Manufacturer.

11. NONSTANDARDS SYMBOLS

Guidewire Compatibility		Sheath Compatibility	
Catheter Working Length		Funnel Size	
Quantity	QTY		
Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.			Rx ONLY

¹ A) Gruntzig A.; Kumpke, D.; Technique of PTA with the Gruntzig Balloon Catheter, Amer Journ of Roent 132:547-552 (1979)

B) Zeitler, E.; Gruntzig A.; Schoop, W.; (Eds), Percutaneous Vascular Recanalization , Springer-Verlag, Heidelberg/New York (1978)

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©2018 Spectranetics Corporation
Manufactured by Spectranetics Corporation