Merknader for operatøren

Lumify diagnostisk ultralydssystem
17. november 2021

4535 621 06811 Rev A

Om Merknader for operatøren

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Baumgartner, Helmut, et al. "Echocardiographic Assessment of Valve Stenosis: EAE/ASE Recommendations for Clinical Practice." *European Journal of Echocardiography*, 10: 1-25, 2009.

Calafiore, P., Stewart, W.J. "Doppler Echocardiographic Quantitation of Volumetric Flow Rate," *Cardiology Clinics*, Vol. 8, No. 2: 191-202, May 1990.

Rudski, Lawrence, et al. "Guidelines for the Echocardiographic Assessment of the Right Heart in Adult: A Report from the American Society of Echocardiography." *Journal of the American Society of Echocardiography*, Vol. 23, No. 7: 685-713, 2010.

Zoghbi, William, et al. "Recommendations for Evaluation of Prosthetic Valves with Echocardiography and Doppler Ultrasound." *Journal of the American Society of Echocardiography*, Vol. 22. No. 9: 975-1014, 2009.

Maks. trykkgradient (forenklet Bernoulli-likning)

Silverman, N. H., Schmidt, K. G. "The Current Role of Doppler Echocardiography in the Diagnosis of Heart Disease in Children." *Cardiology Clinics*, Vol. 7, No. 2: 265–96, mai 1989.

Reynolds, T. *The Echocardiographer's Pocket Reference, Second Edition*. Arizona Heart Institute Foundation, Phoenix, AZ, 2000, p. 382.

Maks. trykkgradient (full Bernoulli-likning)

Silverman, N. H., Schmidt, K. G. "The Current Role of Doppler Echocardiography in the Diagnosis of Heart Disease in Children." *Cardiology Clinics*, Vol. 7, No. 2: 265–96, mai 1989.

Middeltrykkgradient

Reynolds, T. *The Echocardiographer's Pocket Reference, Second Edition*. Arizona Heart Institute Foundation, Phoenix, AZ, 2000, p. 382.

Trykhalveringstid

Hatle, L., Angelsen, B., Tromsal, A. "Noninvasive Assessment of Atrioventricular pressure half-time by Doppler Ultrasound" *Circulation*, Vol. 60, No. 5: 1096-104, November, 1979.

Hastighetstidsintegral (VTI)

Silverman, N. H., Schmidt, K. G. "The Current Role of Doppler Echocardiography in the Diagnosis of Heart Disease in Children." *Cardiology Clinics*, Vol. 7, No. 2: 265–96, mai 1989.

Vaskulære referanser

Dopplerhastighet (HAST)

Krebs, C. A., Giyanani, V. L., Eisenberg, R. L. *Ultrasound Atlas of Vascular Diseases*, Appleton & Lange, Stamford, CT, 1999.

Endediastolisk hastighet (EDV)

Strandness, D. E., Jr. *Duplex Scanning in Vascular Disorders*. Lippincott, Williams & Wilkins, Philadelphia, PA, 2002.

Minimum diastolisk hastighet (MDV)

Evans, D. H., McDicken, W. N. *Doppler Ultrasound Physics, Instrumentation, and Signal Processing, Second Edition*. John Wiley & Sons, Ltd., 2000.

Topp trykkgradient (PG)

Powls, R., Schwartz, R. *Practical Doppler Ultrasound for the Clinician*. Williams & Wilkins, Baltimore, Maryland, 1991.

Topp systolisk hastighet (PSV)

Krebs, C. A., Giyanani, V. L., Eisenberg, R. L. *Ultrasound Atlas of Vascular Diseases*, Appleton & Lange, Stamford, CT, 1999.

Motstandsindeks (RI)

Zwiebel, W. J., ed. *Introduction to Vascular Ultrasonography, Third Edition*. W. B. Saunders Company, Philadelphia, PA 1992.

Systolisk/diastolisk-forhold (S/D)

Zwiebel, W. J., ed. *Introduction to Vascular Ultrasonography, Third Edition*. W. B. Saunders Company, Philadelphia, PA 1992.

Hastighetstidsintegral (VTI)

Reynolds, T. *The Echocardiographer's Pocket Reference, Second Edition*. Arizona Heart Institute Foundation, Phoenix, AZ, 2000, p. 383.

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© 2021 Koninklijke Philips N.V. Med enerett. Publisert i USA.
Philips Ultrasound, Inc.