

PHILIPS

Sleep Position Therapy

NightBalance Lunoa

Specialist portal
user guide



Specialist portal: Setup

1 Account is set up by NightBalance service desk

Philips NightBalance 

Home Patients Physicians Distributors Devices Settings Data Export Inbox About Logout

New Physician


Initials	<input type="text" value="J.R."/>
First name	<input type="text" value="John"/>
Last name	<input type="text" value="Smith"/>
Email	<input type="text" value="jrsmith@hospital.com"/>
Confirm email	<input type="text" value="jrsmith@hospital.com"/>
Phone prefix	<input type="text"/>
Phone number	<input type="text"/>
Mobile prefix	<input type="text"/>
Mobile number	<input type="text"/>
Gender	<input checked="" type="radio"/> Male <input type="radio"/> Female <input type="radio"/> Other
Language	<input type="text" value="English, United States"/>
Country	<input type="text" value="Netherlands"/>
Registration number	<input type="text"/>

Hospitals

« 1 to 1 of 1 »

Specialist portal: **Setup**


2 Service desk confirms registration and sends confirmation email



3 Specialist receives email to verify account




4 Specialist opens email and clicks the link to activate account and set password



Specialist portal: **Setup**


5 After clicking the activation link, the window to change password pops up



The screenshot shows a web form for changing a password. At the top center is the Philips logo, a shield with a cross and waves, with the word "PHILIPS" above it. Below the logo is the text "Philips NightBalance". A horizontal teal line separates the header from the form content. The form contains the following elements:

- A label: "Please enter your new password"
- A text input field labeled "New Password" with a small icon of a password strength indicator on the right.
- A text input field labeled "Confirm Password" with a small icon of a password strength indicator on the right.
- A teal button labeled "Submit New Password".

6 Specialist changes password and receives confirmation




The screenshot shows a confirmation message box with the following elements:

- A title bar: "Information" with a close button (X) on the right.
- A message: "The password has been updated."
- An "OK" button.

Specialist portal: **Setup**

7 Specialist is automatically directed to the login page



The screenshot shows the Philips NightBalance login page. At the top center is the Philips logo and the text "Philips NightBalance". Below this, there are two main sections: "New user?" and "Already have an account?".


New user?

- A dropdown menu showing "English, United States".
- A teal "Create account" button.

Already have an account?

- An input field for "email address (lower case)".
- An input field for "password".
- A teal "Login" button.
- A link for "Forgot your password?".

8 Specialist logs in for the first time




This screenshot shows the "Already have an account?" section of the login page. The email address field contains "jsmith@hospital.com" and the password field is masked with dots. A teal "Login" button and a "Forgot your password?" link are visible at the bottom.

Already have an account?

- Input field: "jsmith@hospital.com"
- Input field: "....."
- Teal "Login" button
- Link: "Forgot your password?"

Specialist portal: Setup

9 Edit and change personal details/ compliance criteria if needed



The 'DETAILS' form allows editing personal information. It includes fields for Email (molly.lewis@philips.com), Initials, First Name (John), Last Name (Doe), Phone (+45), Mobile (+45), Country (Denmark), and Language (English, United States). There are buttons for 'Change Password', 'Change Email', and 'Upload Profile Image'. At the bottom, there are 'Save', 'Cancel', 'Deactivate Account', and 'Delete Account' buttons.


Personal details

The 'COMPLIANCE SETTINGS' form allows adjusting compliance criteria. It includes a field for 'Number of consecutive days over which the device must be used for at least 4 hours to be considered compliant' (30) and a field for 'Percentage of nights during the above period which the device must be used for at least 4 hours to be considered compliant' (70). There are 'Save' and 'Cancel' buttons.

Compliance criteria

Specialist portal: Contents

1 Dashboard



Specialist portal: Contents

2 Patient registration Inbox

Physician will receive a notification that a new patient has selected them as their specialist

From this moment, their data are accessible in the physician portal

Print current page here

The screenshot displays the LUNOA specialist portal interface. On the left is a teal sidebar with navigation links: Home, Patients, My Profile, Inbox, About, and Logout. The NightBalance logo is at the bottom of the sidebar. The main content area is titled 'Inbox' and features a search bar and a refresh button. A notification card is visible with the following text:

LUNOA: NEW PATIENT REGISTRATION 06-12-2018

A patient has selected you as their specialist.

If this patient gave approval to share their data to your Lunoa account, the patient data is now accessible from the Patients page in the Lunoa Portal.

Patient: Molly Lewis

A red line connects the notification text to the explanatory text above. A print icon in the top right corner of the page is connected to the 'Print current page here' text.

Specialist portal: Contents

3 Patient follow-up Patients

Click here to see detailed patient data

Average % spent in supine position during analysis phase

Click here for sorting (ascending or descending)

The screenshot shows the LUNOA Patients portal interface. On the left is a navigation menu with options: Home, Patients, My Profile, Inbox, About, and Logout. The main content area is titled 'Patients' and features a search bar with date filters (25-12-2018 to 23-01-2019) and a 'Show Details' button. Below the search bar is a table with the following data:

Full Name	Baseline % supine (first 2 nights of use)	% supine (weighted average of all days*)	Hours used (average of all days*)	% used at least 4 hours (all days*)
Pedro Caet	55	0	00:00	0
Realtime Aging 1	25	0	00:00	0
Realtime Aging 2	27	0	00:00	0
Realtime Aging 3	27	0	00:00	0
Realtime Aging 4	29	0	00:00	0
Realtime Aging 5	27	0	00:00	0
Realtime Aging 6	25	0	00:00	0
Realtime Aging 7	25	0	00:00	0
Realtime Aging 8	25	0	00:00	0
Saskia VdL	72	0	00:00	0


* Calculated over the selected time period in the calendar boxes

Important:
The first two nights (analysis phase) are excluded from all calculations (% Supine, # Days Used, Compliance)

Day N—starts at noon Day N-1 and finishes at noon on Day N

Specialist portal: Contents

4 Periodic view Patients > Show details



Specialist portal: Contents


5 Day view

Patients > Show details

Adjust the day here

Number of vibrations

Number of reactions



Detailed view of body positions through the night

Specialist portal: Contents

6 Dashboard view single patient

Patients > Show details

The screenshot displays the LUNOA patient dashboard for Sem Mekking. The interface includes a sidebar with navigation options (Home, Patients, My Profile, Inbox, About, Logout) and the NightBalance logo. The main content area shows an overview for the period 24-12-2018 until 23-01-2019. Three key metrics are highlighted with callouts:

- Percent of time slept in supine:** 0%
- Hours used in last session:** 00:00
- Percent of time responded to vibration:** No vibrations

Metric	Value
AVERAGE % SUPINE	0 %
AVERAGE USAGE	00:00
% REACTIONS	No vibrations



©2019 Koninklijke Philips N.V. All rights reserved.
Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

Caution: Federal law restricts this device to sale by or
on the order of a physician.

www.philips.com/respironics

RRDPGH EB 2/15/19 MCI 4108706 PN 1141599
1010 Murry Ridge Lane, Murrysville, PA 15668
800-345-6443 • 724-387-4000