

Notes from today's visit

Brush with toothpaste and a Philips Sonicare power toothbrush designed to help with:

☐ Sensitive teeth ☐ Gum care ☐ Plaque removal ☐ Stain removal

Floss your teeth each day with:

☐ String floss ☐ Floss picks ☐ Power/water flosser

Rinse with a mouthwash that helps with or contains: _____
and use it ____ time(s) a day.

Improve your oral health habits by setting a goal. Focus on:

Create healthy oral care habits at home!

Remember to brush your teeth twice a day for 2 minutes with a Philips Sonicare power toothbrush and don't forget to floss every day.

Thank you for visiting our office today!



PHILIPS

sonicare ZOOM!

Your Path to Health recommendation



DiamondClean Smart 9000 series

- Up to 20x more plaque removal¹
- Up to 15x healthier gums²
- Removes up to 100% more stains in less than 2 days³



ExpertClean 7000 series

- Removes up to 10x more plaque in hard-to-reach areas⁴
- Up to 7x healthier gums in just 2 weeks⁴



ProtectiveClean 4500/4700 series

- Improves gingival health up to 100% more¹
- 4500 removes up to 6x¹ more plaque
- 4700 removes up to 10x¹ more plaque



Philips Sonicare For Kids

- 98% of parents surveyed say it's easier to get their kids to brush longer and better⁵
- Removes up to 75% more plaque⁶



Philips Sonicare Power Flosser/Cordless Power Flosser

- Quad Stream technology creates four wide streams to remove up to 99.9% of plaque in treated areas⁷



Philips Sonicare brush heads



Premium All-in-One (A3)

Removes up to 20x more plaque¹



Premium Plaque Control (C3)

Removes up to 10x more plaque⁸



Premium Gum Care (G3)

Up to 7x healthier gums in just 2 weeks⁹



InterCare (I)

Removes up to 7x more plaque¹



Optimal Gum Care (G2)

Improves gum health up to 4x better¹



DiamondClean (W)

Up to 100% less stains in a week¹⁰



Sensitive (S)

Ultra-soft bristles for gentle, effective cleaning



Optimal Plaque Control (C2)

Removes up to 10x more plaque¹

1 vs. a manual toothbrush.

2 In Clean mode vs. a manual toothbrush; measured by MGI at 6 weeks.

3 In a lab study vs. a manual toothbrush; actual results may vary.

4 vs. a manual toothbrush, when used in Gum Health mode.

5 vs. a manual toothbrush alone.

6 vs. a children's manual toothbrush.

7 In an in-vitro study, actual results may vary.

8 From hard-to-reach areas compared to a manual toothbrush.

9 In Gum Health mode vs. a manual toothbrush; measured by GBI.

10 In White+ mode vs. a manual toothbrush.