

Critical connections. Care for life.

Philips ALS Solutions enable secure data transfer throughout the patient journey. Advanced monitoring includes a broad range of clinical measurements. The unique data capabilities support improved accuracy of records and hand-offs. Data management, monitoring, advisory tools and post-event analysis can improve collaboration and workflow. Scalable to meet changing needs, Philips ALS Solutions partner with customers throughout the life cycle of the relationship to offer seamless connected care.

	 Efficia DFM100 with HeartStart Event Review Pro	 HeartStart Intrepid with Philips Emergency Care Informatics Suite
Size	23.5 cm (H) x 29 cm (W) x 20.5 cm (D)	24.6 cm (H) x 29 cm (W) x 21 cm (D)
Weight (with battery)	~6.1 kg	~6.7 kg
Display size (diagonal)	17.8 cm	21.3 cm
Display resolution	800 x 480 pixels	1024 x 768 pixels
IP rating	IP54	IP54
SpO ₂ , NBP, EtCO ₂	Yes (EtCO ₂ external)	Yes (EtCO ₂ integrated)
Temp	No	Yes
12-Lead ECG	No	Yes
QCPR	No	Yes
Infant/Child AED	Yes	Yes
Patient category button	Yes	Yes
TBI advisory	No	Yes
Battery capacity (per battery)	100 shocks/2.5 hours monitoring	100 shocks/5 hours monitoring
Charge time on battery	≤5 secs to 150J ≤6 secs to 200J	≤5 secs to 150J ≤6 secs to 200J
WiFi	No	Yes
Cellular	No	Yes
Data Management Software	HeartStart Event Review Pro	Emergency Care Informatics Suite