

Premium package

Includes

Ultrasound system

Affiniti CVx Advanced with AutoStrain and stress PLUS X5-1 and xPlane

Philips understand the challenges you face. From a changing patient population to the need to maximise resources, the challenges can be considerable. **Philips Affiniti** helps you overcome these daily challenges so you can provide the best possible care for your patients every day. Affiniti offers a powerful combination of performance and workflow advances to enhance diagnostic confidence making it the choice of clinicians worldwide.



Stress test and equipment

Trackmaster TMX-428 Treadmill



Image management solution

TOMTEC ARENA Network



TOMTEC-ARENA network is a premium solution for private practice. Expand your access and workflow efficiency with a client-server network that offers concurrent user licenses, full DICOM connectivity, automated archiving, and measurement mapping services. At the core is TOMTEC DATACENTER¹: a powerful multi-vendor, multimodality platform managing connectivity and communication between IT systems, data archiving, PACS/EMR integration and much more. TOMTEC-ARENA network offers off-cart image review including STRESS ECHO², advanced quantification and structured reporting. We offer “small business” pricing models tailored to the needs of private practices, with a selection of clinical tools including premium TOMTEC service with an Up-to-Date maintenance contract to ensure you benefit from the latest developments.

1. Is part of TOMTEC-ARENA. TOMTEC-ARENA is a trademark of TOMTEC Imaging Systems GmbH.

2. STRESS ECHO is called ECHO-COM and is part of TOMTEC ARENA. TOMTEC-ARENA is a trademark of TOMTEC Imaging Systems GmbH.

This device is a heavy-duty, medical treadmill specifically designed to be interfaced and controlled by a variety of cardiac and pulmonary stress test systems. Design features such as the handrails have been used to provide an ergonomic and safe experience. The safety tether can be worn by patients and stops the treadmill when pulled.

Stress system – Philips ST80i Stress Testing System (with Tango)

Philips ST80i Stress Testing System opens up new paths for productivity and decision-making. Its wireless design, bi-directional connectivity, and advanced decision-making support tools enhance data review, patient care, and efficiency.



Equipment additions

Philips DigiTrak XT 1810 Basic Suite Software Package with 2x24HR and 2x48HR DigiTrak Digital Holter Recorders inc PC installation and training

Philips DigiTrak XT is the lightest Holter recorder on the market with a large on-board display. This sleek, streamlined recorder is designed for patient comfort and for streamlined workflow with all users in mind.



HeartStart Intrepid

Intrepid defibrillator/monitor – 867172 Base unit 3-lead ECG trunk cable, 2.7m 3-lead ECG leadset, manual and AED modes and Cardioversion.

You never know what you will face when you arrive at the scene of an emergency, or when a hospital patient will suffer sudden cardiac arrest. But you do know you need to be ready. **Philips HeartStart Intrepid** is up to the challenge. This easy-to-use, lightweight and rugged monitor/defibrillator gives you the power to respond quickly and act confidently when your patient's life is on the line.




Borsam ABP-06 ABPM

- Motion tolerance filtering
- Good consistency in ambulatory circumstance
- Accurate reading including body position recording
- Gravity sensor inside to determine 4 body position: sit/stand, lying, slight motion, heavy motion
- Hypertension PC software allows for professional analysis of patient results
- PDF export for your reports



PageWriter TC10

The **Philips PageWriter TC10** is an advanced cardiograph that can grow with you as your workflow evolves. With its compact size and lightweight design, the TC10 is a portable solution that supports patient needs in a wide range of care settings.



Furniture

Adjustable Echo Scanning Couch with cut-out

A high quality 2 section unit specifically designed to enhance patient-doctor interaction for echo sonographers and cardiologists.

Features:

The AMC 2540 Malachite Echo Scanning Couch has fully motorised movements of both height and backrest with a gas ram assisted echo drop down flap. The couch is supplied with a foot switch control for hands free operation. Easy access under couch for lifting machine arms.



Adjustable Scanning Chair

- Electric Hi Lo operation by foot switches
- Sleek design with dual motors
- Minimum height 55cm to maximum height 90cm in 30 seconds
- 50mm retractable castors for excellent mobility
- Adjustable backrest featuring double gas struts, angle adjustable from 0 to +60 degrees
- Able to lift up to 250kgs working weight
- Chair goes fully horizontal for people recovering in emergency situations



Adjustable Scanning Couch

- 2 section couch with electric Hi Lo operation by foot switch
- 75mm retractable, dual locking castors for safety and mobility
- Sleek design with rounded frame legs for safety
- Minimum height 50cm to maximum height 90cm in 30 seconds
- Adjustable head section through dual gas operated struts from -20 to +80 degrees
- 250kg loading weight



Consumables

We offer a wide range of consumables to meet your practice needs

- ECG electrodes
- Gel and wipes
- ECG Skin Preparation Paper

