

Standard package

Includes

Ultrasound system

Affiniti CVx Standard with AutoStrain LV, Stress and S5-1 transducer

Philips understand the challenges you face. From a changing patient population to the need to maximise resources, the challenges can be considerable. **Philips Affiniti** helps you overcome these daily challenges so you can provide the best possible care for your patients every day. Affiniti offers a powerful combination of performance and workflow advances to enhance diagnostic confidence making it the choice of clinicians worldwide.



Stress test and equipment


Trackmaster TMX-428 Treadmill



Image management solution

TOMTEC-ARENA

Single Seat Workstation for Image Review



The **TOMTEC-ARENA single seat workstation** provides the basis for off-cart image review and advanced quantification of ultrasound images and 3D datasets. Expand the capabilities of your single seat workstation by adding **REPORTING¹**: a comprehensive, customisable solution for structured reporting in Adult, Paediatric and Vascular ultrasound labs. The TOMTEC-ARENA single seat workstation is a complete IT solution with full DICOM connectivity. Licensing options are tailored to meet your specific practice needs: affordable, flexible and scalable. TOMTEC Up-to-Date Support and Maintenance ensures you benefit from the premium TOMTEC service and latest product developments and releases.

To learn more, please visit www.tomtec.de or ask to speak with a TOMTEC Sales Solution Specialist today.

1. Is part of TOMTEC-ARENA. TOMTEC-ARENA is a trademark of TOMTEC Imaging Systems GmbH.

With a 227kg weight limit this device is a heavy-duty, medical treadmill specifically designed to be interfaced and controlled by a variety of cardiac and pulmonary stress test systems. With a zero-start speed and unique Master-Track belt alignment system the Trackmaster is the pinnacle of medical treadmills.

Design features include inbuilt hand rails and both safety tether and manual emergency stop for your utmost patient safety.

Philips ST80i Stress Testing System

Philips ST80i Stress Testing System opens up new paths for productivity and decision-making. Its wireless design, bi-directional connectivity, and advanced decision-making support tools enhance data review, patient care, and efficiency.



Equipment additions

Philips DigiTrak XT 1810 Basic Suite Software Package with 2x24HR and 2x48HR DigiTrak Digital Holter Recorders inc PC installation and training

Philips DigiTrak XT is the lightest Holter recorder on the market with a large on-board display. This sleek, streamlined recorder is designed for patient comfort and for streamlined workflow with all users in mind.



HeartStart FRx AED Defibrillator (includes case, battery and pads)

The **Philips HeartStart FRx defibrillator** features intuitive, step-by-step voice instructions, including CPR guidance, and an audible metronome to help guide basic life support (BLS) responders while treating a suspected sudden cardiac arrest (SCA). Pre-connected SMART Pads II can be used for both adults and children. Rugged, lightweight and reliable, it can withstand rough handling and extreme temperatures. When every minute counts.




Borsam ABP-06 ABPM

- Motion tolerance filtering
- Good consistency in ambulatory circumstance
- Accurate reading including body position recording
- Gravity sensor inside to determine 4 body position: sit/stand, lying, slight motion, heavy motion
- Hypertension PC software allows for professional analysis of patient results
- PDF export for your reports



PageWriter TC10

The **Philips PageWriter TC10** is an advanced cardiograph that can grow with you as your workflow evolves. With its compact size and lightweight design, the TC10 is a portable solution that supports patient needs in a wide range of care settings.



Furniture

Adjustable Echo Scanning Couch with cut-out

A high quality 2 section unit specifically designed to enhance patient-doctor interaction for echo sonographers and cardiologists.

Features:

The AMC 2540 Malachite Echo Scanning Couch has fully motorised movements of both height and backrest with a gas ram assisted echo drop down flap. The couch is supplied with a foot switch control for hands free operation. Easy access under couch for lifting machine arms.



Adjustable Scanning Chair

- Electric Hi Lo operation by foot switches
- Sleek design with dual motors
- Minimum height 55cm to maximum height 90cm in 30 seconds
- 50mm retractable castors for excellent mobility
- Adjustable backrest featuring double gas struts, angle adjustable from 0 to +60 degrees
- Able to lift up to 250kgs working weight
- Chair goes fully horizontal for people recovering in emergency situations



Adjustable Scanning Couch

- 2 section couch with electric Hi Lo operation by foot switch
- 75mm retractable, dual locking castors for safety and mobility
- Sleek design with rounded frame legs for safety
- Minimum height 50cm to maximum height 90cm in 30 seconds
- Adjustable head section through dual gas operated struts from -20 to +80 degrees
- 250kg loading weight



Consumables

We offer a wide range of consumables to meet your practice needs

- ECG electrodes
- Gel and wipes
- ECG Skin Preparation Paper

