

Ingenia Elition X

MR Systems

A revolutionary breakthrough in diagnostic quality and speed

A revolutionary breakthrough in diagnostic quality and speed

Every day, healthcare moves forward with innovations in clinical pathways and supporting technologies. For radiology the necessity for high productivity, an improved patient experience while ensuring excellence in imaging can be daunting. The perception is often that MR represents a trade-off between productivity and image quality. The Philips Ingenia Elition X solution offers cutting-edge MR imaging techniques, while setting new directions for clinical research in 3.0T imaging based on gradient- and RF designs.

The Ingenia Elition X delivers on superb image quality, and performs MRI exams up to 50% faster¹. Fast overall exam-time is achieved by improving patient handling setup time at the bore with the touchless guided patient setup, combined with accelerations in both 2D- and 3D scanning. This has been made possible by gradient- and RF designs as well as Compressed SENSE. Furthermore, the Ingenia Elition X offers an immersive audio-visual experience to calm patients and guide them through MR exams. In a study, with the use of the in-bore solution, Herlev Gentofte University Hospital in Denmark managed to reduce the number of rescans by up to $70\%^2$, allowing radiologists to handle more patients per day.

Ingenia Elition X

A confident d

Up to 50% fa

Patient-cen

Dramaticall

Enhance the

iagnosis	6
ster MRI exams with virtually equal image quality ¹	14
ered productivity	22
mprove patient experience	32
value of your MR investment	40





1 Compared to Philips scans without Compressed SENSE.

A **confident** diagnosis

- Achieve **up to 60%** higher resolution¹
- ▶ Scan your DWI images **up to 30%** faster²
- Enhance your diagnostic confidence in neuro oncology
- Unlock new territories in neurofunctional MRI

Up to 50% faster MRI exams¹

- ► Virtually equal image quality
- ► Add extra patient slots in your schedule
- ► Easily fit in unplanned patients
- ▶ Reduce overtime



1 Compared to Philips scans without Compressed SENSE 2 Compared to Ingenia 3.0T Omega HP R5.3

Patient-centered productivity

- ▶ Reduce operator workload
- ▶ Guided exam set-up
- Automated workflow
- Superior image quality

Dramatically improve patient experience

- ► Provide an immersive visual experience
- Guide your patients through the examination
- Comfort in every detail
- ▶ Reduce acoustic noise

Enhance the value of your MR investment

- ▶ Prevent issues before they occur
- ▶ Protect your MR equipment
- ► Standardize your MR fleet
- ► Tailored financing solutions



A confident diagnosis


Ingenia Elition X supports confident diagnosis by innovating on all fronts. The new high-end performance gradient and RF design, combined with innovative imaging solutions such as Compressed SENSE and 3D APT help you reach new levels of precision in anatomical and functional clinical imaging. As a result, you can diagnose the most challenging clinical indications with confidence.

Up to 60% higher resolution

Fully redesigned gradients combined with Compressed SENSE acceleration technology allow up to 60% higher spatial resolution¹ in the same scan time, revealing more details. For example, speeding up isotropic 3D MSK imaging enables a switch from multi-orientation 2D imaging to a single high-resolution 3D efficient scan. It's precision made efficient.

Up to 30% faster DWI images

Ingenia Elition X's high-performance Vega HP gradients allow the most advanced imaging techniques on a 3.0T system. Diffusion scans are up 30% faster while appearing sharper². An average of 70% higher contrast resolution can be achieved in diffusion imaging². Due to a TE up to 15% shorter in diffusion imaging, SNR is further improved or used to generate higher resolution with similar scan time². All this makes Ingenia Elition X an ideal choice in clinical routine and oncology applications, even for challenging anatomies.



Enhanced diagnostic confidence in neuro oncology

3D APT (Amide Proton Transfer) is a unique, contrast-free, brain MR imaging method that addresses the need for more confident diagnosis in neuro oncology. 3D APT uses the presence of endogenous cellular proteins to produce an MR signal that directly correlates with cell proliferation, an indicator of tumoral activity. 3D APT can support trained medical professionals in differentiating low grade from high grade brain gliomas and, in differentiating tumor progression from treatment efffect³.

Unlock new territories in neurofunctional MRI


Ingenia Elition X Vega HP gradients deliver up to 23% higher temporal resolution in fMRI studies as well as 30% shorter TR in diffusion imaging for excellent functional imaging at 3.0T². Ingenia Elition X lets you unlock new territories in the field of neurofunctional MRI, unraveling the connections and functional set-up of the brain. This new level of precision and our next-generation wide bore system can attract an array of new funding with human connectome style protocols, and open up research opportunities.


One processing platform for making your diagnosis

IntelliSpace Portal offers a comprehensive set of over 70 clinical applications in multiple clinical domains, including neurology, cardiology, vascular, oncology, and more. This all-round workstation provides you with the right tools when you need them.


1 In isotropic 3D MSK VIEW scans, compared to Philips scans without Compressed SENSE. 2 Compared to Ingenia 3.0T Omega HP R5.3 3 Togao et al. (2014) Neuro-Oncology. Park KJ et al. (2016) Eur Radiol

















NeuroScience


Resting State fMRI - ABCD research protocol











































Coronal T1w TSE mDIXON XD with gado (In Phase + Water only 0.4 x 0.4 x 2.0 mm, 3:46 min Courtesy: Academic Medical Center, Amsterdam, The Netherlands







Sagittal PDw TSE 0.18 x 0.18 x 1.5 mm, 5:05 min



Axial PDw TSE 0.23 x 0.35 x 2.5 mm, 1:53 min




Axial 3D PDw SPAIR 0.8 x 0.8 x 1.0 mm, 6:54 min



oronal PDw TSE 9.3 x 0.3 x 1.5 mm, 5:05 min




Sagittal PDw SPAIR 0.3 x 0.4 x 2.5 mm, 4:06 mir




Axial T2w Axial DWI Courtesy: Te




Axial T1w Axial DWI Courtesy: Te






SE, 0.46 x 0.57 x 3.0 mm, 5:07 min 00 / b500 / b1000 / b1400, 2.5 x 2.4 x 3.0 mm, 3:44 min Inrical University Munich, Germany



D mDIXON, 1.5 x 1.5 x 2.0 mm, 2:50 min o0 / b50 / b300 / b600, 3.0 x 3.0 x 4.0 mm, 4:03 min hnical University Munich, Germany



Coronal T2w TSE – Two stations .4 x 1.6 x 5.0 mm, 1:48 min /station



Up to 50% faster MRI exams with virtually equal image quality¹

Time is one of the most precious commodities you have in your MR department.

What if we told you there was a way to recover time you have been losing during your MR examinations? And use the time you do have more wisely? Imagine how that could help you make better use of your scarce resources and better meet the demands of referring physicians.

That's exactly what Compressed SENSE can do for your MR department. It accelerates your existing MR scans by up to 50% with virtually equal image quality, frees up time to improve your patient experience and can provide up to 60% higher resolution to enhance diagnostic confidence.¹

Compressed SENSE is suitable for all anatomies and can be used for all anatomical contrasts, in both 2D- and 3D scanning.

Up to 50% faster MRI exams with virtually equal IQ¹ Brain ExamCard





Total Exam time **16:52 min**

/ SE	0.9 x 1.3 x 5.0 mm	2:35
v TSE	0.6 x 0.7 x 5.0 mm	2:04
/I (b1000)	1.8 x 1.4 x 5.0 mm	0:45
v FLAIR	1.0 x 1.2 x 5.0 mm	2:12
v FFE	0.9 x 1.1 x 5.0 mm	1:15 r
Inflow	0.7 x 1.3 x 1.4 mm	7:11 r

Total Exam time **9:06 min** Compressed SENSE

w SE	0.9 x 1.3 x 5.0 mm	1:14 min
2w TSE	0.6 x 0.7 x 5.0 mm	1:30 min
WI (b1000)	1.8 x 1.4 x 5.0 mm	0:45 min
2w FLAIR	1.0 x 1.2 x 5.0 mm	1:30 min
2w FFE	0.9 x 1.1 x 5.0 mm	0:49 min
D Inflow	0.7 x 1.3 x 1.4 mm	3:06 min





PDw SPAIR sag PDw TSE sag T2w TSE ax

Knee ExamCard

Total Exam time **17:32 min**


0.4 x 0.4 x 2.5 mm	4:24 m
0.3 x 0.4 x 1.5 mm	3:55 m
0.4 x 0.4 x 1.5 mm	3:49 n
0.4 x 0.4 x 1.5 mm	2:23 m
0.4 x 0.5 x 2.5 mm	3:01 m

Total Exam time **10:16 min** Compressed SENSE

PDw SPAIR sag	0.4 x 0.4 x 2.5 mm	3:02 mir
PDw TSE sag	0.3 x 0.4 x 1.5 mm	2:20 mir
T2w TSE ax	0.4 x 0.4 x 1.5 mm	1:52 min
T1w TSE cor	0.4 x 0.4 x 1.5 mm	1:13 min
PDw SPAIR cor	0.4 x 0.5 x 2.5 mm	1:49 min

Add extra patient slots to your daily MRI schedule

Many radiology departments and imaging centers are looking for ways to increase the utilization of their MR equipment to meet the rising demand for MRI services. A full MRI exam performed with Compressed SENSE, for example, can save minutes compared to a conventional MRI exam. This could free up one or two extra exam slots in your daily schedule, which can result in much higher productivity and shorter waitlists without adding more operator hours.



Radiologie Dr Wagner in Gottingen , Germany can accommodate > 5 more patients per day, within the same scanning hours, after the introduction of Compressed SENSE.

Easily fit in unplanned patients

Do unscheduled patients disrupt your daily schedule and put extra stress on your staff? With Compressed SENSE you can create a buffer to easily handle emergency cases or urgent patients that are referred on the same day. This extra capacity can help you serve patients and referring physicians faster and make daily workflow go smoother.

⁶⁶ We can now provide a more flexible and faster MRI service to our patients and referring physicians. For instance, when a referring physician is requesting it, we can now quite smoothly insert an additional MRI examination without previous appointment on the same day."


Hideki Koyasu, MD, Neurosurgical Clinic in Kanagawa, Japan






Reduction in overtime, while maintaining same patient throughput per day

Having to work overtime is a recurring issue for many radiology departments and imaging centers that can impact staff satisfaction and run-up operational costs. The stress caused by heavy workloads and overtime hours greatly contribute to burnout among radiology technologists, not to mention long-term mental and physical health issues.¹ By reducing MRI scan times and improving scheduling flexibility, Compressed SENSE helps patients and staff to get home on time. This can improve the experience for all involved.



ComputerTomography Institut in Innsbruck, Austria has been able to reduce overtime by more than one hour, keeping the same patient throughput per day, after the introduction of Compressed SENSE.

1 Vinu, Raj. Occupational stress and Radiography. NCBI. Nov-Dec 2006. https://www.ncbi.nlm.nih.gov/pubmed/17119177 Results from case studies are not predictive of results in other cases. Results in other cases may vary.



Patient-centered productivity


With a growth in the elderly population and constant demands to do more with less, the pressure on healthcare providers is immense. This pressure is also evident in radiology departments and imaging centers. The increasing use of MR to diagnose a variety of conditions and illnesses has led to demands for greater efficiency, even as departments try to manage a shortage of MR operators and variability in staff expertise. Too often, it seems that productivity is at odds with giving patients the time and attention they desire.

With SmartWorkflow, you can achieve high productivity while enabling your staff to focus on patients. It reduces and simplifies the number of steps needed in a conventional MR exam workflow, using technology to guide and coach where required, and automate where possible. An end-to-end workflow solution that directly boosts efficiency through reduced variability and task automation, while supporting a better patient and staff experience, resulting in patient-centered productivity.





Auto patient centering Put your patients at ease, while manual steps in the workflow are automated



In-room exam start Increase productivity and free up time for other tasks²

Guided exam set-up A virtual coach guiding ex set-up; allowing your staff



Touchless respiratory-triggering Remove the hassle of respiratory belt placemen while providing superior image quality¹



SmartWorkflow in **the exam room**

Guided exam set-up and automation, to increase productivity and free up time to focus on the patient

In the exam room, SmartWorkflow provides guided exam set-up and automation, to increase productivity and free up time to focus on the patient. Even new operators who have never worked with the scanner can proceed with confidence. Allow your staff to focus less on technology, and fully engage with patients. Enjoy reduced variability in patient positioning and quality of respiratory triggering, supporting consistent, high-quality studies. SmartWorkflow decreases patient set-up to less than a minute³, and allows operators to initiate the start of the exam with a single touch directly at the MR scanner, starting immediately after closing the door.

⁶⁶ The entire workflow is smooth: Patient positioning and setup; launching the scan as soon as we leave the exam room; the intuitive touchscreen on the gantry; Touchless patient sensing... All of these things are much better than on our old system."

Lauro Barlow, RTMR, MRI Technologist Supervisor at the University of British Columbia

1 Compared to respiratory belt based signal. Requires an unobstructed line-of-sight.

2 With VitalScreen, initiating the exam at the patient side

3 Based on in-house testing.

SmartWorkflow in the control room

Reduce operator workload, with standardized results and increased throughput


In the control room, SmartWorkflow automates exam planning, scanning and processing, improving staff experience and driving efficiency by freeing time to check imaging results or prepare for the next exam. Decreased exam variability results in imaging excellence supporting confident diagnoses while automated patient coaching enhances patient experience. Furthermore, with SmartWorkflow you can confidently offer imaging to patients with MR Conditional implants. The guided and automated workflow supports the staff, which gives them more time for the patient, resulting in patientcentered productivity.

⁶⁶ We don't have to manually direct the patient to breathing and not breathing. We can go ahead and let the machine do the work of the breathing instructions while we continue our planning of the exam."

Carlos Avila, RT, technologist at Miami Cardiac & Vascular Institute



and scanning



Automated patient coaching




Confidence for MR Conditional implants





post-processing





ngenia Elition X

A virtual coach guiding exam set-up

Increase staff confidence and speed up patient set-up through automated real-time guidance and insights on the details of the current patient study. Achieve high quality results, independent from staff's expertise level. VitalScreen provides guidance at your staff's fingertips. Two 12-inch interactive touchscreens on the scanner provide coaching and visual guidance on recommended patient position, study laterality, coil and accessory placement. Moreover, feedback is provided on important exam details, including physiology signals (both VCG and respiratory) and – if applicable– contrast usage and breath– hold guidance.

Put your patients at ease, while manual steps in the workflow are automated

Free up your staff from monotonous, manual steps and enable them to focus on the patient through automatic placement of the region of interest in the scanner iso-center. The manual use of a laser light visor for iso-center positioning has become obsolete. VitalScreen automatically detects landmarks for selected anatomies and places the region of interest in the iso-center of the magnet. Once the patient is positioned on the table, only the push of a button is required to position the patient in the center of the bore.

Increase productivity and free up time for other tasks


Start exams as soon as possible, eliminating extra steps for your staff and decreasing the time the patient has to spent in the magnet, resulting in a more positive patient experience. VitalScreen allows staff to initiate the exam with a single touch of a button at the patient side. The exam starts immediately after the operator has closed the exam room door, so no time is wasted.

Continuous and robust respiratory signal providing **superior image quality**

Relieve your staff from the burden of positioning – and re-positioning – a respiratory belt. Positioning a belt shifts the operator's focus from the patient to the technology at a moment when it is critical that the patient is comfortable and reassured. Enjoy optical sensing and Al¹ to automatically detect patient respiratory patterns. VitalEye touchless patient sensing provides a fast detection of patient's breathing without any operator interaction. With VitalEye, the technologist no longer needs to set up an old-fashioned respiratory belt but receives a continuous and robust respiratory signal without any interaction. This revolution in touchless patient sensing helps your staff to keep a caring eye on your patient. The quality of the physiology signal detected by VitalEye is better than a belt-based approach providing superior image quality, for a broad range of patient sizes.


⁶⁶ It always works, and it's always there."

MR operator, University of Bonn, Germany



Touchless patient sensing

Superior image quality with VitalEye¹, consistently



1 Compared to Philips belt-based signal. Requires an unobstructed line of sight Results from case studies are not predictive of results in other cases. Results in other cases may vary.



Dramatically improve the **patient experience**

Your patients are at the heart of Ingenia Elition – which includes an MR experience that enhances comfort and compliance. With up to 80% acoustic noise reduction1, voice guidance, Immersive in-bore visuals and a comfortable table mattress, Ingenia Elition helps your patients feel at ease, resulting in smooth, fast exams.

Provide an immersive visual experience

Your patients' scanning experience is significantly enhanced with Ingenia Elition. Designed to offer a relaxing sensory experience, Ambient Experience provides positive distractions for patients by incorporating dynamic lighting, projection and sound, contributing to a positive, engaging environment to benefit quality of care. From the moment a patient is moved into the scanner (the point at which people report the most stress), through completion of the scan, the In-Bore Connect solution can help patients to relax, follow directions and minimize motion. In a study, conducted using our in-bore solution, Herlev Gentofte University Hospital in Denmark managed to reduce the number of rescans by up to 70%¹. A case study at Radiologisches Zentrum am Kaufhof, Lübeck, Germany showed that the number of patients needing sedation was reduced by 80%².

⁶⁶ We've had a lot of patients provide compliments on the environment. We have the Ambient solution in there that creates a soothing environment."

Carol Melvin, MD, Miami Cardiac and Vascular Institute

1 Compared to the average of the other 5 Philips Ingenia MR scanners without Ambient Experience and In-Bore Connect. Results from case studies are not predictive of results in other cases. Results in other cases. Results in other cases may vary. 2 Results from case studies are not predictive of results in other cases. Results in other cases may vary. *The tranquilizer referred to is a valium-based derivative called "Diazepam".





Comfort in every detail

Because no detail is too small when it comes to helping your patients feel comfortable, Ingenia Elition includes the ComfortPlus mattress. On average, 90% of patients in severe discomfort find it easy to lie still on the ComfortPlus mattress. Overall comfort for this group of patients can increase by up to 36%.¹

Put your patients at ease and guide them through the examination

Be confident that your patients know what to do and what to expect through automated, consistent instructions, relieving some of the anxiety of an MR exam. AutoVoice supports exam compliance by guiding your patients through the MR examination. Including automatically announcing scan duration and table movements in your choice of 30 languages and dialects. In addition, providing breath hold instructions, with either manual timing or timing synchronized to fit the patient's respiratory cycle.


⁶⁶ The most frequent comment we are getting from our technologists, is that for patients who have had scans on other Philips scanners, this new mattress is really significantly more comfortable." ⁶⁶ AutoVoice enables us to shift our focus from having to manually give the breathing instructions to the patient to now planning the exam."

Carlos Avila, RT, Miami Cardiac and Vascular Institute, USA

Dr. Oswood, Hennepin County Medical Center

Reduce acoustic noise for your patient

No matter how short the exam, a noisy MR scanner can make it seem unbearably long. The Philips unique ComforTone solution achieves up to 80% reduction in acoustic noise¹ with similar image quality and contrast within the same time slot. You can use ComforTone in routine exams such as brain, spine and MSK. Thanks to our ready-to-use ExamCard protocols, ComforTone is simple to implement and use, requiring just a few clicks to get started.



1 Compared to scanning without ComforTone Results from case studies are not predictive of results in other cases. Results in other cases may vary.



Enhance the value of your MR investment


Imaging is both a clinical and an economic challenge. You need to manage a host of financial obligations and opportunities, all while keeping your focus on your patients. We can help, by putting together a package of offerings that keep total cost of ownership in check while providing you with tailored solutions for maintenance, fleet management, cybersecurity, education and financing.



500 parameters are monitored on an MR


+90,000remote connections across 25,000healthcare facilities in

139 countries³



connected Philips MR service cases² are resolved before they cause downtime, due to proactive monitoring





1 For the Philips diagnostic imaging installed base

2 Based on data collected between July '18 and July '19 on all service events registered on remotely connected Philips MR systems [globally]. Downtime does not include time due to planned maintenance 3 Based on global Philips only data.

Prevent issues before they occur

Scanner downtime can disrupt your schedule and delay patient care. We offer maintenance agreements that are suited to your needs, enabled by the latest service innovations and including an uptime guarantee. We prevent issues before they occur through proactive remote monitoring, remote diagnostics and remote and field service support. With e-Alerts and other remote data, we monitor more than 500 parameters of your MR system from a distance, detecting and resolving issues without impacting your department's operations. In fact, more than 50% of MR service cases are resolved remotely.¹ Our Philips-qualified service experts can also proactively resolve issues on-site, fix your system before it causes any disruption, and provide reliable and knowledgeable support.²

Protecting your MR equipment from patient data breaches and cyber-attacks

Protecting patient health information requires constant vigilance. To keep health information and medical devices secure, we employ best practices in medical device security. Our multi-layered defense barriers include security policies, procedures, access controls, technical measures, training, and risk assessments. The Technology Maximizer Plus subscription program conveniently keeps your MR systems up-to-date through access to the latest cybersecurity patches and mandatory safety fixes via regular and ongoing software upgrades and hardware refreshes.

Standardize your MR fleet at a fixed annual cost

If you own more than one Philips scanner, standardizing under the same software release can enhance efficiency through one user interface for operators to learn and use the same ExamCards across multiple scanners. The Ingenia Elition is delivered with the latest available software release, providing a perfect opportunity to upgrade your fleet to this release and enter into a Technology Maximizer Plus subscription program.¹ Under the program, your Ingenia Elition and the rest of your fleet will receive software updates whenever available, giving you the benefits of software improvements and cyber-security advances while maintaining all your MR systems on the same level.

⁶⁶ It was consistently a challenge to plan for annual upgrades and predict their costs. Thanks to Technology Maximizer, we can now continue to have the latest versions of software for all of our MRI systems."

> Eliseo Vañó Galván, , MD, Cardiovascular radiologist, Chairman of the CT & MR Department at Hospital Nuestra Senora del Rosario, Madrid, Spain


in previous years




Cost of maintenance and upgrades



Reduction in accumulated cost of maintenance and upgrades over 4 years Before vs with Technology Maximizer





✓ Technology Maximizer saves cost and provides more frequent updates

Achieve excellence through **ongoing education**

Delivering consistent healthcare day-in and day-out is a challenge, particularly when faced with staff shortages and the need to cross-train department personnel. Our Philips MR Healthcare Education can help unlock the full potential of your staff, technology, and organization through innovative and meaningful healthcare education, delivered on-site or as e-Learning. For example, the Philips MR Technologist Development Program at Burjeel Hospital for Advanced Surgery (BHAS), a leading orthopedic and joint care center in Dubai, UAE, resulted in an average of 30% improvement in image quality across all procedures.¹ Team knowledge increased 30-40% in the key areas of patient care¹, imaging procedures, data acquisition and physics of image formation. The comprehensive, clinicallyrelevant courses, programs, and learning paths are designed to support clinical excellence, enhance operational efficiency and provide high-quality patient care.

Tailored financing solutions in line with your cash flow needs, budgets, and business strategy

Providing access to best-in-class healthcare is a leading priority for facilities like yours around the globe. At the same time, financial security and protecting your assets over time are also high on the agenda. To manage your financial challenges, you need to know whether your healthcare investments are sustainable – and how to get the most from your equipment. Financing your Ingenia Elition helps you exchange variability and unpredictability for visibility and certainty. This helps you avoid the burden and risk of upfront expenditures and benefit from transparent, predictable cost structures. As a result, you can manage and plan budgets more efficiently and free up capital that would otherwise be tied up in fixed assets.





© 2021 Koninklijke Philips N.V. All rights reserved. Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

4522 991 66081 * JAN 2021

How to reach us Please visit www.philips.com healthcare@philips.com