

Portable oxygen solutions

Veterans Affairs

SimplyGo



SimplyGo Mini



Oxygen on the go

Philips SimplyGo and SimplyGo Mini give veterans confidence on the go, 24/7/365

If you use portable oxygen, Philips offers veterans powerful choices to meet your needs, wherever and whenever you need them. With both SimplyGo and SimplyGo Mini, you can have the confidence of knowing you have a durable, lightweight and powerful portable oxygen concentrator at your side day and night.

SimplyGo Mini



Our smallest and lightest oxygen concentrator ever

But, don't let its size fool you. SimplyGo Mini is tough. At just five pounds, it's built to withstand almost anything you can throw at it – from drops, to falls, to rain showers and extreme temperatures.

Count on it night and day

SimplyGo Mini gives you the oxygen therapy you need, wherever and whenever you need it – 24 hours a day, 7 days a week. It has been shown to be as effective as continuous flow liquid oxygen in keeping oxygen levels at 90% or above.¹

Up to 9 hours of battery² operation

From a walk in the park, to a car ride – or even traveling by airplane – SimplyGo Mini offers lithium Ion battery options with standard AC/DC power cords that add convenience and flexibility for when and where you need it the most.



¹When used in accordance with all instructions and conditions outlined in the SimplyGo Mini user manual.

²Battery duration times are based on new, fully charged batteries used with a new SimplyGo Mini system. Battery duration times will degrade with battery age, environmental operating conditions, use over time, and operational condition of the concentrator.

SimplyGo



The all-in-one oxygen concentrator

SimplyGo is a portable and stationary oxygen concentrator. It offers pulse flow and continuous flow in a small and compact unit for easy storage. And, if you need higher pulse settings at night, the Sleep mode does it while providing exceptional comfort.

Transportable and powerful oxygen therapy

SimplyGo's easy-to-carry design makes it the right oxygen concentrator to take anywhere you go. And SimplyGo is the only continuous flow oxygen concentrator to weigh in at 10 pounds – more than twice the oxygen output of any POC weighing 10 pounds or less.

Our new standard in durability

SimplyGo was built with a long-life compressor, quality parts, oversized cart wheels and a high-impact-resistant design that sets our new standard in ruggedness.



Air travel

The manufacturer of these POCs has determined these devices conform to all applicable FAA requirements for POC carriage and use on board aircraft.



Learn more

Contact your provider today to learn more about these oxygen products or visit www.philips.com/oxygen

