


**GS1-128 Barcode Format**



- 1 Application Identifier, AI (01), indicates that the GS1 AI data field contains a GTIN
- 2 GTIN-14
- 3 AI (17) indicates that the GS1 AI data field contains expiration date
- 4 Expiration Date (YYMMDD, example: 160910)
- 5 AI(10) indicates that the GS1 AI data field contains Lot number
- 6 Lot number (XXXXXXXXXX, 9 characters)

**GTIN-14 Barcode Format**



















- 1 Indicator Digit for different levels of product hierarchy (For product with package of 1 unit per box: the 0 indicator is for both pouch and box; For product with package of multiple units: 0 is for pouch and 2 is for box)
- 2 GS1 Company Prefix (SPNC Prefix: 081313202)
- 3 Item Reference Number (example: 074)
- 4 Calculated Check Digit Reference (calculation based on previous 13 digits to ensure global uniqueness).

## Lead Management Product GS1 Barcodes

Model #	Description	GTIN Portion of the GS1 BOX Barcode
<b>Bridge</b>		
590-001	Bridge 80mm	 00813132024666
<b>Bridge (-B Model)</b>		
590-001-B	Bridge 80mm	 00813132029739
<b>Laser Sheath (SLS™II and Glidelight™)</b>		
500-001	SLS II Laser Sheath Kit 12F	 00813132020293
500-012	SLS II Laser Sheath Kit 14F	 00813132020309
500-013	SLS II Laser Sheath Kit 16F	 00813132020316
500-301	Glidelight 12F	 00813132020323

Model #	Description	GTIN Portion of the GS1 BOX Barcode
500-302	Glidelight 14F	 00813132020330
500-303	Glidelight 16F	 00813132020347
<b>GlideLight Laser Sheath (-B Models)</b>		
500-301-B	Glidelight 12F	 00813132027926
500-302-B	Glidelight 14F	 00813132027933
500-303-B	Glidelight 16F	 00813132027940
<b>GlideLight Laser Sheath (-C Models)</b>		
500-301-C	Glidelight 12F	 00810060550595

Model #	Description	GTIN Portion of the GS1 BOX Barcode
500-302-C	Glidelight 14F	 00810060550601
500-303-C	Glidelight 16F	 00810060550618
<b>Lead Locking Device (LLD™)</b>		
518-062	LLD EZ (Box of 3)	 20813132023076
518-067	LLD EZ (Box of 1)	 00813132023065
518-039	LLD E (Box of 3)	 20813132023052
518-018	LLD #1 (Box of 3)	 20813132022994
518-021	LLD #1 (Box of 1)	 00813132023003
518-019	LLD #2 (Box of 3)	 20813132023014

Model #	Description	GTIN Portion of the GS1 BOX Barcode
518-022	LLD #2 (Box of 1)	 00813132023027
518-020	LLD #3 (Box of 3)	 20813132023038
518-023	LLD #3 (Box of 1)	 00813132023041
<b>LLD Accessories</b>		
518-024	Lead Cutter	 00813132020804
518-027	Accessory Kit	 00813132020811
<b>Lead Locking Device (-B Models)</b>		
518-062-B	518-062-B	 20813132029924
518-067-B	518-067-B	 00813132029647
518-039-B	518-039-B	 20813132029917








Model #	Description	GTIN Portion of the GS1 BOX Barcode
518-018-B	518-018-B	 20813132029887
518-021-B	518-021-B	 00813132029616
518-019-B	518-019-B	 20813132029894
518-022-B	518-022-B	 00813132029623
518-020-B	518-020-B	 20813132029900
518-023-B	518-023-B	 00813132029630
<b>LLD Accessories (-B Models)</b>		
518-024-B	518-024-B	 00813132029654

Model #	Description	GTIN Portion of the GS1 BOX Barcode
518-027-B	518-027-B	 <p>00813132029661</p>
<b>Lead Locking Device (-C Models)</b>		
518-062-C	518-062-C	 <p>20810060550711</p>
518-067-C	518-067-C	 <p>00810060550724</p>
518-039-C	518-039-C	 <p>20810060550704</p>
518-018-C	518-018-C	 <p>20810060550629</p>
518-021-C	518-021-C	 <p>00810060550656</p>
518-019-C	518-019-C	 <p>20810060550636</p>





Model #	Description	GTIN Portion of the GS1 BOX Barcode
518-022-C	518-022-C	 00810060550663
518-020-C	518-020-C	 20810060550643
518-023-C	518-023-C	 00810060550670
<b>LLD Accessories (-C Models)</b>		
518-024-C	518-024-C	 00810060550687
518-027-C	518-027-C	 00810060550694
<b>VisiSheath™ Dilator Sheath</b>		
501-012	VisiSheath S 43cm	 00813132020354
501-014	VisiSheath M 43cm	 00813132020361



Model #	Description	GTIN Portion of the GS1 BOX Barcode
501-016	VisiSheath L 43cm	 00813132020378
501-112	VisiSheath S 33cm	 00813132020385
501-114	VisiSheath M 33cm	 00813132020392
501-116	VisiSheath L 33cm	 00813132020408
501-212	VisiSheath S 23cm	 00813132020415
501-214	VisiSheath M 23cm	 00813132020422
501-216	VisiSheath L 23cm	 00813132020439
<b>TorqMax™ Sheath Grip Accessory</b>		
502-001	TorqMax Grip Accessory	 20813132020440





Model #	Description	GTIN Portion of the GS1 BOX Barcode
<b>TightRail™ Rotating Dilator Sheath</b>		
545-509	TightRail 9F	 00813132020873
545-511	TightRail 11F	 00813132020880
545-513	TightRail 13F	 00813132020897
<b>TightRail™ Sub-C</b>		
560-009	<b>TightRail 9F</b>	 00813132022785
560-011	<b>TightRail 11F</b>	 00813132022792
560-013	<b>TightRail 13F</b>	 00813132022808
<b>TightRail Mini™ Rotating Dilator Sheath</b>		
540-009	TightRail Mini 9F	 00813132020859





Model #	Description	GTIN Portion of the GS1 BOX Barcode
540-011	TightRail Mini 11F	 00813132020866
<b>TightRail Guardian Motorized Dilator Sheath</b>		
575-009	TightRail Guardian 9F	 00813132027421
575-011	TightRail Guardian 11F	 00813132027438
575-013	TightRail Guardian 13F	 00813132027445
<b>SightRail™ Manual Dilator Sheath</b>		
550-008	Dilator Sheath 8.5F	 00813132020903
550-010	Dilator Sheath 10F	 00813132020910
550-011	Dilator Sheath 11.5F	 00813132020927
<b>SightRail™ Manual Dilator Sheath (Continued)</b>		
550-013	Dilator Sheath 13F	 00813132020934

Model #	Description	GTIN Portion of the GS1 BOX Barcode
555-508	Long Dilator Sheath 8.5F	 00813132020941
555-510	Long Dilator Sheath 10F	 00813132020958
555-511	Long Dilator Sheath 11.5F	 00813132020965
555-513	Long Dilator Sheath 13F	 00813132020972

Model #	Configuration	GTIN Portion of the GS1 POUCH Barcode
<b>Bridge</b>		
590-001	Bridge 80mm	 00813132024666
<b>Bridge (-B Model)</b>		
590-001-B	Bridge 80mm	 00813132029739
<b>Laser Sheath (SLS™II and Glidelight™)</b>		
500-001	SLS II Laser Sheath Kit 12F	 00813132020293
500-012	SLS II Laser Sheath Kit 14F	 00813132020309
500-013	SLS II Laser Sheath Kit 16F	 00813132020316
500-301	Glidelight 12F	 00813132020323

Model #	Configuration	GTIN Portion of the GS1 POUCH Barcode
500-302	Glidelight 14F	 00813132020330
500-303	Glidelight 16F	 00813132020347
<b>GlideLight Laser Sheath (-B Models)</b>		
500-301-B	Glidelight 12F	 00813132027926
500-302-B	Glidelight 14F	 00813132027933
500-303-B	Glidelight 16F	 00813132027940
<b>GlideLight Laser Sheath (-C Models)</b>		
500-301-C	Glidelight 12F	 00810060550595

Model #	Configuration	GTIN Portion of the GS1 POUCH Barcode
500-302-C	Glidelight 14F	 00810060550601
500-303-C	Glidelight 16F	 00810060550618
<b>Lead Locking Device (LLD™)</b>		
518-062	LLD EZ (Box of 3)	 00813132023072
518-067	LLD EZ (Box of 1)	 00813132023065
518-039	LLD E (Box of 3)	 00813132023058
518-018	LLD #1 (Box of 3)	 00813132022990
518-021	LLD #1 (Box of 1)	 00813132023003
518-019	LLD #2 (Box of 3)	 00813132023010

Model #	Configuration	GTIN Portion of the GS1 POUCH Barcode
518-022	LLD #2 (Box of 1)	 00813132023027
518-020	LLD #3 (Box of 3)	 00813132023034
518-023	LLD #3 (Box of 1)	 00813132023041
<b>LLD Accessories</b>		
518-024	Lead Cutter	 00813132020804
518-027	Accessory Kit	 00813132020811
<b>Lead Locking Device (-B Models)</b>		
518-062-B	LLD EZ (Box of 3)	 00813132029920
518-067-B	LLD EZ (Box of 1)	 00813132029647
518-039-B	LLD E (Box of 3)	 00813132029913









Model #	Configuration	GTIN Portion of the GS1 POUCH Barcode
518-018-B	LLD #1 (Box of 3)	 00813132029883
518-021-B	LLD #1 (Box of 1)	 00813132029616
518-019-B	LLD #2 (Box of 3)	 00813132029890
518-022-B	LLD #2 (Box of 1)	 00813132029623
518-020-B	LLD #3 (Box of 3)	 00813132029906
518-023-B	LLD #3 (Box of 1)	 00813132029630
<b>LLD Accessories (-B Models)</b>		
518-024-B	Lead Cutter	 00813132029654
518-027-B	Accessory Kit	 00813132029661
<b>Lead Locking Device (-C Models)</b>		

Model #	Configuration	GTIN Portion of the GS1 POUCH Barcode
518-062-C	518-062-C	 00810060550717
518-067-C	518-067-C	 00810060550724
518-039-C	518-039-C	 00810060550700
518-018-C	518-018-C	 00810060550625
518-021-C	518-021-C	 00810060550656
518-019-C	518-019-C	 00810060550632
518-022-C	518-022-C	 00810060550663
518-020-C	518-020-C	 00810060550649

Model #	Configuration	GTIN Portion of the GS1 POUCH Barcode
518-023-C	518-023-C	 00810060550670
<b>LLD Accessories (-C Models)</b>		
518-024-C	518-024-C	 00810060550687
518-027-C	518-027-C	 00810060550694
<b>VisiSheath™ Dilator Sheath</b>		
501-012	VisiSheath S 43cm	 00813132020354
501-014	VisiSheath M 43cm	 00813132020361
501-016	VisiSheath L 43cm	 00813132020378
501-112	VisiSheath S 33cm	 00813132020385
501-114	VisiSheath M 33cm	 00813132020392

Model #	Configuration	GTIN Portion of the GS1 POUCH Barcode
501-116	VisiSheath L 33cm	 00813132020408
501-212	VisiSheath S 23cm	 00813132020415
501-214	VisiSheath M 23cm	 00813132020422
501-216	VisiSheath L 23cm	 00813132020439
<b>TorqMax™ Sheath Grip Accessory</b>		
502-001	TorqMax Grip Accessory	 00813132020446
502-001	TorqMax Grip Accessory (Insert Label Barcode)	 00813132020446
<b>TightRail™ Rotating Dilator Sheath</b>		
545-509	TightRail 9F	 00813132020873
545-511	TightRail 11F	 00813132020880

Model #	Configuration	GTIN Portion of the GS1 POUCH Barcode
545-513	TightRail 13F	 00813132020897
<b>TightRail™ Sub-C</b>		
560-009	TightRail 9F	 00813132022785
560-011	TightRail 11F	 00813132022792
560-013	TightRail 13F	 00813132022808
<b>TightRail Mini™ Rotating Dilator Sheath</b>		
540-009	TightRail Mini 9F	 00813132020859
540-011	TightRail Mini 11F	 00813132020866
<b>SightRail™ Manual Dilator Sheath</b>		
550-008	Dilator Sheath 8.5F	 00813132020903
550-010	Dilator Sheath 10F	 00813132020910

Model #	Configuration	GTIN Portion of the GS1 POUCH Barcode
550-011	Dilator Sheath 11.5F	 00813132020927
550-013	Dilator Sheath 13F	 00813132020934
555-508	Long Dilator Sheath 8.5F	 00813132020941
555-510	Long Dilator Sheath 10F	 00813132020958
555-511	Long Dilator Sheath 11.5F	 00813132020965
555-513	Long Dilator Sheath 13F	 00813132020972