

PHILIPS

RESPIRONICS

SimplyGo Mini

Portable oxygen
concentrator

Be simply confident
on the go

SimplyGo Mini is a reliable, powerful choice in lightweight portable oxygen concentrators

Are you living your life? Taking trips to the store, museum or a musical? Spending weekends away with friends and family? If you use portable oxygen, SimplyGo Mini is a powerful choice for a reliable and lightweight portable oxygen concentrator that can help you do what you want, when you want.

Reliability you can count on

We built SimplyGo Mini to be tough. During development, SimplyGo Mini was rigorously tested in real world environments including impacts, vibrations, and temperatures to deliver reliable performance day in and day out.

Simply Portable

From a walk in the park to an overnight trip, SimplyGo Mini may help you remain active and provide you with flexibility while on the go.

The manufacturer of this POC has determined this device conforms to all applicable US Federal Aviation Administration (FAA) requirements for POC carriage and use on board aircraft.



Complies with
RTCA/DO-160G
section 21, category M



You have the **power**

SimplyGo Mini offers easy-to-change lithium ion battery options with standard AC and DC power cords that add convenience and flexibility.

Choose the power option that works best for your lifestyle needs

Available battery options:

- Standard battery for battery duration up to 4.5 hours* (pulse setting of 2 at 20 BPM). Resulting in total SimplyGo Mini weight of 2.3 kg (5 lbs)
- Extended battery for battery duration up to 9 hours* (pulse setting of 2 at 20 BPM). Resulting in total SimplyGo Mini weight of 2.7 kg (6 lbs)

*Battery duration times are based on new, fully charged batteries used with a new SimplyGo Mini system. Battery duration times will degrade with battery age, environmental operating conditions, use over time, and operational condition of the concentrator.



SimplyGo Mini at a glance


Easy-to-use
and read LCD
touch screen

Five pulse dose
settings designed
to keep you
on the go

Lightweight and
easy to carry

Rugged,
impact-resistant
design

Simple-to-
replace/remove
rechargeable
battery



Letting you be you

With its sleek, lightweight, and non-clinical design, SimplyGo Mini doesn't draw attention to itself. A stylish case is included for a comfortable way for you to carry by hand or over your shoulder.

Want to go hands-free? With an optional backpack specifically designed for SimplyGo Mini, you have a comfortable way to do it.





© 2016 Koninklijke Philips N.V. All rights reserved. Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

www.philips.com/respironics

PN449743