

# Lumify Android Device Compatibility List

## Supported Devices\*

Items in **RED** are new to this document revision

Items in **Blue** are compatible, but no longer available

Device Name	Screen Size	Connectivity Options	Internal Storage	Approximate Continuous Scan Time	USB Type	Supports Reacts	Supported OS
<b>Getac ZX10</b>	<b>10 in</b>	<b>Wi-Fi / LTE</b>	<b>128 GB</b>	<b>4hrs 39mins</b>	<b>C</b>	<b>Yes</b>	<b>12</b>
Samsung Galaxy S24 Ultra	6.8 in	Wi-Fi / 5G	256 GB	4 h 51 min	C	Yes	14
Samsung Galaxy S24+	6.7 in	Wi-Fi / 5G	256 GB	6 h	C	Yes	14
Samsung Galaxy S24	6.1 in	Wi-Fi / 5G	256 GB	4 h 15 min	C	Yes	14
Samsung Galaxy S23	6.1 in	Wi-Fi / 5G	256 GB	4 h 41 min	C	Yes	14
Samsung Galaxy S9 Tab	11.0 in	Wi-Fi / 5G	128 GB	5 h 55 min	C	Yes	13
Samsung Galaxy S9+ Tab	12.4 in	Wi-Fi / 5G	256 GB	7 h 20 min	C	Yes	13
Samsung Galaxy S9 Ultra Tab	14.6 in	Wi-Fi / 5G	256 GB	7 h 40 min	C	Yes	13
Samsung Galaxy A9+	11.0 in	Wi-Fi / 5G	128 GB	5 h 5 min	C	Yes	13
Samsung Galaxy Tab Active 5	8 in	Wi-Fi / 5G	128 GB	4 h 15 min	C	Yes	14
Samsung Galaxy Tab Active4 Pro (SM-T630)	10.1	Wi-Fi / 5G	64GB	4 h 16 min	C	Yes	Out of the box
Samsung Galaxy Tab Active 3	8 in	Wi-Fi / 5G	64 GB	4 h 14 min	C	Yes	Out of the box
Samsung Galaxy Tab A8	10.5 in	Wi-Fi / 5G	32 GB	3 h 57 min	C	Yes	Out of the box
Samsung Galaxy Tab S7 FE	12.4 in	Wi-Fi / 5G	64 GB	5 h 54 min	C	Yes	Out of the box
Samsung Galaxy Z Flip 3 5G	6.7 in	Wi-Fi / 5G	256 GB	2 h 51 min	C	Yes	Out of the box
Samsung Galaxy Z Fold 3 5G	7.6 in & 6.2 in	Wi-Fi / 5G	256 GB	3 h 44 min	C	Yes	Out of the box
Samsung Galaxy Tab A7 (SM-T505)	10.4 in	Wi-Fi	32 GB	3 h 58 min	C	Yes	Out of the box
Samsung A7 Lite	8.7 in	Wi-Fi	32 GB	3 h 20 min	C	Yes	Out of the box
Samsung Tab A 10.1 T 515	10.1 in	Wi-Fi	32 GB	4 h 40 min	C	Yes	Out of the box



**Philips Ultrasound LLC**

22100 Bothell Everett Hwy

Bothell, WA 98021-8431, USA

[www.healthcare.philips.com/ultrasound](http://www.healthcare.philips.com/ultrasound)

© Koninklijke Philips N.V. 2024. All rights reserved. Printed in USA.

Reproduction or transmission in whole or in part, in any form or by any means, electronic, mechanical, or otherwise, is prohibited without the prior written consent of the copyright owner.  
239034 Rev AL November, 2024

## Lumify Android Device Compatibility List

Huawei Mate 30 Pro 5G	6.53 in	Wi-Fi / 5G	256 GB	3 h 57 min	C	Yes	Out of the box
Huawei Changxiang 2	10.1 in	Wi-Fi	128 GB	3 h 31 min	C	Yes	Out of the box
Panasonic Toughbook (FZ-S1)	7 in	Wi-Fi / LTE	64 GB	2 h 16 min	C	Yes	Out of the box
Panasonic Toughbook (FZ-A3)	10.1 in	Wi-Fi / LTE	128 GB	2 h 19 min	C	Yes	Out of the box
Samsung Galaxy Tab S7 (SM-T870)	11 in	Wi-Fi	128 GB	5 h 55 min	C	Yes	Out of the box
Samsung Galaxy Tab S7+ (SM-T970)	12.4 in	Wi-Fi	128 GB	5 h 23 min	C	Yes	Out of the box
Samsung Galaxy Tab S7 5G (SM-T878U)	11 in	Wi-Fi / 5G	128 GB	5 h 44 min	C	Yes	Out of the box
Samsung Galaxy Tab S7+ 5G (SM-T978U)	12.4 in	Wi-Fi / 5G	128 GB	6 h 15 min	C	Yes	Out of the box
Samsung Galaxy Tab S6 (Model SM-T860, SM-T865)	10.5 in	Wi-Fi / LTE	128 GB	6 h 22 min	C	Yes	Out of the box
Samsung Galaxy S20 Ultra (SM-G988U)	6.9 in	Wi-Fi / LTE	128 GB	4 h 41 min	C	Yes	Out of the box
Samsung Galaxy S20+ (SM-G986U)	6.7 in	Wi-Fi / LTE	128 GB	5 h 10 min	C	Yes	Out of the box
Samsung Galaxy S20 (SM-G981U)	6.2 in	Wi-Fi / LTE	128 GB	5 h 23 min	C	Yes	Out of the box
Concept FPZ10	10 in	Wi-Fi/ LTE	64 GB	8 hrs +	C	Yes	Out of the box
Concept FPZ08	8 in	Wi-Fi	64 GB	4 h 27 min	C	Yes	Out of the box
Samsung Galaxy Tab S2 8" (Only Model SM-T713)	8 in	Wi-Fi	32 GB	2 h 30 min	B	Yes	Out of the box
<b>Samsung XCover Pro</b>	6.3 in	Wi-Fi / LTE	64 GB	2 h 37 min	C	Yes	Out of the box
<b>Samsung Z-Flip (SM-F700U)</b>	6.7 in	Wi-Fi / LTE	256 GB	4 h 2 min	C	Yes	Out of the box
<b>Samsung Galaxy S10e (SM-G970U, SM-G970F)</b>	5.8 in	Wi-Fi / LTE	128 GB	4 h 7 min	C	Yes	Out of the box
<b>Samsung Galaxy S10 (SM-G973U, SM-G973F)</b>	6.1 in	Wi-Fi / LTE	128 GB	4 h 3 min	C	Yes	Out of the box
<b>Samsung Galaxy S10+ (SM-G975U)</b>	6.4 in	Wi-Fi / LTE	128 GB	3 h 34 min	C	Yes	Out of the box
<b>Samsung Galaxy Note20 Ultra 5G (SM-N986U)</b>	6.7 in	Wi-Fi/LTE/5G	128 GB	3 h 40 min	C	Yes	Out of the box
<b>Samsung Galaxy Note10+ (SM-N975U, SM-N976V)</b>	6.8 in	Wi-Fi / LTE / 5G	256 GB	4 h 52 min	C	Yes	Out of the box



## Lumify Android Device Compatibility List

<b>Samsung Galaxy Note10 (SM-970U)</b>	6.3 in	Wi-Fi / LTE	256 GB	TBD	C	Yes	Out of the box
<b>Samsung Galaxy Note9 (SM-N960U)</b>	6.25 in	Wi-Fi / LTE	128 GB	3 h 26 min	C	Yes	Out of the box
<b>Google Pixel 4A</b>	5.7 in	Wi-Fi / LTE	64 GB	2 h 56 min	C	Yes	Out of the box
<b>Google Pixel 4</b>	5.7 in	Wi-Fi / LTE	64 GB	3 h 5 min	C	Yes	Out of the box
<b>Lenovo YT-X705F</b>	10.1 in	Wi-Fi	64 GB	4 h 21 min	C	Yes	Out of the box

\*Disclaimer:

Devices not included in the preceding table have not undergone compatibility testing with Lumify. Philips cannot guarantee Lumify system performance with non-compatible devices as compatibility testing with all devices on the market is not feasible. The devices have been evaluated using the operating system (OS) versions they were originally equipped with.

Subsequent OS releases, known as maintenance releases, that are higher than the version installed at the time of device launch has not been explicitly verified. Conducting explicit tests for newer major OS releases on older devices, which were initially equipped with a previous OS version, is not feasible. Lumify system functionality and performance may be impacted if significant operating system changes have been applied



# Lumify Android Device Compatibility List



**Philips Ultrasound LLC**  
22100 Bothell Everett Hwy  
Bothell, WA 98021-8431, USA  
[www.healthcare.philips.com/ultrasound](http://www.healthcare.philips.com/ultrasound)

© Koninklijke Philips N.V. 2024. All rights reserved. Printed in USA.  
Reproduction or transmission in whole or in part, in any form or by any means, electronic, mechanical, or otherwise, is prohibited without the prior written consent of the copyright owner.  
239034 Rev AL November, 2024