

Lumify iOS Device Compatibility List

Supported Devices*

Items in **RED** are new to this document revision
 Items in **Blue** are compatible, but no longer available

| Device Name | Year Released | Screen Size (in.) | Connectivity Options | Internal Storage (in GB) | Philips Case Available | LPM Required | Connector Type |
|----------------------------|---------------|-------------------|----------------------|--------------------------|------------------------|--------------|----------------|
| iPad Air 11-inch M2 | 2024 | 11 | Wi-Fi / LTE | 128 | No | Yes | Type C |
| iPad Air 13-inch M2 | 2024 | 13 | Wi-Fi / LTE | 128 | No | Yes | Type C |
| iPad Pro 12.9-inch 6th Gen | 2022 | 12.9 | Wi-Fi / LTE | 128 | No | Yes | Type C |
| iPad Pro 11- inch 4th Gen | 2020 | 11 | Wi-Fi / LTE | 128 | No | Yes | Type C |
| iPad Pro 11- inch (M4) | 2024 | 11 | Wi-Fi / LTE | 256 | No | Yes | Type C |
| iPad Pro 12.9-inch 5th Gen | 2021 | 12.9 | Wi-Fi / LTE | 128 | No | Yes | Type C |
| iPad Air 5th Gen M1 | 2022 | 10.9 | Wi-Fi / LTE | 64 | No | Yes | Type C |
| iPad Pro 13- inch M4 | 2024 | 13 | Wi-Fi / LTE | 512 | No | Yes | Type C |
| iPhone 15 | 2023 | 6.12 | Wi-Fi / LTE | 128 | No | Yes | Type C |
| iPhone 15 Plus | 2023 | 6.69 | Wi-Fi / LTE | 128 | No | Yes | Type C |
| iPhone 15 Pro | 2023 | 6.12 | Wi-Fi / LTE | 128 | No | Yes | Type C |
| iPhone 14 | 2022 | 6.06 | Wi-Fi / LTE | 128 | No | Yes | Lightning |
| iPhone 14 Pro | 2022 | 6.1 | Wi-Fi / LTE | 128 | No | Yes | Lightning |
| iPad 10th Gen | 2022 | 10.9 | Wi-Fi / LTE | 64/256 | No | Yes | Type C |
| iPad Pro 6th Gen | 2022 | 12.9 | Wi-Fi / LTE | 128/256/512/1TB/2TB | No | Yes | Type C |
| iPad Air 5th Gen | 2022 | 10.9 | Wi-Fi / LTE | 64/256 | No | Yes | Type C |



Lumify iOS Device Compatibility List

| | | | | | | | |
|----------------------------|------|------|-------------|-----------------|-----|-----|-----------|
| iPad Mini 6th Gen | 2021 | 8.3 | Wi-Fi / LTE | 64/256 | No | Yes | Type C |
| iPhone 13 mini | 2021 | 5.4 | Wi-Fi / LTE | 128/256/512 | No | Yes | Lightning |
| iPhone 13 | 2021 | 6.1 | Wi-Fi / LTE | 128/256/512 | No | Yes | Lightning |
| iPhone 13 Pro | 2021 | 6.1 | Wi-Fi / LTE | 128/256/512/1TB | No | Yes | Lightning |
| iPhone 13 Pro Max | 2021 | 6.7 | Wi-Fi / LTE | 128/256/512/1TB | No | Yes | Lightning |
| iPad (9th gen) | 2021 | 10.2 | Wi-Fi / LTE | 64/256 | Yes | Yes | Lightning |
| iPhone 12 Pro Max | 2020 | 6.7 | Wi-Fi / LTE | 128/256/512 | No | Yes | Lightning |
| iPhone 12 Mini | 2020 | 5.4 | Wi-Fi / LTE | 64/128/256 | No | Yes | Lightning |
| iPhone 12 | 2020 | 6.1 | Wi-Fi / LTE | 64/128/256 | No | Yes | Lightning |
| iPhone 12 Pro | 2020 | 6.1 | Wi-Fi / LTE | 128/256/512 | No | Yes | Lightning |
| iPad (8 th gen) | 2020 | 10.2 | Wi-Fi / LTE | 32/128 | Yes | Yes | Lightning |
| iPhone 12 Pro | 2020 | 6.1 | Wi-Fi / LTE | 128/256/512 | No | Yes | Lightning |
| iPhone 12 | 2020 | 6.1 | Wi-Fi / LTE | 64/128/256 | No | Yes | Lightning |
| iPhone 11 Pro Max | 2019 | 6.5 | Wi-Fi / LTE | 64/256/512 | No | Yes | Lightning |
| iPhone 11 Pro | 2019 | 5.8 | Wi-Fi / LTE | 64/256/512 | No | Yes | Lightning |
| iPhone 11 | 2019 | 6.1 | Wi-Fi / LTE | 64/128/256 | Yes | Yes | Lightning |
| iPhone XS Max | 2018 | 6.5 | Wi-Fi / LTE | 64/256/512 | No | Yes | Lightning |
| iPhone XS | 2018 | 5.8 | Wi-Fi / LTE | 64/256/512 | Yes | Yes | Lightning |
| iPhone XR | 2018 | 6.1 | Wi-Fi / LTE | 64/256/512 | Yes | Yes | Lightning |
| iPhone X | 2017 | 5.8 | Wi-Fi / LTE | 64/256 | Yes | Yes | Lightning |
| iPhone 8 Plus | 2017 | 5.5 | Wi-Fi / LTE | 64/256 | No | Yes | Lightning |
| iPhone 8 | 2017 | 4.7 | Wi-Fi / LTE | 64/256 | Yes | Yes | Lightning |
| iPhone 7 Plus | 2016 | 5.5 | Wi-Fi / LTE | 32/128/256 | No | Yes | Lightning |



Lumify iOS Device Compatibility List

| | | | | | | | |
|--------------------------------|------|------|-------------|------------|-----|-----|-----------|
| iPhone 7 | 2016 | 4.7 | Wi-Fi / LTE | 32/128/256 | Yes | Yes | Lightning |
| iPad (7 th gen) | 2019 | 10.2 | Wi-Fi / LTE | 32/128 | Yes | Yes | Lightning |
| iPad Air (3 rd gen) | 2019 | 10.5 | Wi-Fi / LTE | 64 / 256 | No | Yes | Lightning |
| iPad Mini 5 | 2019 | 7.9 | Wi-Fi / LTE | 64/256 | Yes | Yes | Lightning |
| iPad (6 th gen) | 2018 | 9.7 | Wi-Fi / LTE | 32/128 | Yes | Yes | Lightning |
| iPad (6 th gen) | 2017 | 9.7 | Wi-Fi / LTE | 32/128 | Yes | Yes | Lightning |
| iPad Pro (2 nd gen) | 2017 | 12.9 | Wi-Fi / LTE | 64/256/512 | No | Yes | Lightning |
| iPad Pro (2 nd gen) | 2017 | 10.5 | Wi-Fi / LTE | 64/256/512 | No | Yes | Lightning |



Lumify iOS Device Compatibility List

***Disclaimer:**

Devices not included in the preceding table have not undergone compatibility testing with Lumify. Philips cannot guarantee Lumify system performance with non-compatible devices as compatibility testing with all devices on the market is not feasible. The devices have been evaluated using the operating system (OS) versions they were originally equipped with.

Subsequent OS releases, known as maintenance releases, that are higher than the version installed at the time of device launch has not been explicitly verified. Conducting explicit tests for newer major OS releases on older devices, which were initially equipped with a previous OS version, is not feasible. Lumify system functionality and performance may be impacted if significant operating system changes have been applied.

Please note that Reacts tele-ultrasound solution is not available on Lumify when using the above iOS devices. All iOS devices currently compatible require use of Lumify Power Module and utilize USB-C transducer cables and will all continuously scan for 2 to 2.5 hrs.

Devices on the supported list satisfy all the following criteria: 1. Passed software testing on the latest Lumify software, 2. Passed Philips hardware testing, 3. Is readily available from the device manufacturer.

Philips cannot test every device available on the market, so this list represents a range of devices to satisfy a variety of user needs with respect to form factor, price, battery life and connectivity options.

