

# PHILIPS

## Ultrasound Systems DICOM SR template support

	<b>OB/Gyn (TID 5000)</b>	<b>Vascular (TID 5100)</b>	<b>Adult Echo (TID 5200)</b>	<b>Pediatric (TID 5220)</b>	<b>Fetal Echo (TID 5220)</b>	<b>Small Parts (TID 995400)</b>
<b>ClearVue</b>	X	X	X			
<b>CX50/CX30</b>	X	X	X	X		
<b>HD11</b>	X	via serial only	X			
<b>HD15 R1.5.1</b>	X	via serial only	X			
<b>HD15 R3</b>	X	X	X	X		
<b>HD6</b>	X		X			
<b>HD7</b>	X		X			
<b>iE33</b>		X	X	X		
<b>iU22</b>	X	X	X			
<b>Lumify</b>	No SR export	No SR export	No SR export	No SR export	No SR export	No SR export
<b>QLAB</b>			X <sup>1</sup>			
<b>Sparq</b>	X	X	X			
<b>Xperius</b>	No SR export	No SR export	No SR export	No SR export	No SR export	No SR export

