

Lumify Android Device Compatibility List

Supported Devices*

Items in **RED** are new to this document revision
 Items in **Blue** are compatible, but no longer available

Device Name	Screen Size	Connectivity Options	Internal Storage	Approximate Continuous Scan Time	USB Type	Supports Reacts
Samsung Galaxy Tab Active4 Pro (SM-T630)	10.1	Wi-Fi / 5G	64GB	4 h 16 min	C	Yes
Lenovo Tab K10	10.3 in	Wi-Fi	64 GB	6 h 44 min	C	Yes
Samsung Galaxy Tab Active 3	8 in	Wi-Fi / 5G	64 GB	4 h 14 min	C	Yes
Samsung Galaxy Tab A8	10.5 in	Wi-Fi / 5G	32 GB	3 h 57 min	C	Yes
Samsung Galaxy Tab S7 FE	12.4 in	Wi-Fi / 5G	64 GB	5 h 54 min	C	Yes
Samsung Galaxy Z Flip 3 5G	6.7 in	Wi-Fi / 5G	256 GB	2 h 51 min	C	Yes
Samsung Galaxy Z Fold 3 5G	7.6 in & 6.2 in	Wi-Fi / 5G	256 GB	3 h 44 min	C	Yes
Samsung Galaxy Tab A7 (SM-T505)	10.4 in	Wi-Fi	32 GB	3 h 58 min	C	Yes
Samsung A7 Lite	8.7 in	Wi-Fi	32 GB	3 h 20 min	C	Yes
Samsung Tab A 10.1 T 515	10.1 in	Wi-Fi	32 GB	4 h 40 min	C	Yes
Huawei Mate 30 Pro 5G	6.53 in	Wi-Fi / 5G	256 GB	3 h 57 min	C	Yes
Huawei Changxiang 2	10.1 in	Wi-Fi	128 GB	3 h 31 min	C	Yes
Panasonic Toughbook (FZ-S1)	7 in	Wi-Fi / LTE	64 GB	2 h 16 min	C	Yes
Panasonic Toughbook (FZ-A3)	10.1 in	Wi-Fi / LTE	128 GB	2 h 19 min	C	Yes
Samsung Galaxy Tab S7 (SM-T870)	11 in	Wi-Fi	128 GB	5 h 55 min	C	Yes
Samsung Galaxy Tab S7+ (SM-T970)	12.4 in	Wi-Fi	128 GB	5 h 23 min	C	Yes
Samsung Galaxy Tab S7 5G (SM-T878U)	11 in	Wi-Fi / 5G	128 GB	5 h 44 min	C	Yes
Samsung Galaxy Tab S7+ 5G (SM-T978U)	12.4 in	Wi-Fi / 5G	128 GB	6 h 15 min	C	Yes
Samsung Galaxy Tab S6 (Model SM-T860, SM-T865)	10.5 in	Wi-Fi / LTE	128 GB	6 h 22 min	C	Yes



Lumify Android Device Compatibility List

Items in **RED** are new to this document revision

Items in **Blue** are compatible, but no longer available

Device Name	Screen Size	Connectivity Options	Internal Storage	Approximate Continuous Scan Time	USB Type	Supports Reacts
Samsung Galaxy S20 Ultra (SM-G988U)	6.9 in	Wi-Fi / LTE	128 GB	4 h 41 min	C	Yes
Samsung Galaxy S20+ (SM-G986U)	6.7 in	Wi-Fi / LTE	128 GB	5 h 10 min	C	Yes
Samsung Galaxy S20 (SM-G981U)	6.2 in	Wi-Fi / LTE	128 GB	5 h 23 min	C	Yes
Concept FPZ10	10 in	Wi-Fi/ LTE	64 GB	8 hrs +	C	Yes
Concept FPZ08	8 in	Wi-Fi	64 GB	4 h 27 min	C	Yes
Samsung Galaxy Tab S2 8" (Only Model SM-T713)	8 in	Wi-Fi	32 GB	2 h 30 min	B	Yes
Samsung XCover Pro	6.3 in	Wi-Fi / LTE	64 GB	2 h 37 min	C	Yes
Samsung Z-Flip (SM-F700U)	6.7 in	Wi-Fi / LTE	256 GB	4 h 2 min	C	Yes
Samsung Galaxy S10e (SM-G970U, SM-G970F)	5.8 in	Wi-Fi / LTE	128 GB	4 h 7 min	C	Yes
Samsung Galaxy S10 (SM-G973U, SM-G973F)	6.1 in	Wi-Fi / LTE	128 GB	4 h 3 min	C	Yes
Samsung Galaxy S10+ (SM-G975U)	6.4 in	Wi-Fi / LTE	128 GB	3 h 34 min	C	Yes
Samsung Galaxy Note20 Ultra 5G (SM-N986U)	6.7 in	Wi-Fi/LTE/5G	128 GB	3 h 40 min	C	Yes
Samsung Galaxy Note10+ (SM-N975U, SM-N976V)	6.8 in	Wi-Fi / LTE / 5G	256 GB	4 h 52 min	C	Yes
Samsung Galaxy Note10 (SM-970U)	6.3 in	Wi-Fi / LTE	256 GB	TBD	C	Yes
Samsung Galaxy Note9 (SM-N960U)	6.25 in	Wi-Fi / LTE	128 GB	3 h 26 min	C	Yes
Google Pixel 4A	5.7 in	Wi-Fi / LTE	64 GB	2 h 56 min	C	Yes
Google Pixel 4	5.7 in	Wi-Fi / LTE	64 GB	3 h 5 min	C	Yes
Lenovo TB-8705F	8.0 in	Wi-Fi	32 GB	4 h 52 min	B	No
Lenovo YT-X705F	10.1 in	Wi-Fi	64 GB	4 h 21 min	C	Yes



Philips Ultrasound, Inc.

22100 Bothell Everett Hwy

Bothell, WA 98021-8431, USA

www.healthcare.philips.com/ultrasound

© Koninklijke Philips N.V. 2017. All rights reserved. Printed in USA.

Reproduction or transmission in whole or in part, in any form or by any means, electronic, mechanical, or otherwise, is prohibited without the prior written consent of the copyright owner.

239034 Rev AL April, 2023

Lumify Android Device Compatibility List

*Devices on the supported list satisfy all of the following criteria: 1. Passed software testing on the latest Lumify software, 2. Passed Philips hardware testing, 3. Device is readily available from the device manufacturer. Philips cannot test every device available on the market, so this list represents a range of devices to satisfy a variety of user needs with respect to form factor, price, battery life and connectivity options.



Philips Ultrasound, Inc.

22100 Bothell Everett Hwy

Bothell, WA 98021-8431, USA

www.healthcare.philips.com/ultrasound

© Koninklijke Philips N.V. 2017. All rights reserved. Printed in USA.

Reproduction or transmission in whole or in part, in any form or by any means, electronic, mechanical, or otherwise, is prohibited without the prior written consent of the copyright owner.

239034 Rev AL April, 2023

Lumify Android Device Compatibility List

No Longer Supported**

Items in **BLUE** have been removed in this document revision

Device Name	Screen Size	Connectivity Options	Internal Storage	Approximate Continuous Scan Time	USB Type
Samsung Galaxy S8 Ultra	14.6 in	Wi-Fi / 5G	128 GB	5 h 47 min	C
Samsung Galaxy S22 Ultra	6.8 in	Wi-Fi / 5G	128 GB	4 h 48 min	C
Samsung Galaxy Note20 5G (SM-N981U)	6.7 in	Wi-Fi/LTE/5G	128 GB	3 h 45 min	C
Samsung Galaxy Tab S4 (Model SM-T830, SM-T835, SM-T837)	10.5 in	Wi-Fi or Wi-Fi/LTE	64 GB	5 h 30 min	C
Samsung Galaxy Tab A- 2018 (Model SM-T387)	8.0 in	Wi-Fi /LTE	32 GB	3 h 29 min	B
Samsung Galaxy S9 (SM-G960U)	5.8 in	Wi-Fi / LTE	32 GB	3 h 45 min	C
Samsung Galaxy S9+ (SM-G965U)	6.2 in	Wi-Fi / LTE	64 GB	4 h	C
Samsung Galaxy Note8 (SM-N950U)	6.3 in	Wi-Fi / LTE	64 GB	3 h 30 min	C
Samsung Galaxy S8 (SM-G950U)	5.8 in	Wi-Fi / LTE	32 GB	3 h 30 min	C



Philips Ultrasound, Inc.

22100 Bothell Everett Hwy

Bothell, WA 98021-8431, USA

www.healthcare.philips.com/ultrasound

© Koninklijke Philips N.V. 2017. All rights reserved. Printed in USA.

Reproduction or transmission in whole or in part, in any form or by any means, electronic, mechanical, or otherwise, is prohibited without the prior written consent of the copyright owner.

239034 Rev AL April, 2023

Lumify Android Device Compatibility List

Samsung Galaxy S8+ (SM-G955U)	6.2 in	Wi-Fi / LTE	64 GB	4 h	C
Samsung Galaxy Tab A – 2019 (Model SM-T510)	10.1 in	Wi-Fi	64 GB	4 h 37 min	C
Samsung Galaxy Tab A w/ S-Pen (SM-P200, SM-P205)	8.0 in	Wi-Fi / LTE	32 GB	3 h 16 min	C
Samsung Galaxy Tab A -2018 (Model SM-T590, SM-T595)	10.5 in	Wi-Fi /LTE	32 GB	3 h 38 min	C
Huawei Mediapad M5 (Model SHT-W09, SHT-AL09) *removed for reason other than compatibility	8.4 in	Wi-Fi or Wi-Fi/LTE	64 GB	3 h 30 min	C
Huawei Mediapad M3 (Model BTV-W09, BTV-DL09) *removed for reason other than compatibility	8.4 in	Wi-Fi or Wi-Fi/LTE	32 GB	3 h 30 min	B
Samsung Galaxy S7 Edge (SM-G935)	5.5 in	Wi-Fi/LTE	32 GB	2 h 15 min	B
Samsung J7 (Model SM-J710)	5.5 in	Wi-Fi/LTE	16 GB	4 h 15 min	B
Samsung Galaxy Tab S3 (Model SM-T820,SM-T825)	9.7 in	Wi-Fi/LTE	32 GB	4 h 30 min	C
Samsung Galaxy View	18.4 in	Wi-Fi or Wi-Fi/LTE	32 GB	5 h	B
Google Nexus 6 (not Nexus 6P)	6 in	Wi-Fi/Cellular	32/64 GB	1 h 45 min	B
Google Pixel 3a / Pixel 3a XL	5.6 in/ 6.0 in	Wi-Fi / LTE	64 GB	3 h 26 min	C
Google Pixel 3/ 3 XL	5.5 / 6.3 in	Wi-Fi/LTE	64 GB	3 h	C
Google Pixel 2 / 2 XL	5.5 /6.3 in	Wi-Fi/LTE	64 GB	3 h	C
NVIDIA SHIELD Tablet (V1)	8 in	Wi-Fi	16 GB	1 h 30 min	B
Samsung Galaxy Tab Active	8 in	Wi-Fi	16 GB	3 h 30 min	B
Samsung Galaxy Tab S	8.5 in	Wi-Fi	16 GB	2 h 15 min	B
Nvidia Shield Tablet K1	8 in	Wi-Fi	16 GB	2 h 30 min	B



Philips Ultrasound, Inc.

22100 Bothell Everett Hwy

Bothell, WA 98021-8431, USA

www.healthcare.philips.com/ultrasound

© Koninklijke Philips N.V. 2017. All rights reserved. Printed in USA.

Reproduction or transmission in whole or in part, in any form or by any means, electronic, mechanical, or otherwise, is prohibited without the prior written consent of the copyright owner.
239034 Rev AL April, 2023

Lumify Android Device Compatibility List

Samsung Galaxy Tab S2 8" (Model SM-T710) – Not Compatible with S4-1	8 in	Wi-Fi	32 GB	2 h 30 min	B
Samsung Galaxy Tab S2 9.7" (all models)	9.7 in	Wi-Fi	32 GB	3 h	B
Nokia 6 (Model TA-1003)	5.5 in	Wi-Fi/LTE	32 GB	2 h 30 min	B
Samsung Galaxy Note5	5.7 in	Wi-Fi/Cellular	32/64/128 GB	3 h 15 min	B
Samsung Galaxy Tab S	8.5 in	Wi-Fi	16 GB	2 h 15 min	B
Google Nexus 9	9 in	Wi-Fi/LTE	16/32 GB	2 h	B
Lenovo Tab4 8	8.0 in	Wi-Fi / LTE	16 GB	4 h 26 min	B
Lenovo Smart Tab P10 (Model TB-X705F)	10.1 in	Wi-Fi	32 GB	4 h 53 min	C
Lenovo Smart Tab M10 (Model TB-X605F)	10.1 in	Wi-Fi	32 GB	3 h 43 min	C
Lenovo Yoga Tab 3 (Model YT3-850F, P/N#ZA090094US)	8 in	Wi-Fi	16 GB	4 h	B
Lenovo Yoga Tab 3 (Model YT3-X50F P/N#ZA0H0022US, ZA0H0064US)	10 in	Wi-Fi	16 GB	4 h	B
Panasonic Toughpad FZ-A2	10.1 in	Wi-Fi/LTE	32 GB	3 h 30 min	C
Panasonic Toughpad FZ-B2 (USB C models only)	7 in	Wi-Fi/LTE	32 GB	3 h 30 min	C

**The devices on this list represent devices that at one point supported the Lumify system but no longer do for one or more of the following reasons: 1. Did not pass software testing on the latest Lumify software, 2. Did not pass Philips hardware testing, or 3. Is no longer available from the manufacturer. Lumify may still run on these devices, but due to the above criteria, Philips can no longer verify Lumify compatibility

Tablet Compatibility is subject to change, and Philips reserves the right to make changes without notice. Devices listed under “No Longer Supported” are not tested for use with the current version of Lumify by Philips.

