

PHILIPS

Ultrasound

Lumify



The ultrasound you trust.
On the devices you love.

Lumify ultrasound on Apple iOS

Fast diagnosis and treatment.

Our app-based ultrasound makes it easy to get started. Three simple steps, and you're ready to go.

1



Download
the Lumify app to
your compatible
smart device.*

2



Connect
the Lumify transducer
to your smart device.

3



Scan
with the exceptional
quality of Philips
imaging technology.

Continuously scan for more than 2.25 hours
without risk of overheating.



2 *See the Lumify website for the list of compatible smart device options:
philips.com/lumify-compatible-devices.

Exceptional imaging performance



Lumify
C5-2
broadband
curved array
transducer

High-resolution imaging for abdomen,
gallbladder, lung and vascular applications




Lumify
L12-4
broadband
linear array
transducer

High-resolution imaging for lung, musculoskeletal,
soft tissue, superficial and vascular applications




Lumify
S4-1
broadband
phased array
transducer


High-resolution imaging for cardiac, lung,
abdominal, FAST and Ob/Gyn applications



Thyroid nodule using L12-4 transducer




Carotid artery with color
using L12-4 transducer




Brachial plexus using L12-4 transducer




Pancreas using C5-2 transducer




IVC and M-mode using S4-1 transducer



Portal vein using C5-2 transducer



Jugular vein and carotid artery
with color using L12-4 transducer



Rotator cuff using L12-4 transducer

Keep the power where it belongs. With you and your device.

The power to turn an Apple device into handheld ultrasound

The Lumify power module provides the necessary bridge between an Apple device and an ultrasound transducer. Lumify offers the longest scan time among handheld ultrasound devices.

The Lumify swappable power module is available for Apple devices that utilize a lightning connector, enabling lightweight, high-quality Philips transducers to scan all day long on Apple iPhone and iPad.



Your data.
Your way.

With Lumify, your data and images stay with you and are stored however you choose. Lumify is HIPAA-compliant and is in use in hospitals worldwide.

Decades of imaging expertise. Quality you can trust.

Crystal-clear diagnostic insights depend on high-resolution imaging performance across shallow and deeper applications.



For a full list of compatible devices go to
philips.com/lumify-compatible-devices.

Accessible pricing

Because Lumify works on compatible Apple devices, it can relieve you of the burden of paying for a custom ultrasound machine.

When you purchase Lumify, you'll receive everything you'll need to get started.

- Lumify iOS power module and cables
- A Philips Lumify transducer
- Lumify ultrasound app
- Ongoing software upgrades
- Manufacturer's warranty on the hardware
- Access to Philips service and education

Call **1-844-MYLUMIFY**
or visit **philips.com/lumify**



