

Accessories and Supplies



Everything you need.

Airway adapters, sensors, masks and more.

Latex Free.

A wide range of innovative Accessories and Supplies from Respironics.

Accessories & Supplies – Simple. Convenient. Reliable.

As a recognized global leader in respiratory monitoring, Respironics has been providing innovative and cost effective solutions to our partners for over 20 years. With an entire family of patented products, including airway adapters, nasal cannulas, flow sensors and more, we offer a wide array of solutions to meet your patient's needs in any clinical environment. Respironics provides comprehensive technical, clinical, and marketing support to assure optimal utilization of our accessories and supplies.

PRODUCT FEATURES:

- Reusable & disposable versions available
- Unique combined CO₂/flow sensors
- Clinically reliable, cost effective solutions
- Proprietary sample cell protects internal NDIR components from contamination
- Easy to use for all patient populations
- Private label option

A complete monitoring solution



Accessories and supplies.
Just one more way it all comes together when you work with Respironics for OEM solutions.

LoFlo Sidestream CO₂ Sampling Kits

A complete line of CO₂ sampling kits and adapters for the LoFlo sidestream CO₂ sensor which samples at 50 mL/minute

- All in one, single patient use sample cell, filter and cannula
- Adult, pediatric and neonatal cannulas and adapters
- Cannulas available for nasal CO₂ sampling with O₂ delivery
- Cannulas available for high humidity environments

Flow Sensors

Flow sensors for proximal measurement of flow, pressure and volume on mechanically ventilated patients. Designed for use with the FloTrak Elite system.

- Single patient use
- Patented, automatic sensor identification
- No calibration required
- Adult, pediatric and neonatal sizes

Combined CO₂/Flow Sensors

Combined CO₂ /Flow sensors that allow simultaneous measurement of CO₂ flow, pressure and volume for volumetric capnography measurements. Includes VCO₂, airway deadspace, alveolar ventilation, Vd/Vt and more.

- Single patient use
- No calibration required
- Adult, pediatric and neonatal sizes



Mainstream CO₂ Airway Adapters

A variety of single patient use and reusable mainstream CO₂ airway adapters designed for use with the CAPNOSTAT sensor:

- Single patient use and reusable airway adapters available
- No calibration required
- Adult/pediatric and infant/neonatal sizes

CAPNO₂mask

Non-invasive mask that simultaneously delivers oxygen and measures mainstream end-tidal CO₂. Designed for use with the CAPNOSTAT sensor:

- Single patient use
- Delivers O₂ concentrations equivalent to that of a traditional O₂ mask
- Large-adult, adult and pediatric sizes

Removable Sample Cell and Filter

Patented disposable, external sample cell and filter are integrated into each sample line to protect internal NDIR components from contamination.

RESPIRONICS
Envisioning tomorrow. Improving today.

For more information, visit <http://oem.respironics.com> or call the OEM team at 1.800.243.3444, Option 3, or 203.697.6488

Specifications subject to change without notice. Customer is responsible for all regulatory approvals and market clearance. CAUTION: US law restricts this device to sale by or on the order of a licensed medical practitioner. Respironics, CAPNOSTAT, LoFlo, FloTrak Elite, CAPNO₂mask, and Envisioning tomorrow. Improving today, are registered trademarks belonging to Koninklijke Philips Electronics N.V.
©2010 Respironics, Inc. and its affiliates.