



# Inspiring innovation

Philips is committed to enhancing patients' lives by delivering cost-effective solutions that drive better awareness, diagnosis, treatment and monitoring.

Our Respiratory Drug Delivery product portfolio provides both patients and caregivers with a wide range of high quality products. These have been designed to deliver drugs effectively and enable patients with Asthma, Chronic Obstructive Pulmonary Disease (COPD), Cystic Fibrosis and other respiratory diseases to successfully manage their conditions.



# Be inspired

Philips Respironics solutions catalog

© 2017 Koninklijke Philips N.V. All rights reserved. Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. (Royal Philips) or their respective owners.

[www.philips.com/respironics](http://www.philips.com/respironics)  
JS 3/20/17 PN470050



RRP: \$19.95

## OptiChamber Diamond

Designed for ease of use to help maximize delivery of medication into the lungs.

- ▶ Anti-static material – suspends the aerosol longer, giving patients more time to take their medication
- ▶ Low resistance exhalation valve – allows user to breathe easily through the chamber and visually confirms exhalation
- ▶ Available with detachable LiteTouch Masks



RRP: \$15.00

### LiteTouch Masks with soft seal technology

Soft contoured cushion that gently rests on the face to minimise leakage. Available in 3 sizes:

- ▶ Small – 0 to 18 months
- ▶ Medium – 1 to 5 years
- ▶ Large – 5 years to adult

#### Part number OptiChamber Diamond

1079820	OptiChamber Diamond
1083785	LiteTouch small mask
1083784	LiteTouch medium mask
1083783	LiteTouch large mask



RRP: \$26.00

## PersonalBest

Enables user to monitor lung function. Can help patients keep their asthma under control when used with an asthma action plan.

- ▶ Three zone system for easy and objective way of monitoring asthma control
- ▶ Fold out handle keep the hand away from the airstream, aiding accurate peak flow readings
- ▶ Compliant with international ISO standard

#### Part number PersonalBest

1022973	PersonalBest EU Universal Range Peak Flow Meter
HS714-200	Plastic mouthpiece with one way valve (Pack of 200)



RRP: \$179.00

## InnoSpire Elegance

### Robust and powerful

Designed to cope with frequent use, the InnoSpire Elegance compressor is combined with the SideStream Reusable nebulizer and can be used continuously, at home or in a hospital environment.

#### Full system package includes

- ▶ Reusable SideStream with mouthpiece and tubing
- ▶ Adult (1) and Pediatric mask (1)
- ▶ Replacement air inlet filters (4)
- ▶ Unit carrying case (1)

#### Part number InnoSpire Elegance

1099970	InnoSpire Elegance - 230V 50Hz AU plug
1101452	InnoSpire Elegance carrying case
<b>Filters</b>	
1102088	Air inlet filters (All InnoSpire Essence, Elegance models)



RRP: \$290.00

## InnoSpire Deluxe

### Therapy you can trust

Designed to easily manage prolonged use in a clinical environment, the InnoSpire Deluxe is suitable for both professional use and frequent home use. The InnoSpire Deluxe is combined with either SideStream Plus or SideStream reusable.

#### Full system package includes

- ▶ Reusable SideStream Plus with mouthpiece and tubing
- ▶ Replacement air inlet filter (4)

#### Part number InnoSpire Deluxe

<b>Compressor nebulizer systems - including SideStream Plus kit</b>	
1110066	InnoSpire Deluxe – 230V/50Hz
<b>Filters</b>	
1113525	Air inlet filters (pack of 4) – (for 230V/50Hz Models)



- ▶ Durable lid design with 180° hinge



- ▶ Nebulizer handset holder in integral storage area



- ▶ Handle for easy transport



RRP: \$325.00

## InnoSpire Mini

### Lightweight and battery powered

InnoSpire Mini allows patients to get their therapy whenever they want and wherever they may be.

#### Full system package includes

- ▶ SideStream Plus with mouthpiece and tubing
- ▶ Rechargeable battery
- ▶ 12V DC adapter
- ▶ AC adapter
- ▶ Replacement air inlet filter (4)
- ▶ Unit carrying case

#### Part number InnoSpire Mini

1111452	InnoSpire Mini with SideStream Plus and battery
<b>Filters and accessories</b>	
1109488	Air inlet filters (Pack of 4)
1113183	AC power/ battery charger
1109455	12V DC (car) adapter
1110726	Rechargeable lithium ion battery

## SideStream Nebulizers

A venturi nebulizer that can be used to reliably and rapidly nebulize the majority of commonly prescribed bronchodilators, antibiotics or steroids



RRP: \$19.95

### SideStream Reusable – To be replaced after 12 months

#### Part number SideStream reusable

1100311	SideStream Patient Pack Contents include: Sidestream Resuable nebulizer, mouthpiece, tubing, adult and child mask, 4 X filters for Portaneb, Inspiration Elite and InnoSpire Elegance
---------	--

### SideStream Plus – Breath enhanced nebulizer –

To be replaced after 12 months

- ▶ Low resistance valve to deliver fast treatments with less waste

#### Part number SideStream Plus

1121451	SideStream Plus Patient Pack Contents include: Sidestream Plus nebulizer, mouthpiece, tubing, adult and child mask, 4 X filters for InnoSpire Deluxe and InnoSpire Mini
---------	--

#### Optional SideStream plus accessories:



- ▶ Filter attachment



- ▶ Aerosol hose

