


How to choose the right breast pump

When women are committed to breastfeeding, it's good to know which pump will suit her needs best



Quality pumps that satisfy the high expectations of moms and midwives¹



Proven clinically and in everyday use, Philips Avent's lightweight pumps use Natural Motion Technology to comfortably stimulate nipples by combining gentle massage with ideal suction levels. For even more ease-of-use and comfort, our soft, silicone cushion is designed to adapt to 99.98% of nipple sizes. And with sound levels are just 50dB – as quiet as a peaceful office or home – it's no wonder our pumps scored 97% in overall effectiveness with lactating moms.²



References:

1. Tested with 100 midwives in Germany, 05/2020, www.hebammen-testen.de
2. Based on product clinical trial results with 20 participants (Netherlands, 2019), average score for single and double electric breast pump; 95% of participants agree that our breast pumps are effective (single electric); 100% of participants agree that our breast pumps are effective (single electric).