

CT CROSS REFERENCE

26/11/20

Subject to change without notice. All products listed may be trademarked by the associated OEM.

Brand	CT Scanner Model	Dunlee Replacement Tube	Dunlee Part #
GE	BrightSpeed™ Edge	DA200ULTRA	9896-055-88704
GE	BrightSpeed™ Elite	DA200ULTRA	9896-055-88704
GE	BrightSpeed™ Excel	DA200ULTRA	9896-055-88704
GE	BrightSpeed Edge™ (Select / Lite)	DA165NP	9896-055-83901
GE	BrightSpeed Elite™ (Select / Lite)	DA165NP	9896-055-83901
GE	BrightSpeed Excel™ (Select / Lite)	DA165NP	9896-055-83901
GE	Brivo™ CT315/325	DA135 CT/e	9896-055-80800
GE	Brivo™ CT385	DA135 CT/e	9896-055-80800
GE	Discovery™ NM/CT570c	REEVO240G	9896-055-96401
GE	Discovery™ CT590 RT	REEVO240G	9896-055-96401
GE	Discovery™ NM/CT 670	DA200ULTRA	9896-055-88704
GE	Discovery™ NM/CT 670 Pro	DA200ULTRA	9896-055-88704
GE	Discovery™ PET/CT 610	DA200ULTRA	9896-055-88704
GE	Discovery™ PET/CT 710	DA200ULTRA	9896-055-88704
GE	Discovery™ RX	DA200ULTRA	9896-055-88704
GE	Discovery™ ST (PET/CT)	DA200ULTRA	9896-055-88704
GE	Discovery™ STE (PET/CT)	DA200ULTRA	9896-055-88704
GE	Discovery™ VCT	REEVO240G	9896-055-96401
GE	HiSpeed CT/e, CT Dual™	DA135 CT/e	9896-055-80800
GE	HiSpeed DX/i	DA165NP	9896-055-83901
GE	HiSpeed FX/i	DA165NP	9896-055-83901
GE	HiSpeed LX/i	DA165NP	9896-055-83901
GE	LightSpeed™ ¹⁶	DA200ULTRA	9896-055-88704
GE	LightSpeed™ Plus	DA200ULTRA	9896-055-88704
GE	LightSpeed™ Pro ³²	REEVO240G	9896-055-96401



Brand	CT Scanner Model	Dunlee Replacement Tube	Dunlee Part #
GE	LightSpeed™ QX/i	DA200ULTRA	9896-055-88704
GE	LightSpeed™ RT ¹⁶	REEVO240G	9896-055-96401
GE	LightSpeed™ Ultra	DA200ULTRA	9896-055-88704
GE	LightSpeed™ VCT	REEVO240G	9896-055-96401
GE	LightSpeed™ VCT Select	REEVO240G	9896-055-96401
GE	LightSpeed™ VCT XT	REEVO240G	9896-055-96401
GE	LightSpeed™ VCT XTe	REEVO240G	9896-055-96401
GE	LightSpeed™ Xtra	REEVO240G	9896-055-96401
GE	Optima™ CT520	DA165NP	9896-055-83901
GE	Optima™ CT520 Pro	DA165NP	9896-055-83901
GE	Optima™ CT540	DA200ULTRA	9896-055-88704
GE	Optima™ CT580	REEVO240G	9896-055-96401
GE	Optima™ CT580 RT	REEVO240G	9896-055-96401
GE	Optima™ CT580 W	REEVO240G	9896-055-96401
GE	Optima™ CT660	DA200P40 DA200P40+LMB	4598-006-53451 9890-000-87251
GE	Optima™ PET/CT 560	DA200ULTRA	9896-055-88704
GE	Optima™ PET/CT 560 FX	DA200ULTRA	9896-055-88704
GE	Revolution™ EVO	DA200P40+LMB	9890-000-87251
Siemens	Somatom Sensation 4	S532B/Akron B	9896-055-95702
Siemens	Somatom Sensation Volume Zoom	S532B/Akron B	9896-055-95702
Siemens	Somatom Sensation Volume Access	S532B/Akron B	9896-055-95702
Siemens	Somatom Sensation 10	S532Q/Akron Q	9896-055-95902
Siemens	Somatom Sensation 16	S532Q/Akron Q	9896-055-95902
Siemens	Somatom Sensation Cardiac	S532Q/Akron Q	9896-055-95902

All rights are reserved. Reproduction in whole or in part is prohibited without the prior written consent of the copyright holder.

Dunlee reserves the right to make changes in specification and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Dunlee is a brand of the Philips Company Group

How to reach us
www.dunlee.com
marketing.dunlee@philips.com

Product information
www.dunlee.com